Human Serum Albumin (HSA) 100mg/mL in Normal Saline

Catalog No. 9988

100 mL, 12 x 5 mL Kits

- For assisted reproductive procedures.
- Für assistierte Reproduktionsverfahren.
- Per tecniche di riproduzione assistita.
- Para utilización en técnicas de reproducción asistida.
- Pour les techniques de procréation médicalement assistée.
- Para técnicas de reprodução assistida.
- Για διαδικασίες υποβοηθούμενης αναπαραγωγής.
- Pro postupy asistované reprodukce.
- Til assisteret reproduktionsbehandling.
- Til assisteeret reproduktionsbehandling.
- Ar palīgīdzekļiem veicamām reproduktīvām procedūram.
- Voor geassisteerde voortplantingsprocedures.
- Do procedur wspomaganego rozrodu.
- Pentru proceduri de reproducere asistată.
- För procedurer för assisterad befruktning
- Kasutamiseks abistatud viljastamisprotseduurides.
- Asszisztált reprodukciós eljárásokhoz.
- Skirta pagalbinio apvaisinimo procedūroms.
- Yardımcı üreme işlemleri içindir.
- Na postupy asistovanej reprodukcie.
- За процедури за асистирана репродукција.
- Za postupke potpomognute oplodnje.
- Għal proċeduri ta' riproduzzjoni assistita.
- Za postopke asistirane reprodukcije.

Glossary of Symbols*:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature 2-8°C
	Do not re-use
	Do not re-sterilize
	Do not use if package is damaged
	Manufacturer
	U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	CE Mark
	Emergo Europe - Prinsessegracht 20 2514 AP The Hague The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

ENGLISH

EU Caution: For Professional Use Only.

INDICATION FOR USE

HSA is intended for use in assisted reproductive procedures which include gamete and embryo manipulation, as a supplement for culture media. Multipurpose Handling Medium (MHM) Catalog #90163 and Continuous Single Culture (CSCM) Catalog #90164.

DEVICE DESCRIPTION

HSA consists of human serum albumin from therapeutic-grade source material (100 mg/mL) in a saline solution.

COMPOSITION

<u>Salts and Ions</u>	<u>Water</u>	<u>Protein Source</u>
Sodium Chloride	WFI Quality	Human Serum Albumin

QUALITY ASSURANCE

HSA is a protein supplement which is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10⁻⁶.

Each lot of HSA is tested for:

- Endotoxin by LAL methodology
- Biocompatibility by Mouse Embryo Assay
- Sterility by the current USP Sterility Test <71>

All results are reported on a lot-specific Certificate of Analysis which is available upon request.

DIRECTIONS FOR USE

General laboratory practice includes protein supplementation when using Multipurpose Handling Medium (MHM) Catalog #90163 and Continuous Single Culture (CSCM) Catalog #90164 for gamete and embryo procedures. The amount of protein supplementation can vary among laboratories and is dependent on the phase of processing/growing the gametes and embryos.

The following are recommendations for protein supplementation based upon the corresponding product indications for use:

Warm the protein supplemented media to 37°C and equilibrate to desired pH in a CO₂ incubator prior to use.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

Multipurpose Handling Medium (MHM) Catalog #90163

For Sperm Washing:

The Human Serum Albumin (HSA) 100mg/ml solution is supplemented into Multipurpose Handling Medium (MHM) at 5 mg/ml according to Table 1 below.

Table 1: Supplementation of Multipurpose Handling Medium (MHM) at 5 mg/mL

Medium (mL)	HSA – Solution (mL)	Final Volume (mL)
9.5	0.5	10.0
19.0	1.0	20.0
28.5	1.5	30.0
38.0	2.0	40.0
47.5	2.5	50.0
57.0	3.0	60.0
66.5	3.5	70.0
76.0	4.0	80.0
85.5	4.5	90.0
95.0	5.0	100.0

For Oocyte Retrieval:

When using the HSA 100 mg/mL solution supplement the Multipurpose Handling Medium (MHM) at 5mg/mL according to Table 1 above.

For Embryo Transfer:

The HSA 100 mg/mL solution is supplemented into Multipurpose Handling Medium (MHM) at 5 mg/mL according to Table 1 above.

Continuous Single Culture (CSCM) Catalog #90164:

For Fertilization and Embryo Culture:

The HSA 100 mg/mL solution is supplemented into Continuous Single Culture (CSCM) at 5mg/mL according to Table 1 above.

STORAGE INSTRUCTIONS AND STABILITY

Store the unopened bottles refrigerated at 2° to 8°C,

Do not freeze or expose to temperatures greater than 39°C.

HSA is stable until the expiration date shown on the vial and kit label when stored as directed.

Remove desired volume using aseptic procedures.

Do not return any volume to the original container.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any bottle of medium which shows evidence of particulate matter, cloudiness or is not clear or pale yellow in color or shows any evidence of microbial contamination.

Do not use any bottle in which the sterile packaging has been compromised.

HSA is not for use as an injectable product and is antibiotic free.

To avoid problems with contamination, handle using aseptic techniques and discard any excess medium that remains in the bottle after the procedure is completed.

Information on known characteristics and technical factors that could pose a risk if the product were to be re-used have not been identified therefore the product is not to be used following the initial use of the container.

EU: Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. It is strongly recommended that every time an FUJIFILM Irvine Scientific, Inc. culture media are administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

US: Human Serum Albumin (HSA) is human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal pre-cautions.

CONTRAINDICATION

Patient is allergic to human serum albumin.

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