# Complete Human Tubal Fluid (HTF) Medium with SSS

---

Catalog No. 9922

2 x 20 mL

---

For assisted reproductive procedures.

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
Irvine**Scientific**<sup>®</sup>

## **INTENDED USE**

Complete HTF Medium with SSS is intended for use as a fertilization and culture medium for human embryos through day 3 of development.

## **PRODUCT DESCRIPTION**

Complete HTF Medium with SSS consists of HTF Medium supplemented with 10% Serum Substitute Supplement (SSS), for a final total protein concentration of 6 mg/mL. The protein component consists of 5 mg/mL therapeutic grade human serum albumin\* and 1 mg/mL alpha globulins\*. This medium contains the antibiotic gentamicin sulfate.

## **QUALITY ASSURANCE**

Complete HTF Medium with SSS is a culture medium which is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of  $10^{-3}$ .

Each lot of Complete HTF Medium with SSS is tested for:

- Endotoxin by LAL methodology
- Biocompatibility by Mouse Embryo Assay
- Sterility by the Current USP Sterility Test <71>

All results are reported on a lot specific Certificate of Analysis which is available upon request.

## **BUFFER SYSTEM**

Complete HTF Medium with SSS uses Sodium Bicarbonate as a buffering system. This is specifically designed for use in a CO<sub>2</sub> incubator.

## **DIRECTIONS FOR USE**

Complete HTF Medium with SSS should be warmed to 37°C and equilibrated to desired pH in a CO<sub>2</sub> incubator prior to use. When equilibrating in a CO<sub>2</sub> incubator, the bottle of Complete HTF Medium with SSS should be loosely capped to allow for the exchange of gas and pH equilibration.

After the medium has been pre-warmed and equilibrated in a CO<sub>2</sub> incubator, allow fertilization to occur in a culture dish that contains Complete HTF Medium with SSS. Once fertilization

occurs, the reproductive specialist should transfer the zygote/embryo into a new dish with pre-warmed equilibrated Complete HTF Medium with SSS for the growth phase. Allow embryo to grow until desired developmental stage (up to 3 days).

For additional details on the use of these products, the laboratory should consult its own standard procedures and protocols which have been specifically developed and optimized for the individual medical facility.

### **STORAGE INSTRUCTIONS AND STABILITY**

Store the unopened bottles refrigerated at 2° to 8°C.

Do not freeze or expose to temperatures greater than 39°C.

Complete HTF Medium with SSS is stable until the expiration date shown on the bottle label when stored as directed.

### **PRECAUTIONS AND WARNINGS**

Do not use any bottle of medium which shows evidence of particulate matter, cloudiness or is not reddish-orange in color.

To avoid problems with contamination, handle using aseptic techniques, and discard any excess medium that remains in the bottle after the procedure is completed.

**Complete HTF Medium with SSS contains the antibiotic Gentamicin. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.**

**\*Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies for Hepatitis B surface antigen (HBsAg), antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are non-infectious. Handle all human source material as if it were capable of transmitting infection, using universal precautions. Donors of the source material have also been screened for CJD.**

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## REFERENCES

Pool, TB, Atiee, SH and Martin, JE: Oocyte and embryo culture: Basic concepts and recent advances. In: Assisted reproductive laboratory considerations. Infertility and reproductive medicine, clinics of North America. (J. May, ed.) 9(2):181-203, 1998.

Weathersbee, PS, Pool, TB and Ord, T: Synthetic Serum Substitute (SSS): A globulin-enriched protein supplement for human embryo culture. *J Asst. Reprod and Genet.* 12:354-360, 1995.

Irvine**Scientific**<sup>®</sup>

2511 Daimler Street, Santa Ana, California 92705-5588

Telephone: 1 949 261 7800 • 1 800 437 5706

Fax: 1 949 261 6522 • [www.irvinesci.com](http://www.irvinesci.com)

PN 40636 Rev. 4