Complete P-1[®] (Preimplantation Stage One) Medium with SSS

Catalog No. 9926

2 x 20 mL

For assisted reproductive procedures.

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
Irvine**Scientific**[®]

INTENDED USE

Complete P-1[®] Medium with SSS is intended for use as a fertilization and culture medium for human embryos through day 3 of development.

PRODUCT DESCRIPTION

Complete P-1 Medium with SSS consists of P-1[®] Medium supplemented with 10% Serum Substitute Supplement (SSS), for a final total protein concentration of 6 mg/mL. The protein component consists of 5 mg/mL therapeutic grade human serum albumin* and 1 mg/mL alpha globulins*. This medium contains the antibiotic gentamicin sulfate.

QUALITY ASSURANCE

Complete P-1 Medium with SSS is a culture medium which is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10^{-3} .

Each lot of Complete P-1 Medium with SSS is tested for:

- Endotoxin by LAL methodology
- Biocompatibility by Mouse Embryo Assay
- Sterility by the Current USP Sterility Test <71>

All results are reported on a lot specific Certificate of Analysis which is available upon request.

BUFFER SYSTEM

Complete P-1 Medium with SSS uses Sodium Bicarbonate as a buffering system. This is specifically designed for use in a CO₂ incubator.

DIRECTIONS FOR USE

Complete P-1 Medium with SSS should be warmed to 37°C and equilibrated to desired pH in a CO₂ incubator prior to use. When equilibrating in a CO₂ incubator, the bottle of Complete P-1 Medium with SSS should be loosely capped to allow for the exchange of gas and pH equilibration.

After the medium has been pre-warmed and equilibrated in a CO₂ incubator, allow fertilization to occur in a culture dish that

contains Complete P-1 Medium with SSS. Once fertilization occurs, the reproductive specialist should transfer the zygote/embryo into a new dish with pre-warmed equilibrated Complete P-1 Medium with SSS for the growth phase. Allow embryo to grow until desired developmental stage (up to 3 days).

For additional details on the use of these products, the laboratory should consult its own standard procedures and protocols which have been specifically developed and optimized for the individual medical facility.

STORAGE INSTRUCTIONS AND STABILITY

Store the unopened bottles refrigerated at 2° to 8°C.

Do not freeze or expose to temperatures greater than 39°C.

Complete P-1 Medium with SSS is stable until the expiration date shown on the bottle label when stored as directed.

PRECAUTIONS AND WARNINGS

Do not use any bottle of medium which shows evidence of particulate matter, cloudiness or is not reddish-orange in color.

To avoid problems with contamination, handle using aseptic techniques, and discard any excess medium that remains in the bottle after the procedure is completed.

Complete P-1 Medium with SSS contains the antibiotic Gentamicin. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

***Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies for Hepatitis B surface antigen (HBsAg), antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are non-infectious. Handle all human source material as if it were capable of transmitting infection, using universal precautions. Donors of the source material have also been screened for CJD.**

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REFERENCES

- Conaghan, J. et. al.: Effects of pyruvate and glucose on the development of human preimplantation embryos in vitro. *Journal of Reproduction & Fertility* (1993), 87-95.
- Quinn, P. et. al.: Successful human in vitro fertilization using a modified human tubal fluid medium lacking glucose and phosphate ions. *Fertility and Sterility*, 1995; 63: 922-4.
- Carrillo, A. et al.: Improved clinical outcomes for in vitro fertilization with delay of embryo transfer from 48 to 72 after oocyte retrieval: Use of glucose-and phosphate-free media. *Fertil Steril* 1998; 69: 329-34.
- Pool, T.B. et. al.: Oocyte and Embryo Culture: Basic Concepts and Advances. *Infertility and Reproductive Medicine, Clinics of North America* 9(2): 181-203 Assisted Reproduction Laboratory Considerations (J. May, ed.) 1998.

Irvine**Scientific**[®]

2511 Daimler Street, Santa Ana, California 92705-5588

Telephone: 1 949 261 7800 • 1 800 437 5706

Fax: 1 949 261 6522 • www.irvinesci.com

PN 40637 Rev. 4