FUJIFILM

Complete MultiBlast Medium with Gentamicin and DSS

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REFERENCES

Biggers, J.D. and Racowsky, C. *The development of fertilized human ova to the blastocyst stage in KSOM[®] medium: is a two-step protocol necessary?* RBMOnline, 5:133-140, 2002.

Pool, T.B. *Recent advances in the production of viable human embryos in vitro.* RBMOnline, 4:294-302, 2002.

Biggers, J.D. *Thoughts on embryo culture conditions.* RBMOnline, 4 (suppl.1):30-38, 2001.

Lane, M., Hooper, K., and Gardner, D.K. *Effect of essential amino acids on mouse embryo viability and ammonium production.* J. Asst. Reprod. Genet. 18: 519-525, 2001

Devreker, F., Hardy, K., Van den Bergh, M., Vanin, A.S., Emiliani, S., and Englert, Y. *Amino acids promote human blastocyst development in vitro.* Hum. Reprod. 16:749-756, 2001.

Devreker, F., Van den Bergh, M., Biramane, J., Winston, R.M.L., Englert, Y., and Hardy, K. *Effects of taurine on human embryo development in vitro.* Hum. Reprod. 14: 2350-2356, 1999.

Pool, T.B., Atiee, S.H. and Martin J.E. *Oocyte and Embryo Culture: Basic Concepts and Recent Advances. Infertility and Reproductive Medicines,* Clinics of North America 9(2):181-203. Assisted Reproduction-Laboratory Considerations. (J. May, ed.) 1998.

FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706 Fax: 1 949 261 6522 • www.irvinesci.com

© 2019 FUJIFILM Irvine Scientific, Inc. All rights reserved. The FUJIFILM Irvine Scientific logo, and MultiBlast Medium are trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions. PN 40644 Rev.8

ENGLISH

EU CAUTION: For Professional Use Only.

INDICATION FOR USE

Complete MultiBlast Medium with Gentamicin and DSS[®] is intended for use in assisted reproductive procedures, which include the culture of human embryos from day 3 to the blastocyst stage of development.

DEVICE DESCRIPTION

Complete MultiBlast Medium with Gentamicin and DSS consists of MultiBlast Medium supplemented with 10% Dextran Serum Supplement (DSS), for a total macromolecule concentration of 7 mg/mL. The macromolecule components consist of 5 mg/mL therapeutic grade Human Serum Albumin (HSA) and 2 mg/mL DSS. This medium contains the antibiotic Gentamicin Sulfate (9 µg/mL).

COMPOSITION:

Salts and Ions	Amino Acids
Sodium Chloride	Alanine
Potassium Chloride	Arginine
Magnesium Sulfate	Asparagine
Calcium Chloride	Aspartic Acid
Potassium Phosphate	Glutamic Acid
Dipeptide	Glycine
Alanyl-Glutamine	Histidine
Antibiotic	Isoleucine
Gentamicin Sulfate	Leucine
pH Indicator	Lysine
Phenol Red	Methionine
Buffer	Phenylalanine
Sodium Bicarbonate	Proline
Energy Substrate	Serine
Dextran	Taurine
Glucose	Threonine
Sodium Lactate	Tryptophan
Sodium Pyruvate	Tyrosine
Protein Source	Valine
Human Serum Albumin	Cysteine
Antioxidant	Water
Sodium Citrate	WFI Quality

QUALITY ASSURANCE

Complete MultiBlast Medium is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10⁻³.

Each lot of Complete MultiBlast Medium with Gentamicin and DSS is tested for:

- Endotoxin by Limulus Amebocyte Lysate (LAL) methodology
- Biocompatibility by Mouse Embryo Assay (one-cell)
- Sterility by the current USP Sterility Test <71>

All results are reported on a lot specific Certificate of Analysis which is available upon request.

BUFFER SYSTEM

Complete MultiBlast Medium uses sodium bicarbonate as a buffering system. This is specifically designed for use in a CO₂ incubator.

DIRECTIONS FOR USE

EQUILIBRATION

Complete MultiBlast Medium should be warmed to 37°C and equilibrated to the desired pH in a CO₂ incubator prior to use.

After achieving desired developmental stage at day 3, transfer embryos to a culture dish containing pre-warmed, equilibrated Complete MultiBlast Medium. The reproductive laboratory specialist may continue growth until transfer on day 4 to 6.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

STORAGE INSTRUCTIONS AND STABILITY

Store the medium refrigerated at 2° to 8° C.

Do not freeze or expose to temperatures greater than 39°C.

For unopened bottles:

The shelf life of the product remains as the expiration date stated on bottle when stored under the recommended conditions of 2° to 8°C.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any bottle of medium which shows evidence of particulate matter, cloudiness or is not reddish-orange in color.

To avoid problems with contamination, handle using aseptic techniques and discard any excess medium that remains in the bottle or vial after the procedure is completed.

Not for injection use.

Do not use any bottle in which the sterile packaging has been compromised.

Information on known characteristics and technical factors that could pose a risk if the product were to be re-used have not been identified therefore the product is not to be used following the initial use of the container.

EU: Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. It is strongly recommended that every time FUJIFILM Irvine Scientific, Inc. Reproductive Media Products culture media are administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

US: This product contains Human Serum Albumin (HSA). Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal pre-cautions. Donors of the source material have also been screened for CJD.

CONTRAINDICATION

Product contains Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

Glossary of Symbols*:

Catalog Number

Lot Number

Sterilized using aseptic processing techniques (filtration)

Expiration: Year - Month - Day

Caution, consult accompanying documents

Consult instructions for use

Storage Temperature 2-8°C

Do not resterilize

Do not use if package is damaged

Do not use if package is damaged

Manufacturer

U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

CE Mark 0050

Emergo Europe - Prinsessegracht 20 2514 AP The Hague The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

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