## CHANG Medium MF Mitogen-Free

Catalog No. 91005

100 mL, 500 mL

For *in vitro* diagnostic use.

Zur *In-vitro*-Diagnostik.

Solo per uso diagnostico *in vitro*.

Para uso diagnóstico *in vitro*.

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
Za diagnostično uporabo *in vitro*.

### Glossary of Symbols\*:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature below -10°C
	Do not resterilize
	Do not use if package is damaged
	Manufacturer
	CE Mark
	Emergo Europe - Prinsessegracht 20 2514 AP The Hague, The Netherlands

\*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

 **FUJIFILM Irvine Scientific, Inc.**

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706 • Fax: 1 949 261 6522 • [www.irvinesci.com](http://www.irvinesci.com)

## REFERENCES

1. Arakaki, D.T. and Sparkes, R.S. (1963), *Cytogenetics*, 2, 57-60.
2. Li, J.G. et Osgood, E.E. (1949). *Blood*, 4, 670
3. Maluish, A.E. and Strong, D. M. (1986). In *Manual of Clinical Laboratory Immunology*, pp. 274-281, 3<sup>rd</sup> Edition, Eds. N.R. Rose, H. Friedman and J.L. Fahey. *American Society for Microbiology*, Washington, D.C.
4. Nowell, P. C. (1960). *Cancer Research* 20, 462-466.
5. Waithe, W. I. And Hirschron, K. (1978). In *Handbook of Experimental Immunology*, 3<sup>rd</sup> Edition, D. M. Weir, Ed. Blackwell Scientific Publications, Oxford, Chapter 26- Lymphocyte responses to activators.
6. Watt, J.L. and Stephen G.S. (1986). In *Human Cytogenetics: a practical approach*, pp. 39-55. Eds.D. E. Rooney and B.H. Czepulkowski, IRL Press Ltd., Oxford.
7. A proposed Standard System of Nomenclature of Human Mitotic Chromosomes. (1960) *Lancet*, *I*, 1063.

**INDICATION FOR USE**

CHANG Medium MF is intended for use in culturing peripheral blood and other specimens for purposes of cytogenetic analysis.

**DEVICE DESCRIPTION**

CHANG Medium MF is a mitogen-free, ready to use medium consisting of RPMI containing 20% FBS, 2 mM glutamine, 20 mM HEPES buffer and the antibiotic Gentamicin Sulfate (35 µg/mL). CHANG Medium MF may require the addition of mitogenic agents, such as phytohemagglutinin (PHA) to optimize the growth of peripheral blood and other cells. The required concentration of PHA (or other mitogens) should be determined by the individual laboratory.

**COMPONENTS**

<u>Amino Acid</u>	<u>Water</u>	<u>Salts &amp; Ions</u>
Arginine	WFI Quality	Sodium chloride
Asparagine	<u>Proteins, Hormones, and Growth Factors</u>	Choline chloride
Aspartic Acid	Fetal bovine serum (FBS)	Calcium nitrate
Cystine	<u>pH Indicator</u>	Potassium chloride
Glutamic Acid	Phenol red	Magnesium sulfate
Glutamine	<u>Other</u>	Sodium phosphate
Glycine	Biotin	<u>Antibiotic</u>
Histidine	Hydroxyproline	Gentamicin Sulfate
Isoleucine	Glutathione	<u>Vitamins and trace elements</u>
Leucine	<u>Buffers</u>	Folic acid
Lysine	Sodium bicarbonate	Nicotinamide
Methionine	HEPES	Riboflavin
Phenylalanine	<u>Energy Substrates</u>	Thiamine
Proline	Glucose	Pantothenic acid
Serine	Inositol	Cobalamin
Threonine		Pyridoxine
Tryptophan		Aminobenzoic acid
Tyrosine		
Valine		

**QUALITY ASSURANCE**

**STERILITY**

Serum used in the production of CHANG Medium MF has been tested for viral contamination per CFR Title 9 Part 113.53. It has also been screened for mycoplasma contamination. CHANG Medium MF is sterilized by filtration through a 0.1 micron filter. Samples of CHANG Medium MF are tested for possible bacteriological contamination following the sterility testing protocol described in the current USP Sterility test <71>.

Several factors including source of sample, culture conditions and selection of reagents can influence the result obtained. Users are advised to run each new batch of reagent in parallel with reference material of known suitable activity before adoption in routine use.

Each Lot of CHANG Medium MF is evaluated using peripheral blood for mitotic index, chromosomal length and quality compared to a control medium. Results are reported on a lot specific Certificate of Analysis.

**MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED**

1. Phytohemagglutinin (PHA) Solution (9 mg/mL), Sterile
2. Phenol-Free Heparin
3. Colcemid 25 µg/mL Solution
4. Potassium Chloride Solution, 75 mM
5. Acetic Alcohol, 1 part Glacial Acetic Acid: 3 parts Methanol (Analytical Reagent Grade)
6. Giemsa or 2% Acetic Acid-orcein
7. Chromic acid
8. Mountant
9. Glass Slides and Coverslips
10. Plastic Centrifuge Tubes
11. CO<sub>2</sub> Incubator
12. Bench Centrifuge
13. Vortex Mixer
14. Light Microscope

## **SPECIMEN COLLECTION AND HANDLING**

The success of whole blood cultures for cytogenetic analysis is influenced by the level of lymphocytes with normal function at the time of sampling. As this level can be affected by infection and drugs, wherever possible subjects for cytogenetic studies should not have taken drugs for 7 days prior to collection of blood for tests. Similarly, the mitotic index may be greatly reduced during the anergic phases of certain diseases (e.g. Hodgkin's disease, sarcoidosis, etc.) and, to a lesser degree, in normal individuals during the later stages of pregnancy. Preservative must not be added to blood samples for lymphocyte culture. Aseptic techniques are essential.

Blood samples should be tested without delay whenever possible. If absolutely necessary, they may be stored at 2°C to 8°C for no longer than 48 hours.

## **PREPARATION FOR USE**

CHANG Medium MF should be thawed overnight in the refrigerator (2-8°C) then gently mixed to assure homogeneity. Aseptically dispense 10 mL of medium into sterile culture flasks and equilibrate to 37°C for immediate use for bone marrow cultures.

Note: Calcium Carbonate crystals commonly form in CHANG Medium MF. The presence of these crystals has not been shown to cause any detrimental effect on product performance.

## **DIRECTIONS FOR USE**

Whole blood or separated leukocytes have been used for cytogenetic studies, but the former is the simplest and most widely used in routine studies. As with any cell culture procedure, optimal results are dependent on establishing adequate culture conditions. Since the relative content of active PHA may vary slightly between different batches, it may be beneficial to test two concentrations of PHA.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

### **A. Preparation of Whole Blood Microcultures**

1. Collect 5-20 mL of fresh blood in phenol-free heparin and mix by inversion.
2. Reconstitute PHA by adding 5 mL sterile distilled water using a sterile syringe.
3. Using aseptic technique, prepare required volume of CHANG Medium MF (microculture medium) allowing 5 mL for each of two centrifuge tubes per blood sample.
4. Dispense CHANG Medium MF to appropriately labeled sterile screw-capped centrifuge tubes and aseptically add 0.1 mL reconstituted PHA.
5. Immediately before culture, add 0.4 mL of heparinized blood using a sterile 1 mL pipette. Use 0.3 mL of blood for infants and children and 0.5 mL for pregnant women.
6. Incubate the cultures at 37°C for 72 hours. For infant and children samples, incubate one of the two cultures at 37°C for 48 hours. For pregnant women, one of the two cultures can be incubated at 37°C for 96 hours. Mix each tube by inversion daily.
7. Add 0.05 mL (50 µL)/5 mL of working solution of methotrexate to each culture 16-18 hours before addition of thymidine.
8. Add 0.1 mL (100 µL)/5 mL of working solution of thymidine to each culture after completion of methotrexate synchronization (5-6 hours before harvest).

### **B. Harvesting the Cultures**

1. Pre-warm the hypotonic solution (75 mM potassium chloride) in a 37°C water bath prior to harvest.
2. Add 0.5 µL of a solution of Colcemid to each culture.
3. Mix gently and return to the 37°C incubator.
4. Remove the cultures from the incubator after 30 minutes.
5. Centrifuge the cultures in a bench centrifuge at 1,000 rpm for 10 minutes.
6. Remove most of the supernatant fluid and discard.
7. Resuspend the pellet in 6 to 8 mL of 75 mM potassium chloride solution pre-warmed to 37°C for 10 minutes.
8. Centrifuge as in step 5 and discard supernatant.
9. Using a Pasteur pipette slowly add 6 to 8 mL of freshly prepared modified Carnoy's fixative (3 parts absolute methanol : 1 part glacial acetic acid) to the pellet while agitating constantly on a vortex mixer. Add the fixative dropwise at first, followed by a slow trickle to minimize cellular damage and formation of lumps.
10. Leave at room temperature for 10 minutes.
11. Centrifuge as in step 5 and remove the supernatant fluid as before and slowly add an additional 5 mL of fixative to resuspend pellet.
12. Repeat step 11 twice more, resuspending finally in 0.5 mL fixative. Use this cell suspension to prepare slides for examination. Care must be taken to avoid agitation of the cells.

### **C. Preparation of Slides**

1. Slides must be scrupulously clean. A suitable cleaning procedure is to soak the slides in chromic acid overnight, after which they should be washed in running water for at least half an hour and polished with a glass cloth.
2. Apply 1 or 2 drops of the resuspended cell preparation to the center of a glass slide from 3-4 inches above top of the slide and allow to spread.
3. Wipe excess fixative from the edges of the slide with filter paper.
4. At the first appearance of Newton's Rings, blow gently to speed final drying of the slide.
5. Stain with Giemsa and mount.

## **CLINICAL APPLICATIONS**

The modal human chromosome number is 46. Human chromosomes have been classified according to their length and the position of the centromere (Denver Classification). Aberrations of chromosomal constitution have been associated with a number of congenital disorders such as Down syndrome (typically with an additional small autosome) and syndromes associated with indeterminate sexuality (Turner's syndrome, Klinefelter's syndrome and others, where the sex chromosomes are found to be abnormal). An acquired chromosomal abnormality can be detected in a proportion of the leukocytes in chronic myelocytic leukemia (the "Philadelphia" chromosome) and the progress of treatment may be assessed by following this marker. As the limit of tolerance is approached in radiation therapy, there is a marked increase in the proportion of cells with bizarre chromosomal constitution and the appearance of these cells has been used as a guide to dosage.

## **STORAGE AND STABILITY**

Store CHANG Medium MF frozen, below  $-10^{\circ}\text{C}$ , until ready to use. Chang Medium MF is stable until the expiration date shown on the bottle label when stored frozen. After thawing, any unused product can be dispensed into working aliquots and refrozen (maximum of two times) for later use, or tightly capped and stored at  $2-8^{\circ}\text{C}$  for up to 30 days. Protect from fluorescent light.

## **PRECAUTIONS AND WARNINGS**

This device is intended to be used by staff trained in procedures that include the indicated application for which the device is intended. CHANG Medium MF contains FBS and should be handled with universal laboratory precautions. The medium contains an antibiotic (gentamicin sulfate), to reduce the potential of bacterial contamination, but aseptic techniques should always be used when dispensing the media. Do not use any medium that is not red in color.

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