# FUJIFILM

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### REFERENCES

Biggers, JD. and Racowsky, C. *The development of fertilized human ova to the blastocyst stage in KSOM<sup>AA</sup> medium: is a two-step protocol necessary?* RBMOnline, 5:133-140, 2002.

Pool, TB. *Recent advances in the production of viable human embryos in vitro.* RBMOnline, 4:294-302, 2002.

Biggers, JD. *Thoughts on embryo culture conditions.* RBMOnline, 4 (suppl.1):30-38, 2001.

Lane, M., Hooper, K., and Gardner, DK. *Effect of essential amino acids on mouse embryo viability and ammonium production.* J. Asst. Reprod. Genet. 18: 519-525, 2001

Biggers, JD. and McGinnis, LK. *Evidence that glucose is not always an inhibitor of mouse preimplantation development in vitro.* Hum. Reprod 16:153-163, 2001.

Devreker, F., Van den Bergh, M., Biramane, J., Winston, RML., Englert, Y., and Hardy, K. *Effects of taurine on human embryo development in vitro.* Hum. Reprod. 14: 2350-2356, 1999.

FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706 Fax: 1 949 261 6522 • www.irvinesci.com

© 2019 FUJIFILM Irvine Scientific, Inc.. All rights reserved. The FUJIFILM Irvine Scientific logo, Early Cleavage Medium, ECM and P-1 are trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions. PN 40762 Rev.8

### ENGLISH

**EU Caution:** For Professional Use Only

#### INDICATION FOR USE

Early Cleavage Medium (ECM) is intended for use in assisted reproductive procedures which include human gamete and embryo manipulation. These procedures include the use of ECM as a culture medium through day 3 of development.

#### DEVICE DESCRIPTION

ECM is a modification of P-1 Medium that contains glucose, alanyl-glutamine and EDTA, components shown to benefit human embryo development. ECM is bicarbonate based and is designed for procedures utilizing a CO<sub>2</sub> incubator. Product contains the antibiotic Gentamicin Sulfate (10 µg/mL).

#### COMPOSITION

<u>Salts and Ions</u>	<u>pH Indicator</u>
Sodium Chloride	Phenol Red
Potassium Chloride	<u>Buffer</u>
Magnesium Sulfate	Sodium Bicarbonate
Calcium Chloride	<u>Antibiotic</u>
<u>Antioxidant</u>	Gentamicin Sulfate
Sodium Citrate	<u>Dipeptide</u>
EDTA	Alanyl-Glutamine
<u>Amino Acids</u>	<u>Water</u>
Taurine	WFI Quality
<u>Energy Substrate</u>	
Sodium Pyruvate	
Sodium Lactate	
Glucose	

#### QUALITY ASSURANCE

Early Cleavage Medium is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10<sup>-3</sup>.

Each lot of ECM is tested for:

- Endotoxin by Limulus Amebocyte Lysate (LAL) methodology
- Biocompatibility by Mouse Embryo Assay (one-cell)
- Sterility by the current USP Sterility Test <71>

All results are reported on a lot specific Certificate of Analysis which is available upon request.

#### BUFFER SYSTEM

ECM uses sodium bicarbonate as a buffering system. This is specifically designed for use in a CO<sub>2</sub> incubator.

#### DIRECTIONS FOR USE

##### PROTEIN SUPPLEMENTATION

ECM does not contain protein components. General laboratory practice includes protein supplementation when using this medium for gamete and embryo culture. The amount of protein supplementation may vary among laboratories and is dependent on the phase of processing/growing the gametes and embryos. Consult your individual laboratory protocols.

##### EQUILIBRATION

ECM (supplemented with protein) should be warmed to 37°C and equilibrated to the desired pH in a CO<sub>2</sub> incubator prior to use.

After the medium has been pre-warmed and equilibrated in a CO<sub>2</sub> incubator, allow fertilization to occur in a culture dish that contains ECM supplemented with protein. Once fertilization occurs, the reproductive specialist should transfer the zygote/embryo into a new dish with pre-warmed equilibrated Early Cleavage Medium containing the desired concentration of protein for the growth phase. Allow embryo to grow until desired developmental stage (up to 3 days).

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

#### STORAGE INSTRUCTIONS AND STABILITY

Store the unopened bottles refrigerated at 2° to 8°C.

Do not freeze or expose to temperatures greater than 39°C.

Duration Following Bottle Opening:

Product should be used within 8 weeks from opening when stored under the recommended conditions of 2° to 8°C.

#### PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any bottle of medium which shows evidence of particulate matter, cloudiness or is not pale orange in color.

CAUTION: Ammonium has been shown to induce retarded embryo development in vitro, specifically in the development of the inner cell mass.

To avoid problems with contamination, handle using aseptic techniques and discard any excess medium that shows any evidence of contamination after opening.

**Not for injection use.**

#### CONTRAINDICATION

Product contains Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

#### Glossary of Symbols\*:

REF Catalog Number

LOT Lot Number

STERILE A Sterilized using aseptic processing techniques (filtration)

Expiration:  
Year - Month - Day

Caution, consult accompanying documents

Consult instructions for use

Storage Temperature  
2-8°C

Do not resterilize

Do not use if package is damaged

Manufacturer

Rx Only U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

CE Mark  
0050

EC REP Emurgo Europe - Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

\*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

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