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
Glossary of Symbols*:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature 2-8°C
	Do not resterilize
	Do not use if package is damaged
	Manufacturer
	U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	CE Mark
	Emergo Europe - Prinsessegracht 20 2514 AP The Hague The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

EU CAUTION: For Professional Use Only

INDICATION FOR USE

MHM with Gentamicin Sulfate is intended for use in assisted reproductive procedures which involve the manipulation of gametes or embryos. Specifically, MHM is indicated for use as an oocyte retrieval medium during ovarian follicle aspiration procedures (not for flushing ovarian follicles), washing sperm prior to IVF and ICSI fertilization procedures, and for transport of the embryo to the uterus during embryo transfer procedures.

DEVICE DESCRIPTION

MHM is dual buffered solution (HEPES and MOPS) that provides a safe and secure environment to maintain viability of gametes and embryos during manipulations under ambient conditions. It is a versatile solution for swim up preparation, sperm washing, oocyte retrieval and rinsing, IUI, ICSI, and embryo transfer. Product needs proteins supplement. MHM contains 10 µg/mL of the antibiotic Gentamicin Sulfate.

QUALITY ASSURANCE

MHM is a handling medium which is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10⁻³.

Each lot of MHM is tested for:

Endotoxin by Limulus Amebocyte Lysate (LAL) methodology (≤ 0.25 EU/mL)
Biocompatibility by Mouse Embryo Assay (one-cell at ≥80% expanded blastocyst 96h) .
Sterility by the current USP Sterility Test <71>
Human Sperm Survival Assay (HSSA) (≥70% motility at 24h).

All results are reported on a lot specific Certificate of Analysis which is available upon request.

COMPOSITION:

<u>Salts and Ions</u>	<u>pH Indicator</u>
Sodium Chloride	Phenol Red
Potassium Chloride	
Magnesium Sulfate	<u>Buffer</u>
Potassium Phosphate	Sodium Bicarbonate
Calcium Chloride	HEPES
	MOPS
<u>Amino Acids</u>	<u>Energy Substrate</u>
Glycine	Sodium Lactate
Taurine	Glucose
	Sodium Pyruvate
<u>Antibiotic</u>	<u>Water</u>
Gentamicin Sulfate	WFI Quality

BUFFER SYSTEM

MHM uses a buffering system composed of a HEPES (N-2-Hydroxyethylpiperazine-N'-2-ethanesulfonic acid), MOPS(3 Morpholinopropane-1-sulfonic acid) and Sodium Bicarbonate combination. This buffering system provides pH maintenance over the physiologic range (7.2 to 7.4) and does not require the use of a CO₂ incubator.

PROTEIN SUPPLEMENTATION

MHM does not contain protein components. The amount of protein supplementation may vary between laboratories and is dependent on the phase of processing/growing of gametes and embryos. Consult your individual laboratory protocols.

The following are recommendation for protein supplementation based upon the indications for use of the MHM:

For Sperm Washing:

When using FUJIFILM Irvine Scientific Inc. Human Serum Albumin (HSA) a 100 mg/mL solution use at 5 mg/mL. For 10 mL of medium, add 0.5 mL of HSA solution to 9.5 mL of the medium. When using FUJIFILM Irvine Scientific, Inc. Serum Substitute Supplement (SSS) a 50mg/mL protein solution use at 10% (v/v). For 10 mL of medium, add 1.0 mL SSS to 9.0 mL of medium.

For Oocyte Retrieval:

When using HSA a 100 mg/mL solution use at 5mg/mL. For 10 mL of medium, add 0.5 mL of HSA solution to 9.5 mL of the medium. When using SSS a 50mg/mL protein solution use at 10% (v/v). For 10 mL of medium, add 1.0 mL SSS to 9.0 mL of medium.

For Embryo Transfer:

When using HSA a 100 mg/mL solution use at 30 mg/mL. For 10 mL of medium, add 3.0 mL of HSA solution to 7.0 mL of the medium. When using SSS a 50mg/mL protein solution use at 50% (v/v). For 10 mL of medium, add 5.0 mL SSS to 5.0 mL of medium.

DIRECTIONS FOR USE

The following are general procedures for the indications of use for the MHM.

Sperm Washing:

The general procedure for washing sperm from its surrounding seminal fluid includes:

- Bring medium to room temperature or 37°C.
- Allow the semen to liquefy at room temperature for 20 to 30 minutes.
- Using aseptic techniques, transfer the liquefied semen into a sterile 10 mL conical centrifuge tube add 2 to 3 volumes of room temperature MHM (for example, a 2 mL semen sample requires 4 to 6 mL of medium). Should the volume of the sperm medium mixture be greater than 5 mL, divide into two sterile conical centrifuge tubes, minimizing the volume per tube to 4 – 6 mL, the recovery of sperm will be maximized. Samples having high viscosity may require a further processing to ensure total sperm recovery. (See the Special Processing Considerations section).
- Centrifuge the tubes at ambient temperature for 10 minutes using a force of 200 - 300 x g.
- Using a sterile pipette, remove and discard the supernatant above the "sperm pellet" by aspiration. The sperm should then be resuspended by gently flicking the tube externally with the index finger. (Note: Do not use a vortex mixer for this step). Resuspend the sperm in 1 to 2 mL of fresh medium, recap and gently mix by inversion. Samples which were fractionated for the first centrifugation step should now be recombined into one tube.
- Recentrifuge as in Step 4.
- Using a sterile pipette, remove and discard the supernatant and resuspend the sperm pellet gently by manual agitation. Add fresh medium to a final volume of 0.5 mL. The sperm are ready for assisted reproductive procedures. (Note: The total volume of the nonravid uterus is 15 - 56 mL).

SPECIAL PROCESSING CONSIDERATIONS

Processing the highly viscous semen sample:

Some samples are naturally highly viscous even after liquefaction. These samples have the consistency of heavy syrup and may be among the most difficult to process.

- After the medium is added to an ejaculate, aspirate and expel the mixture gently using an 18 gauge needle and syringe. This will "shear" some of the viscous mucous.
- Limit the amount of medium-sperm mixture from Step 1 to 5 mL per centrifuge tube for the first centrifugation step.
- If after preprocessing the sample with the needle and syringe (Step 1), the sperm do not "pellet" in a normal manner (the sperm will appear as a "cloudy fiber" attached to the bottom of the centrifuge tube), carefully aspirate as much of the supernatant as possible without disrupting the "cloudy sperm fiber" using a sterile needle and syringe. This can be done by keeping the beveled edge of the needle firmly against the wall of the centrifuge tube and slowly start aspiration from the top of the tube downward. When as much of the supernatant as possible has been removed, add 2 or 3 mL of fresh medium. Repeat the process of drawing the mixture through the 18 gauge needle and syringe. Recentrifuge the mixture. The sperm should pellet normally after the second processing.

- On subsequent sample collections, the patient should be requested to produce a split ejaculate which will minimize the viscosity in the sperm rich portion of the sample.

Oocyte Retrieval (not for flushing ovarian follicles):

MHM may be supplemented with quality tested pharmaceutical grade heparin (2.5 – 10 units/mL) to reduce clotting of the follicular aspirates containing blood.

- Bring protein supplemented medium to room temperature or 37°C.
- The collected follicle aspirates should be transferred into an empty sterile dish.
- Identify the oocytes and remove them from follicular fluid and possible blood contamination using sterile pipettes using pre-rinsed and supplemented MHM.
- Rinse the oocytes in warmed and supplemented MHM.
- Place oocytes into an equilibrated culture medium for further handling.

Embryo Transfer:

Transfer of embryos from culture medium on day 3 or day 5:

- On day 3 or day 5 following assessment of embryos for development, bring protein supplemented medium to room temperature or 37°C.
- Set up one sterile wash dish containing pre-warmed protein supplemented MHM for each set of embryos.
- Place 1.0 mL of the pre-warmed protein supplemented MHM into the well of a sterile 1-well dish.
- Place the wash dish onto a heated stage.
- Wash the embryos in the wash dish by picking up the embryos 2 - 3 times and moving them around in minimal volume of the pre-warmed protein supplemented MHM within the well.
- After washing the embryos are ready for transfer into the patient.

For additional details on the use of MHM each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

STORAGE INSTRUCTIONS AND STABILITY

Store the unopened bottles refrigerated at 2° to 8°C.

Do not freeze or expose to temperatures greater than 39°C.

Duration Following Bottle Opening:

Product should be used within (5) weeks from opening when stored under the recommended conditions of 2° to 8°C.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended. This device is not for use in ovarian follicles flushing procedure. This media is not for use in oocyte flushing procedures.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable

Do not use any bottle of medium which shows evidence of particulate matter or cloudiness.

MHM should be tightly capped when used in a CO₂ incubator to avoid pH levels of 7.0 or less.

To avoid problems with contamination, handle using aseptic techniques and discard any excess medium that remains in the bottle or vial after the procedure is completed.

CONTRAINDICATION

Product contains Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

 FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706 Fax: 1 949 261 6522 • www.irvinesci.com

© 2019 FUJIFILM Irvine Scientific, Inc. All rights reserved. The FUJIFILM Irvine Scientific logo, Multipurpose Handling Medium, MHM, Serum Substitute Supplement, and SSS are trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions. PN 40933 Rev.7

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