## CHANG Amnio with Gentamicin and L-Glutamine

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Za diagnostično uporabo *in vitro*.

### Glossary of Symbols\*:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature store below -10°C
	Do Not Re-Sterilize
	Do not use if package is damaged
	Manufacturer
	CE Mark
	Emergo Europe - Prinsessegracht 20 2514 AP The Hague The Netherlands

\*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

### ENGLISH

#### INDICATION FOR USE

CHANG Amnio with Gentamicin and L-Glutamine may be used for the following applications:

1. the primary culture of amniotic fluid cells
2. growing passaged amniotic fluid cells
3. solid amnion tissue from chorionic villi sampling.

This medium has been designed for use in CO<sub>2</sub> incubators (cultures equilibrated with 5%-8% CO<sub>2</sub> atmosphere).

The final pH must be 6.65-7.44 Please see Directions for Use.

#### DEVICE DESCRIPTION

CHANG Amnio is a complete, ready-to-use medium for the primary culture of human amniotic fluid cells (AFC), chorionic villus sampling (CVS), and products of conception (POC) for use in karyotyping and other prenatal genetic testing. It has been optimized for both flask and *in situ* methodologies. This product contains the antibiotic Gentamicin Sulfate (50 µg/mL)

#### COMPONENTS

<u>Buffers</u> Sodium bicarbonate	<u>Proteins, Hormones, and Growth Factors</u> Fetal bovine serum (FBS) Newborn bovine serum Human transferrin Fibroblast growth factor (FGF) Insulin Progesterone Testosterone Beta estradiol Hydrocortisone	Magnesium sulfate Sodium phosphate
<u>Antioxidant</u> Thioctic acid	<u>Nucleic acids</u> Deoxyadenosine Deoxycytidine Deoxyguanosine Adenosine Cytidine Guanosine Thymidine Uridine	
<u>Antibiotic</u> Gentamicin Sulfate	<u>Other</u> Ethyl alcohol Thyronine	
<u>Amino Acid</u> Alanine Arginine Asparagine Aspartic Acid Cysteine Glutamic Acid Glutamine Glycine Histidine Isoleucine Leucine Lysine Methionine Phenylalanine Proline Serine Threonine Tryptophan Tyrosine Valine	<u>pH Indicator</u> Phenol red	<u>Vitamins and trace elements</u> Ascorbic acid Folic acid Nicotinamide Riboflavin Thiamine Pantothenic acid Cobalamin Pyridoxine Biotin
	<u>Energy Substrates</u> Glucose Pyruvate Inositol	<u>Water</u> WFI Quality

#### QUALITY ASSURANCE

**STERILITY**  
Serum used in the production of CHANG Amnio has been tested for viral contamination per CFR Title 9 Part 113.53. It has also been screened for mycoplasma contamination. CHANG Amnio is sterilized by filtration through a 0.1 micron filter. Samples of CHANG Amnio are tested for possible bacteriological contamination following the sterility testing protocol described in the current USP Sterility test <71>.

#### PREPARATION FOR USE

Thaw rapidly by swirling bottle in a 37°C water bath.

CHANG Amnio contains gentamicin (50 mg/L). Additional antibiotics may be added if desired.

#### ALIQUOTING CHANG Amnio

1. Thaw CHANG Amnio according to the above instructions.
2. Distribute aseptically into convenient sized aliquots and refreeze.
3. Thaw aliquots in 37°C water bath when ready to use.

#### DIRECTIONS FOR USE

Use of CHANG Amnio for Primary Cultures: *in situ* Methodologies

1. Centrifuge amniotic fluid at approx. 1,200 rpm for 10 minutes to concentrate the cells.
2. Aspirate supernatant from the centrifuged tube, leaving approx. 0.5 mL above cell pellet (or about 2x volume of pellet) of spun amniotic fluid. Aliquot supernatant (at least 1 mL, if possible) for alpha-fetoprotein (AFP)

1. acetyl cholinesterase assays, if necessary. If specimen is bloody, prepare an additional aliquot for further testing. Resuspend cell pellet in small volume of patient's own amniotic fluid. Add sufficient CHANG Amnio to the concentrated cell suspension to allow for final plating volume of 0.5 mL per coverslip (total of 4 coverslips, depending on size of cell pellet) or 2 mL per flaskette. If specimen is received from a patient in the third trimester of pregnancy, the pellet may be larger but contain less viable cells, thus requiring heavier seeding (less media than normal).
2. Incubate cultures undisturbed at 37°C 5%-8% CO<sub>2</sub> atmosphere.
3. Flood cultures on day 2 by adding 2 mL of CHANG Amnio.
4. After 4 to 5 days, cultures should be checked for growth. Cultures should be fed once growth has been observed. Feed cultures by removing all of the culture supernatant and replacing with 2 mL fresh CHANG Amnio. It is recommended that cultures be fed every 2 days thereafter. For bloody specimens, cultures may require more frequent media changes.
5. Check cultures for growth on, or after, day 5, and harvest when sufficient colonies are observed.
6. Best results obtained when cultures are fed with CHANG Amnio the day before the harvest.

#### Use of CHANG Amnio for Primary Cultures: Flask Methodologies

1. Centrifuge amniotic fluid at approx. 1,200 rpm for 10 minutes to concentrate the cells.
2. Aspirate supernatant from the centrifuged tube, leaving approx. 0.5 mL above cell pellet (or about 2x volume of pellet) of spun amniotic fluid. Aliquot supernatant (at least 1 mL, if possible) for alpha-fetoprotein (AFP) and acetyl cholinesterase assays, if necessary. If specimen is bloody, prepare an additional aliquot for further testing. Resuspend cell pellet in small volume of patient's own amniotic fluid. Add 4 mL of CHANG Amnio for a total volume of 5 mL per flask. If specimen is received from a patient in the third trimester of pregnancy, the pellet may be larger but contain less viable cells, thus requiring heavier seeding (less media than normal).
3. Incubate cultures undisturbed at 37°C, 5%-8% CO<sub>2</sub> atmosphere.
4. Check for growth on day 5. Change medium with 2 mL of fresh CHANG Amnio and harvest if sufficient cell growth is observed.
5. Check cultures for growth and completely change medium every day thereafter until sufficient colonies are observed and are ready to harvest. For bloody specimens, cultures may require more frequent media changes.
6. Best results obtained when cultures are fed with CHANG Amnio the day before the harvest.

#### Use of CHANG Amnio for Growing Passaged Amniotic Fluid Cells:

To passage the cells, treat cultures with trypsin (or pronase, etc.) as you would normally do when cells are grown in conventional medium. However, protease treatment should be carefully monitored. Amniotic fluid cells grown in CHANG Amnio tend to be more sensitive to protease treatment than when grown in conventional medium. It may be necessary to modify your protocol to take this into account.

Note: The pH of the medium used to feed the cultures must be between 6.65 -7.44 (i.e. the medium must be slightly yellowish salmon color). pH can easily be adjusted by placing the medium in a 5%-8% CO<sub>2</sub> incubator with the cap slightly loosened for about 30 minutes.

#### Use of CHANG Amnio for dissection of CVS:

1. Obtain one 100 mm X 15 mm petri dish divided into 2 partitions and one 60 mm X 15 mm petri dish. Dispense 2 mL of media (without serum) in the smaller petri dish. Dispense 4 mL of the media in one partition of the large petri dish and 6 mL of media in the other partition of the same dish.
2. Carefully aspirate media and villi from the original centrifuge tube containing the patient sample. Dispense media and villi into the partition of the larger

1. petri dish with 4 mL of aspiration media.
2. Using an inverted microscope and two sterile forceps, remove blood clots and any maternal decidua present from the exterior of the chorionic villi. Be careful to avoid damaging the fragile villi. Transfer the cleaned villi to the 6 mL partition of the larger petri dish.
3. Perform a final cleaning, using forceps to grasp villi and agitate gently while immersed in the media to remove any excess decidua, blood clots or debris. Choose villi with visible branches and veins, if possible. Determine the amount of villi present to prepare the optimal number of cultures (5 mg is ideal amount to use per culture; don't set up more than 20 mg of villi).
4. Place cleaned villi into smaller petri dish with 2 mL of media (without serum) using forceps and then transfer villi and media to a 15 mL centrifuge tube.

#### Use of CHANG Amnio for Primary Cultures of CVS

1. Add 4 drops of antibiotic (ie. Gentamicin Sulfate, 50 µg/mL) to the centrifuge tube and let sit for 30 minutes.
2. Centrifuge villi at approx. 1,400 rpm for 5 minutes.
3. Aspirate supernatant from the centrifuged tube, leaving 0.5 mL of media above cell pellet (or about 2x volume of pellet).
4. Gently resuspend pellet. Add 2 mL of CHANG Amnio media to centrifuge tube.
5. Repeat steps 2-3. Add 2 mL of trypsin and incubate culture undisturbed at 37°C, 5%-8% CO<sub>2</sub> atmosphere for 10 minutes. Remove tube from incubator, resuspend pellet and place in incubator for 10 additional minutes.
6. Remove centrifuge tube from incubator, resuspend pellet and centrifuge for 8-10 minutes at 1,400 rpm.
7. Aspirate supernatant from the centrifuged tube. Resuspend pellet, then add 1 mL of collagenase to tube and place in incubator for 5 minutes.
8. Remove from incubator and visually check to see if pellet is cloudy and no distinct individual pieces of villi can be seen. If pellet is not cloudy, place back in incubator for 5 minutes.
9. Repeat step 8 until pellet is cloudy. Add 3 mL of CHANG Amnio to centrifuge tube to stop the action of the collagenase.
10. Centrifuge tube for 8-10 minutes at 1,400 rpm.
11. Aspirate supernatant from the centrifuged tube, leaving 0.5 mL of media above cell pellet (or about 2x volume of pellet). Resuspend pellet before adding CHANG Amnio used for set up.
12. Set up optimal number of cultures (approx. 1 culture per every 5 mg of villi used) using 0.5 mL of CHANG Amnio per culture for each petri dish containing a coverslip.
13. Incubate cultures undisturbed at 37°C 5%-8% CO<sub>2</sub> atmosphere.
14. Flood cultures on day 2 by adding 1.5 mL of CHANG Amnio.
15. At 4 days, cultures should be checked for growth. If growth is observed, remove media and add 2 mL of fresh CHANG Amnio to each coverslip. Cultures should be fed every 2 days thereafter. For bloody specimens, cultures may require more frequent media changes.
16. Check cultures for growth on day 5, and harvest when sufficient colonies are observed.
17. Best results obtained when cultures are fed with CHANG Amnio the day before the harvest.

#### STORAGE AND STABILITY

Store frozen below -10°C. Product is stable until the expiration date on the bottle label when stored frozen. Unused product can be dispensed into working aliquots and refrozen for later use, or tightly capped and stored at 2°C to 8°C for up to 30 days; it may be frozen a maximum of two times. Protect from fluorescent light

#### PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in procedures that include the indicated application for which the device is intended. Do not use any bottle in which the sterile packaging has been compromised. Do not use CHANG Amnio beyond the expiration date indicated on the label.

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