

Insemination Cannula – Standard model Insemination Cannula – Small model Genitor®

Warning

Read the complete directions before use.

- Failure to properly follow the instructions, warnings and cautions may lead to serious surgical consequences or injury to the patient.
- These procedures should only be performed by persons having adequate training and familiarity with these techniques.
- Consult medical literature regarding techniques, complications and hazards prior to performance of these procedures.
- To be used by or under the direction of qualified persons in line with local guidelines governing Assisted Reproductive Techniques, if applicable.

STERILE: Contents sterile unless package has been opened or damaged. Discard if product or packaging is damaged.

Description

Order Code	IUI-S	IUI-G
Designation	Insemination Cannula Standard model	Insemination Cannula Small model Genitor®
Units per box	Box of 25 units	
Internal volume	0.40 ml	0.10 ml
Internal diameter	1.65 mm	0.65 mm
External diameter	2.1 mm	1.60 mm
Total length	17 cm	
Length of the flexible tip	4.5 cm	
Hysterometry guide-marks	5.5 / 6.5 cm	
Specific feature	2 distal openings	
Components	Polyethylene	

Catheter with a very flexible tip in transparent polyethylene.
Standardized connector to connect the tube to a syringe. The catheter and its very flexible blunt-end tip provide access to the uterine cavity without trauma. The sperm preparation is released near the tubal ostia through the two distal openings on the tip.

Class I medical device, complying with Directive 2007/47/CE.
CE 0120 Marking.
U.S. Federal law restricts this device for sale by or on the order of a physician.

Single use: Discard after single use.
Latex free.
Individually sterile package.
Sterilized by irradiation. Do not resterilize.

Indications

- Intrauterine artificial insemination with prepared sperm.
- The catheter helps reach the uterine cavity without trauma and release the sperm preparation close to the tubal ostia.

Contraindications

The cannula should not be used:

- in the presence of chronic cervical infection
- in the presence of or after recent pelvic inflammatory disease
- for intra-fallopian procedures
- for embryo transfers

Caution

- Prior to use, all devices coming in contact with gametes should be checked for integrity and rinsed with appropriate biological media.
- The cannula should never be forced against digitally felt resistance while inside the uterine cavity, as this may result in damage to the endometrial tissue and bleeding.

Instructions for use

- Connect the cannula to a syringe and rinse it with the appropriate biological medium.
- Load with the sperm preparation and purge to expel air.
- Insert a speculum.
- Wipe the cervix carefully.
- If necessary, place the cervix in the uterine axis using Pozzi forceps to help the cannula introduction.
- A slight rotation of the cannula helps overcome certain cervical obstacles.
- Once the endocervical orifice has been cleared, gently push the cannula tip up to the insemination site.
- Expel the sperm preparation very slowly 1 cm from the tubal ostia.
- Remove the cannula.