

# Needle Ovo Art of Prodimed

## Warning

Read the complete directions before use.

- Failure to properly follow the instructions, warnings and cautions may lead to serious surgical consequences or injury to the patient.
- These procedures should only be performed by persons having adequate training and familiarity with these techniques.
- Consult medical literature regarding techniques, complications and hazards prior to performance of these procedures.
- To be used by or under the direction of qualified persons in line with local guidelines governing Assisted Reproductive Techniques, if applicable.

STERILE: Contents sterile unless package has been opened or damaged. Discard if product or packaging is damaged.

## Description

Simple flow, triple bevel needle, with a 5 mm laser marked echogenic zone before the beveled tip. The needle is fitted with a Luer cone to connect a syringe.

Order code	<b>AOVO</b>
Designation	<b>Needle Ovo Art of Prodimed</b>
Units per box	10
Type of needle	Triple bevel - Simple flow
Usable length	30 cm
External diameter	1.4 mm
Internal diameter	1 mm
Echogenic zone	5 mm laser mark
Connection	Luer-lock

Class IIa medical device, complying with Directive 2007/47/CE.

CE 0120 Marking.

U.S. Federal law restricts this device for sale by or on the order of a physician.

Single use: Discard after single use.

Latex free.

Individually sterile package.

Sterilized by irradiation. Do not resterilize.

**Indications**

- In Vitro Fertilization (IVF) procedure: Percoelioscopy or ultra-sound guided (transvesical, transvaginal, per urethral and endovaginal) oocyte retrieval.
- Ovarian cyst aspiration.

**Contraindications**

The device should not be used:

- in the presence of chronic cervical infection
- in the presence of or after recent pelvic inflammatory disease

**Caution**

Prior to use:

- All devices coming in contact with gametes should be checked for integrity and rinsed with appropriate biological media.
- The system must be tested to ensure all connections are sound and flow can be achieved.

**Instructions for use**

- Check the integrity of the device (bevel sharpness, connection part).
- Rinse the needle and syringe set with a biological medium adapted to the collection.
- Depending on the approach chosen, perform the fine needle aspiration using ultrasound guidance or coelioscopy.
- At the end of the procedure, rinse the needle and syringe to make sure there is no oocyte left inside.
- Discard the needle in accordance with local medical hazardous waste practices.