FUJIFILM

Wako

Code No. 121-05063 (20 µg)
125-05061 (20 µg × 5)

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Lysyl Endopeptidase, originally isolated from the soil bacterium by Masaki *et al.*, cleaves specifically the peptides on the carboxy-terminal side of Lysine residues. This enzyme is very useful both in protein sequence analysis and in enzymatic synthesis of Lys-X compounds.

Source	: Bacteria (See Notice.)
Appearance	: Lyophilized form containing 2 mmol/L Tris-HCl buffer, pH 8.0
Activity	: Indicated on the label
Molecular Weight	: 27,000 (Gel filtration), 30,000 (SDS-PAGE)
Solubility	: Soluble in water or buffer solutions.
Stability	: Stable at 4°C for at least 24 hours, when dissolved in Tris buffer at pH 4 - 11 ⁵⁾ . Stable at 30°C in the range of pH 6 - 11, but unstable over 50°C.
Optimal pH	: 9.0 ~ 9.5 (Amidase activity)
Isoelectric point	: 6.9 ~ 7.0
Substrate specificity	: Hydrolysable substrate ···· Tos-Lys-OMe, Bz-Lys-NH ₂ , Bz-Lys-pNA, Lys-pNA Unhydrolysable substrate ···· Bz-Arg-NH ₂ , Bz-Arg- pNA, Arg-pNA
Inhibitors	: DFP, PMSF, TLCK

Assay method

1. Reagents

- A. 0.2 mol/L AMP buffer, pH 9.5
Dissolve 4.2 g of 2-Amino-2-methyl-1,3-propanediol in 150 mL of water, adjust to pH 9.5 with 1 mol/L HCl, and then add water to bring the volume up to 200 mL.
- B. 2.5 mmol/L Substrate solution
Dissolve 22.6 mg of *N*^α-Benzoyl-DL-lysine-*p*-nitroanilide Hydrobromide in 20 mL of water.

- 1/6 -

- C. 2 mmol/L Tris-HCl buffer, pH 8.0

Dissolve 0.24 g of 2-Amino-2-hydroxymethyl-1,3-propanediol in 900 mL of water, adjust to pH 8.0 with 0.1 mol/L HCl, and then add water to bring the volume up to 1 L.

- D. Enzyme solution

Dissolve 1 vial of Lysyl Endopeptidase in 1 mL of Reagent C added directly into the vial.

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2. Procedure

Reagent	Test	Blank
A	2.6 mL	2.6 mL
B	0.3 mL	0.3 mL
Pre-incubate at 30°C for 5 minutes.		
D	0.1 mL	—
C	—	0.1 mL
Immediately mix, and incubate 30°C for exactly 25 minutes		
E	1.0 mL	1.0 mL

Immediately, measure the absorbance at 405 nm of wavelength with water as the control.

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- 2/6 -

3. Unit definition

One amidase unit (AU) is the amount of enzyme, which will produce 1 μ mol of *p*-nitroaniline per minute at 30°C, pH 9.5.

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[Protocols for In-Gel Digestion]

Use a siliconized microcentrifuge tube and Pipette-Tip to prevent trapping any proteins. Use a gel staining kits for Mass Spectrometry, such as Wako's Silver Stain MS Kit (Code No. 299-58901) and Negative Gel Stain MS Kit (Code No. 293-57701).

1. Separate protein samples by electrophoresis.
2. Cut the protein bands from the gel and slice them into a microcentrifuge tube.
3. Destain the gel. (Use Destaining solution which is included in the gel staining kit for Mass Spectrometry)
4. Add 300 μ L Acetonitrile(ACN) to the tube. Shake the gel pieces for 30 minutes with a mixer (Dehydration).
5. Remove the ACN and cover the microcentrifuge tube with Parafilm.
6. Make some hole on the parafilm and vacuum dry for 15 minutes.
7. Reduce the proteins of the gel pieces with 100 μ L of 10 mmol/L Dithiothreitol (DTT) in 100 mmol/L Ammonium Bicarbonate and incubate for 1 hour at 56°C.
8. After cooling to room temperature, replace the DTT solution with the same volume of 50 mmol/L Iodoacetamide in 100 mmol/L Ammonium Bicarbonate and incubate the gel pieces for 45 minutes in the dark with occasional vortexing.
9. Wash the gel pieces for 10 minutes with 100 μ L of 100 mmol/L Ammonium Bicarbonate.

10. Dehydrate the gel pieces for 15 minutes with 300 μ L ACN.
11. Swell the gel pieces for 15 minutes with 100 μ L of 100 mmol/L Ammonium Bicarbonate.
12. Dehydrate again the gel pieces for 15 minutes with 300 μ L ACN.
13. Remove the liquid phase, and vacuum dry the gel pieces for 15 minutes.
14. Swell the gel pieces for 45 minutes with 100 μ L Lysyl Endopeptidase solution* in an ice-cold bath.
※ Dilute the Lysyl Endopeptidase in 50 mmol/L Tris-HCl, pH 8.5 to 10 μ g/mL.
15. Remove the Lysyl Endopeptidase solution, and incubate the gel pieces overnight at 37°C with 10 μ L of 50 mmol/L Tris-HCl, pH 8.5.
16. Shake the gel pieces for 20 minutes with 50 μ L of 20 mmol/L Ammonium Bicarbonate and extract the peptides.
17. Shake the gel pieces three times for 20 minutes with 5% Formic acid in 50% ACN and extract the peptides.
18. If necessary, concentrate the peptides with Speed Vac.
19. Desalt and purify the peptides with ZipTip®
20. If necessary, concentrate the peptides to 2 μ L with weak vacuum.
21. Add the matrix and analyze it by mass spectroscopy.

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Notice

The source of this product had been indicated as “*Achromobacter lyticus*”, based on the physiological and morphological properties of the bacteria. However, by recent bacterial taxonomy, the bacteria were identified as *Lysobacter enzymogenes*.

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* : **Trans** is the Japanese translation.

Lysyl Endopeptidase[®] is registered trademark of FUJIFILM Wako Pure Chemical Corporation.

ZipTip[®] is registered trademark of Millipore.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>