


Simplified Warming Protocol for HSV Device

Oocytes and Embryos

STEPS 9–12 MUST BE PERFORMED AT ROOM TEMPERATURE (22–27°C)

Do not begin warming procedure until you have a pre-equilibrated dish of appropriate culture medium supplemented with SSS or DSS at 20% (v/v), or HSA at 12 mg/mL for final recovery of specimen(s).

Have all necessary materials, tools and equipment ready and easily accessible before starting procedure.

- Set up thawing dishes (as shown in diagram):
 - At 37°C:** Aseptically dispense a minimum volume of 250 μ L of TS and warm to 37°C (incubator without CO₂) or on a heating stage at least 30 minutes prior to starting warming procedure
 -  For oocytes, dispense a minimum of 1 mL of TS.
- Identify HSV Straw(s) to be warmed from LN₂ storage and quickly transfer to a LN₂ filled holding reservoir in preparation for warming procedure.
- Place LN₂ filled holding reservoir in close proximity to the working area and stage of the microscope in order to achieve subsequent rapid manipulation from reservoir to TS.
- Remove TS dish from 37°C incubator or heating stage and place it under focus on top of the microscope stage.
- Lift the straw with forceps enough to expose the colored handling rod. **Make sure the end with the specimen(s) remains immersed in the LN₂.**
- Use a Knipex (or other wire cutter device) to cut the straw at the height of the colored handling rod. The red cut-length guide on the Knipex should be positioned in maximum length position or removed.
 - Alternatively, use fingers and thumb to spin the straw while making cutting movements with scissors, 10 mm under the top of the colored handling rod.
- With one swift but controlled motion, quickly grab the handling rod and extract it out of the straw.
- Immediately plunge the curved spatula (or gutter) of the handling rod into the 37°C TS and gently swirl to detach specimens from device and leave the oocyte or embryo for a total of 1 minute in the TS.


After 30 seconds following the initial plunge, gently pipette the specimen (if floating) and place at the bottom of the TS drop/well.

Steps 9-12 must be performed at room temperature (22-27°C).

- At room temperature:** Aseptically dispense one (1) 50 μ L drop of DS on a sterile Petri dish
- Transfer specimen(s) to DS for 4 minutes. Gently pipette specimens once to ensure complete rinse in DS.
 - During the 4 minute exposure in DS, aseptically dispense two (2) 50 μ L drops of WS (WS1, WS2) as shown in diagram.
 - Transfer specimen(s) to each WS1 then WS2 for 4 minutes each, undisturbed.
 - Transfer warmed OOCYTE(S) to pre-equilibrated culture medium with 20% (v/v) protein supplement or 12 mg/mL for recovery (2–3 hours to allow time for spindle re-formation) prior to subsequent manipulations.



There are two options for warmed EMBRYO(S):

- For immediate transfer to patient: transfer EMBRYO(S) to pre-equilibrated 'transfer' medium containing 20% (v/v) protein supplement or 12 mg/mL.
- For further culture: transfer EMBRYO(S) to pre-equilibrated culture medium containing 20% (v/v) protein supplement or 12 mg/mL for a 4-hour recovery period. After recovery period, transfer EMBRYO(S) to culture medium with 10% (v/v) protein and incubate accordingly until desired developmental stage has been reached for transfer to patient.



Minimum 1 mL of TS for oocytes in a 4-well dish

Transfer to culture medium with 20% (v/v) protein for recovery

KEY	
TS	Thawing Solution
DS	Dilution Solution
WS	Washing Solution
	Transfer specimen to next drop
	HSV Device

Simplified Warming Protocol for HSV Device

Tips

- Have all necessary materials, tools and equipment ready and easily accessible before starting procedure.
 - Refer to Directions for Use or Product Insert that accompany the HSV Device for specific instructions.
 - Pre-equilibrate a dish of appropriate culture medium with 20% (v/v), or 12 mg/mL protein for final recovery of specimen(s).
 - HSV Device must remain submerged in LN₂ until ready to warm. When transferring HSV Device from LN₂ filled holding reservoir, or between LN₂ storage tanks, vitrification devices should always be submerged in a LN₂ filled goblet to prevent uncontrolled/premature warming in air.
 - Set up warming dish with a minimum volume of 250 µL of TS for embryos, or a minimum volume of 1 mL of TS for oocytes, at least 30 minutes prior to beginning procedure.
 - Remove the TS dish from 37°C incubator or heating stage right before the warming procedure.
 - Immediately immerse the curved spatula (or gutter) into the 37°C TS after removed from LN₂ (within 2 seconds).
 - Maintain proper function of your Knipex (or other wire cutter device) with regular monitoring and replacing when needed.
 - Gently swirl the gutter in TS to help detach the specimen(s) from the device.
 - Keep focusing under the microscope while immersing gutter into TS to locate the specimen(s).
 - Limit exposure to light while moving the specimens through the solutions.
 - If you choose to, you may overlay DS and WS drops with 8.0 mL to 8.5 mL of equilibrated mineral oil, at least 45 minutes prior to starting warming procedure.
- 📄 *Following complete recovery (2–4 hours post-warming), oocytes must be fertilized by ICSI due to zona hardening during vitrification.*

www.irvinesci.com

FUJIFILM
Value from Innovation



FUJIFILM IRVINE SCIENTIFIC – CORPORATE

1830 E Warner Avenue, Santa Ana, CA 92705 USA

Phone: 1 (949) 261-7800

Toll Free: 1 (800) 437-5706

Fax: 1 (949) 261-6522

Support: tmrequest@irvinesci.com