

Michael De La Cruz* and Tom Fletcher

Department of Research and Development, FUJIFILM Irvine Scientific, 1830 E. Warner Ave., Santa Ana, CA 92705

*michael.de.lacruz@fujifilm.com

Introduction

Three Manufacturing Science and Technology (MSAT) Service Case Studies

- Company A**
- Simplify hydration method of production medium for continuous process
- Company B**
- Investigate potential sources of precipitate in concentrated feed medium
- Company C**
- Investigate if factors in preparation method cause the aberrant process performance


MSAT Service goals are:

- Develop user-friendly cell culture media and other process solutions for late-stage drug substance manufacturing, and cell and gene therapy process development. (focus shifts from speed to increasing efficiency and reducing risk)
- Improve media process workflows by:
 - Simplifying preparation methods
 - Reducing errors associated with avoidable adjustments, handling, and weighing
- Reduce cost by improving the reliability of the media preparation unit operation increasing the overall efficiency in the facility

Manufacturers' comments specific powder-media associated problems to be avoided:

- Mixing and Handling Problems/Needs**
- I would like to avoid the **long mixing times** required to mix the components.
- Consistency Problems**
- Variation** in powders is normal. We see differences [inconsistencies] between the media we purchase for GMP manufacturing and non-GMP research... **grinding** is likely done at different scale, likely using different mixing equipment.
- Lack of Industry Information and Support**
- Media companies are **not adding value**. All they do is fill our orders. Never any investigation or discussion of the powders we ask them to make. They never suggest how we could be doing things better.
- Need for Closed Systems**
- Regulators want more **closed processing** or **cleaner air** in media prep. If all powder media components were mixed together without **exposure to the environment**, these would be preferred and a big advantage.

MSAT Framework



Applied Design Principles

- Simplicity
- Risk Reduction
- Efficiency
- Risk Mitigation
- Cost Reduction
- Maintained product quality attributes
- Maintained productivity

Challenges faced by End Users

- Powder hydration issues:
- Solubility
 - Component Interactions
 - Complexation
 - Reduction/Oxidation
- Error due to human intervention from:
- Handling and weighing
 - Avoidable adjustments to volume(q.s.), pH, and osmolality

Abstract

Manufacturing Science and Technology (MSAT) Service was created to develop user-friendly cell culture media and other process solutions for the late stages of drug substance and cell and gene therapy process development when customers' efforts shift from speed to market to increasing efficiency and reducing risk. These three case studies presented describe how the MSAT service collaborates with FISI's customers to improve their cell culture media process workflows by simplifying their preparation methods and rebalancing the powder formulas to reduce errors associated with avoidable adjustments, handling, and weighing. Overall, the MSAT service improved the reliability of customer's media and preparation methods, resulting in significant reduction in labor and cost.

Case Study A

Goal: Simplify cell culture media preparation method

A. Complex/Insoluble ABC Blend Simple/Soluble

B.

C.

- Combined components to balance the formula
- Engineered formula for solubility and reliable, easy preparation
- Reduced intervention
- Reduced time and labor
- Designed for scale-up

The total preparation time was reduced from 4 hours to around 1 hour

- Final production medium powder composition:**
- 45% Growth medium
 - 20% Feed medium
 - 35% Supplement
 - 100x liquid supplement
- Simplified medium and preparation method**
- Single powder
 - No adjustments to pH or osmolality during hydration
 - No q.s. step
 - Glutamine and NaHCO₃ added as separate powders




Customer better prepared for adopting automated Oceo Rover™ methods – based on the same basic chemistry.

Received positive customer feedback: "The media preparation was a wonderful experience."

- Customer stated:**
- "We managed to complete 10L preparation and filtering within an hour. In contrast to the 3 to 4 hour process of manually weighing individual component."
 - "The pH and osmolality fall within right range without adjusting."
 - "Metabolite measurement are close to the one manually prepared."
 - Finally, our customer "confirmed that the cells are growing well with the MSAT medium."

Three liter medium preparation




- pH shifts from ~5 to 7 after the addition of the 100x liquid supplement, leading to full powder dissolution
- Targeted pH is achieved after sodium bicarbonate is added
- No q.s. step required by starting with the total required volume


Case Study B

Goal: Identify the cause of precipitation and implement reliable preventative method

Feed medium stability failing the required 12-day storage hold. Observed in three powder qualification batches.



- Batch record data
- Method
- Equipment
- Parameters
- Time
- Environment
- Raw materials
- Sources of materials



The powder was slightly darker near the surface and the storage container emitted sulfur-like odor

Side-by-side EMS powders screen

Filter Membranes

MSAT TEST	MSAT TEST	MSAT TEST	MSAT TEST
981800R210004	981800R210001	981800R210003	981800R210002
205 (No)	204 (Yes)	205 (No)	204 (Yes)
2 (Yes)	3 (Yes)	3 (Yes)	2 (Yes)
6 (Yes)	6 (Yes)	6 (Yes)	6 (Yes)
None (30+)	None (30+)	None (30+)	None (30+)
38	38	38	38
10	10	10	10
Supplier A	Supplier A	Supplier A	Supplier B

Effect of moisture on GMP powder

Preparation #	1	2	3	4	(-) L-Cystine
Powder storage day	30	36	43	67	31
Liquid media storage day	Day 0	Day 1	Day 12	Day 12	Day 12

Conclusions:

- Tighter control of age and exposure of powder
- Powder solubility diminish during storage after
 - 67 days in large lined bucket (5 kg)
 - 27 days in small container (1 kg)
- Hygroscopic components in powder- Choline Chloride
- Strong sulfide like odor indicates a degradation reaction
- Precipitate was observed after more than 30 days of storage

Powder instability was due to exposure: time + moisture during storage and handling

Powder degradation leads to precipitation.

- Changes implemented at FISI**
- Purge mill with dry Nitrogen (currently practiced in cGMP)
 - Control humidity in manufacturing areas (currently practiced in cGMP)
 - Add desiccants in packaging
 - Additional moisture barriers (bags)
 - Identify and then implement a superior custom container/closure system
- Recommendations for Company B**
- Allow packages to warm up to ambient temp prior to opening and use
 - Minimize time containers are open
 - Do not save/ reuse opened containers (use complete contents or discard excess powder if any remaining)
 - Frequent ordering of right-sized batches to minimize storage times

Case Study C

Goal: Investigate if factors in preparation method cause the aberrant cell culture performance

Altered metabolism and growth rate observed with liquid medium supplemented with L-Glutamine solution at point-of-use.

- Slower L-Glutamine utilization.
Lag in growth.
Premature decline in viability.
Increase osmolality due to:
- Elevated lactate level-> increased base demand
 - Reduced uptake of the main feed

Two factors identified for OFAT experiments


Factor Name	Tested Levels
Copper concentration	2x
Glutamine addition	Liquid Powder

Variable space for full factorial design of four factors.

The fractional factorial design has only 8 out of the possible sixteen conditions.

Identified factors for Design-of-Experiment (DOE)

Factor Name	Units	Low	High
Final pH		7.2	7.4
Final osmolality	mOsm/Kg	270	300
Final concentration	%	100	104
Glutamine concentration	m/M	4	8




Fractional factorial experiment designed to evaluate the effect of the four selected factors

Generated nine prototypes for customer to conduct 20 run screening study.

Run	Final pH	Final osmolality (mOsm/Kg)	Final concentration (%)	Glutamine concentration (mM)
1	7.2	270	100	4
2	7.2	270	104	8
3	7.2	300	104	8
4	7.2	300	100	4
5	7.4	270	100	4
6	7.4	270	104	8
7	7.4	300	104	4
8	7.4	300	100	8
9	7.4	270	100	8
10	7.4	270	104	4
11	7.4	300	104	8
12	7.4	300	100	4
13	7.4	270	100	8
14	7.4	270	104	4
15	7.4	300	104	8
16	7.4	300	100	4
17	7.4	270	100	8
18	7.4	270	104	4
19	7.4	300	104	8
20	7.4	300	100	4

Nine unique replicated conditions
Resolution IV, alpha = 0.05

Cell culture data generated using the Ambr[®]15 system and analyzed using Design-Expert[®] software.



Characterized media preparation process to fine tune method's set points and avoid edges of failure

Lowest growth observed in condition with all factors at the high levels of final pH, osmolality, concentration, and Glutamine concentration.

- Higher lactate seen in conditions with lower growth
- Low osmolality leads to lower lactate production.
- High osmolality overall negative effect
- High Final concentration appears to be beneficial to growth
- Interaction between pH and other factors appear to negatively affect growth when tested at their high levels
- High variation seen in midpoint's lactate production
- Copper effect seems to be cell line specific