For Research Use Only. Not for use in diagnostic procedure.

Code No. 291-96701

## LabAssay™ LDL-Cholesterol (LDL-Cholesterol Kit)

### 1. Intended use

This kit is a research reagent for measuring LDL-cholesterol.

One of the lipoproteins is LDL, with the cholesterol carried by LDL termed LDL-cholesterol. LDL-cholesterol carries cholesterol produced in the liver throughout the body, and when present in large amounts in the blood, it deposits and accumulates in blood vessel walls, causing atherosclerosis and resulting in myocardial infarction and cerebral infarction. In humans, the normal range of LDL-cholesterol is less than 140 mg/dL, and more than 140 mg/dL is considered to reflect hyper-LDL cholesterolemia.

### 2. Storage and expiration date

Store at 2°C to 10°C and do not freeze. The expiration date is indicated on the label on the outer box of the kit.

### 3. Kit Component Reagents



Components		Use Status	Amount
(1)	Pretreatment	Solution Ready for use.	20 mL / 1 vial
(2)	Reacting Solution	Solution Ready for use.	8 mL / 1 vial
(3)	LDL-Cholesterol Standard	Freeze-dried. Use after reconstitution.	2 vial
(4)	Standard Diluent	Solution Ready for use.	10 mL / 1 vial

\* For 100 times

### 4. Principle of the method

#### The first reaction (elimination of non-LDL-cholesterol)

The block polymer containing hydrophilic and hydrophobic moiety in Pretreatment selectively binds to LDL in the specimens, and protects it from enzyme (CHE, CO) reactions. CHE and CO react with non-LDL lipoprotein [chylomicron (CM), very low density lipoprotein (VLDL), high density lipoprotein (HDL)]. Hydrogen peroxide produced by the enzyme reactions with non-LDL cholesterol is decomposed to water by catalase in Pretreatment reagent.

#### The second reaction (color reaction of LDL-cholesterol)

When Reacting Solution is added, the cholesterol and its derivatives in LDL produce hydrogen peroxide by CHE and CO. Hydrogen peroxide produced by the enzyme reactions with LDL-cholesterol yields a blue color pigment upon oxidative condensation with N-(3-sulfopropyl)-3-methoxy-5-methylaniline (HMMPS) and 4-aminoantipyrine in the presence of peroxidase (POD). The amount of LDL-cholesterol contained in the sample is determined by measuring the absorbance of the blue color.

### 5. Equipment or supplies required but not provided in the kit

- 96-well microplate (transparent type)
- Micropipette
- Incubator maintained at 37°C
- Plate mixer
- Microplate reader with 600 nm wavelength filter

### 6. Preparation of reagents and standard solutions

- ① Pretreatment : Ready to use. After opening the bottle, store at 2-10°C and use within one month.
- ② Reacting Solution : Ready to use. After opening the bottle, store at 2-10°C and use within one month.
- ③ Dilution series of standard solution :

LDL-Cholesterol Standard with the correct volume\* of deionized water (or distilled water) described in a separate sheet to prepare the original standard solution (300 mg/dL). Then prepare a dilution series of the standard solutions as shown below using the Standard Diluent brought to room temperature.

\*Because the volume of deionized water (or distilled water) to be added to the freeze-dried standard differs depending on the lot, see the specified volume in the separate sheet provided. The standard solution containing purified water should be stored at 2-10°C and used within one month. Standard solutions prepared for each concentration should be used immediately and should not be stored.

#### Example of preparation of dilution series of standard solution

Conc. (mg/dL)	300	150	75.0	37.5	18.8	9.38	0
Standard solution (μL)	50	50*	50*	50*	50*	50*	0
Standard Diluent (μL)	0	50	50	50	50	50	50

\*One rank higher standard.

### 7. Preparation of specimen

#### Serum/Plasma

- Specimen analysis should be performed immediately after collection.
- Use fresh specimens. Do not use specimens after repeated freeze-thawing as this may denature lipoproteins.
- Anticoagulants such as heparin, citrate and EDTA do not significantly influence the assay when used in normal amounts.
- Hemolysis does not significantly influence the assay.
- Dilute specimen with saline and repeat the assay if the measured value exceeds the measurable range, and multiply the result by the dilution factor.
- When triglyceride in the specimen exceeds 1,000 mg/dL, dilute the specimen with saline and multiply the result by the dilution factor.

### 8. Assay procedure

#### Bring reagents to room temperature (20°C- 25°C) before use.

- (1) Dispense 180 μL of the Pretreatment into the 96-well microplate.
- (2) Dispense 5 μL of each concentration of the standard solution into the well in which the standard is to be measured.
- (3) Dispense 5 μL of sample into the well where the sample is to be measured.
- (4) Shake the 96-well microplate on a plate mixer (\*①).
- (5) React in a thermostatic bath at 37°C for 10 minutes.
- (6) Remove the 96-well microplate from the thermostatic bath and dispense 60 μL of reaction reagent into each well.
- (7) Shake the 96wells microplate on a plate mixer (\*①).
- (8) React in a thermostatic bath at 37°C for 10 minutes.
- (9) Measure the absorbance of each well at 600 nm (reference wavelength, 700 nm (\*①)) immediately using a plate reader.

\*① Guideline for shaking : 500 rpm-600 rpm for 10 seconds × 3 times.

### 9. Calculation

- (1) Create a calibration curve for each measurement, with the standard solution concentration (mg/dL) on the X-axis and absorbance on the Y-axis.
- (2) From the calibration curve, read the concentration (mg/dL) corresponding to the absorbance of the sample. If a diluted sample is used the concentration reading is multiplied by the sample dilution ratio to obtain the measurement value.


\*If the absorbance of the sample deviates from the calibration curve absorbance, prepare the sample to an appropriate

dilution with saline and repeat the assay.

\* For arithmetic operations in computer software, we recommend a cubic polynomial, 4- or 5-parameter specification.

## 10. Standard Curve

• Measuring range : 9.38-300 mg/dL

## 11. Performance

### • Sensitivity

- When a sample subject to a standard dilution is measured, the absorbance is less than 0.030.
- When a standard solution of a specific concentration (LDL-cholesterol 300 mg/dL) is measured as a sample, the absorbance ranges from 0.190 to 0.350.

### • Accuracy

- When measuring control serum of known concentration, performance is within  $\pm 15\%$  of the known concentration.

## Notes

- Store the reagents under the specified conditions. Do not use reagents that have passed the expiration date stated on each reagent container label.
- Do not use reagents that were frozen in error. Such reagents may give false results.
- Operate out of direct sunlight.
- Avoid contamination from the micropipette when collecting reagents.
- For Research Use Only. Not for use in diagnostic procedure.
- If any reagent comes into contact with the mouth, eyes, or skin, wash off immediately with a large amount of water. Consult a physician if necessary.
- When discarding the reagents, dispose of them according to local or national regulations.
- All devices including reagents and reagent bottles that come into contact with specimens should be considered potentially infectious.
- The Reacting Solution contains 0.09% sodium azide as a preservative. Sodium azide may react with copper or lead plumbing to form explosive compounds. Even though the reagent contains minute quantity of sodium azide, drains should be flushed copiously with a large amount of water, when discarding these reagents.

## 【References】

- 1) NIH Publication No. 95-3044, Recommendations on Lipoprotein Measurement (1995)
- 2) Japan Atherosclerosis Society : Guidelines for Prevention of Atherosclerotic Cardiovascular Diseases, 5-7 (2002)

LabAssay™ LDL-Cholesterol

[Storage condition] Store the kit at 2°C - 10°C (Do not freeze)

[Term of validity] Expiration date is indicated on the container

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