### CHANG Medium C For Human Amniotic Fluid Cells

#### Catalog No.

**T101-019 kit includes (C100 & C107)**

**T101-059 kit includes (C101 & C109),**

**C100 (CHANG Medium Basal)**

**+C106 (CHANG Medium C Frozen Supplement-14mL),**

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Za diagnostično uporabo *in vitro*.

#### Glossary of Symbols\*:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature 2-8°C for T101-019, T101-059, C100, and C101
	Storage Temperature -10°C for C106 and C108
	Do Not Re-Sterilize
	Do not use if package is damaged
	Manufacturer
	CE Mark
	Emergo Europe - Prinsessegracht 20 2514 AP The Hague The Netherlands

\*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

#### INDICATION FOR USE

CHANG Medium C may be used for the following applications:

1. the primary culture of amniotic fluid cells
  2. growing passaged amniotic fluid cells
  3. the culture of bone marrow cells
  4. solid amnion tissue from chorionic villi sampling.
- THIS MEDIUM HAS BEEN DESIGNED FOR USE IN BOTH OPEN (CULTURES EQUILIBRATED WITH A 5% CO<sub>2</sub> ATMOSPHERE) AND CLOSED CULTURE SYSTEMS (CULTURES FLOODED WITH 5% CO<sub>2</sub> AND INCUBATED WITH THE CULTURE VESSELS TIGHTLY CAPPED).

#### DEVICE DESCRIPTION

CHANG Medium C was developed for the primary culture of human amniotic fluid cells for use in karyotyping and other antenatal genetic testing. THIS FORMULA HAS BEEN OPTIMIZED FOR BOTH OPEN AND CLOSED SYSTEMS.

#### COMPONENTS

Salts & Ions	Proteins	pH Indicator
Sodium Chloride	Hormones, and	Phenol red
Sodium Selenite	Growth Factors	Vitamins and trace elements
Calcium Chloride	Fetal bovine serum (FBS)	Ascorbic acid
Choline Chloride	Newborn bovine serum	Folic acid
Potassium Chloride	Human transferrin	Nicotinamide
Potassium Phosphate	Fibroblast growth factor (FGF)	Thiamine
Magnesium Chloride	insulin	Riboflavin
Magnesium Sulfate	Progesterone	Pantothenic acid
Sodium Phosphate	Testosterone	Cobalamin
Ferrous Sulfate	Beta estradiol	Pyridoxal
Cupric Sulfate	Hydrocortisone	Pyridoxine
Zinc Sulfate	Water	
Buffer	WFI Quality	
Sodium Bicarbonate	Other	
Amino Acids	Hypoxanthine	
Alanine	Biotin	
Arginine	Putrescine	
Asparagine	Ethyl alcohol	
Aspartic Acid	Antioxidant	
Cysteine	Thioctic acid	
Cystine	Nucleic acids	
Glutamic Acid	Deoxyadenosine	
Glutamine	Deoxycytidine	
Glycine	Deoxyguanosine	
Histidine	Adenosine	
Isoleucine	Cytidine	
Leucine	Guanosine	
Lysine	Thymidine	
Methionine	Uridine	
Phenylalanine	Energy Substrates	
Proline	Glucose	
Serine	Pyruvate	
Threonine	Inositol	
Tryptophan		
Tyrosine		
Valine		

#### QUALITY ASSURANCE

##### STERILITY

Serum used in the production of CHANG Medium C Supplement has been tested for viral contamination per CFR Title 9 Part 113.53. It has also been screened for mycoplasma contamination. Both CHANG Medium B Basal and CHANG Medium C Supplement are sterilized by filtration through a 0.1 micron filter. Representative samples of CHANG Medium B Basal and CHANG Medium C Supplement are tested for possible bacteriological contamination following the sterility testing protocol described in the current USP Sterility test <71>.

##### PREPARATION FOR USE

###### Lyophilized Supplement

1. Allow CHANG Medium C Lyophilized Supplement to equilibrate to room temperature.
2. Using aseptic techniques, add 10 mL or 50 mL of sterile distilled water to the lyophilizate. The 10 mL vial is for 90 mL CHANG Medium B Basal; the 50 mL vial is for 450 mL CHANG Medium B Basal.
3. Swirl the vial gently to effect complete solution. (NOTE: It may be necessary to incubate the vial at 37°C for 5 to 10 minutes to completely dissolve the supplement).
4. Aseptically transfer the entire contents of the reconstituted CHANG Medium C Supplement to the

bottle of CHANG Medium B Basal.

5. Mix the complete CHANG Medium C well by swirling the bottle.
6. Add L-Glutamine, 10.0 mL/L (200 mM). Antibiotics may be added if desired.

###### Frozen Supplement

1. Thaw CHANG Medium C Frozen Supplement rapidly by swirling vial in a 37°C water bath.
2. Aseptically transfer the entire contents of the thawed CHANG Medium C Supplement to the bottle of CHANG Medium B Basal.
3. Mix the complete CHANG Medium C well by swirling the bottle.
4. Add L-Glutamine, 10.0 mL/L (200 mM). Antibiotics may be added if desired.

###### ALIQUOTING CHANG MEDIUM C SUPPLEMENT

###### Lyophilized Supplement

If you are unable to use the complete CHANG Medium C within 10 days you may wish to make up smaller amounts by following this modification of procedure:

1. Reconstitute CHANG Medium C Lyophilized Supplement.
2. Distribute aseptically into convenient-sized aliquots and freeze.
3. Also aseptically dispense the CHANG Medium B Basal into a corresponding number of aliquots.
4. To prepare the complete CHANG Medium C, add one aliquot of thawed CHANG Medium C Supplement to one aliquot of CHANG Medium B Basal.
5. Mix well.

###### Frozen Supplement

1. Thaw CHANG Medium C Frozen Supplement.
2. Distribute aseptically into convenient-sized aliquots and refreeze. (NOTE: The supplement vial contains 14 mL or 70 mL. Distribute equally.)
3. Also aseptically dispense the CHANG Medium B Basal into a corresponding number of aliquots.
4. To prepare the complete CHANG Medium C, add one aliquot of thawed CHANG Medium C Supplement to one aliquot of CHANG Medium B Basal.
5. Mix well.

###### DIRECTIONS FOR USE

THE pH OF THE MEDIUM USED TO FEED THE CULTURES MUST BE BETWEEN 6.8 - 7.2 (i.e. the medium must be slightly yellowish-salmon color). pH can easily be adjusted by placing the medium in a 5% - 8% CO<sub>2</sub> incubator with the cap slightly loosened. The final pH must be 6.8 - 7.2.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

###### Use of CHANG Medium C for Primary Cultures: In Situ Methodologies

1. Centrifuge amniotic fluid at low speed to concentrate the cells.
2. Resuspend the cell pellet in a small volume of the patient's own amniotic fluid. For example, aspirate the supernate of 10 mL of spun amniotic fluid to 0.5 mL above the cell pellet and resuspend. Add sufficient CHANG Medium C to the concentrated cell suspension to allow for final plating volume of 0.5 mL per cover slip (total of 4 coverslips) or 2 mL per flaskette.
3. Incubate cultures undisturbed at 37°C 5% CO<sub>2</sub> atmosphere.
4. Flood cultures on day 2 by adding 2 mL of CHANG Medium C.
5. After 4 to 5 days, the cultures should be checked for growth. Cultures should be fed once growth has been observed. Feed cultures by removing all of the culture supernatant and replacing with 2 mL of fresh CHANG Medium C. It is recommended that cultures be fed every 2 days thereafter.
6. Check cultures for growth on/or after day 5, and harvest when sufficient colonies are observed.
7. Best results are obtained when the cultures are fed with CHANG Medium C the day before the harvest.

###### Use of CHANG Medium C for Primary Cultures: Flask Methodologies

1. Centrifuge amniotic fluid at low speed to concentrate cells.
2. Resuspend the cell pellet in a small volume of the patient's own amniotic fluid. For example, aspirate the supernate of 10 mL of spun amniotic fluid to 1 mL above the cell pellet and resuspend. Add 4 mL of CHANG Medium C for a total volume of 5 mL per flask.
3. Incubate cultures undisturbed at 37°C 5% CO<sub>2</sub> atmosphere.
4. Check for growth on day 5. Change medium with fresh CHANG Medium C and harvest if sufficient cell growth is observed.
5. Check cultures for growth and completely change medium every other day thereafter until sufficient colonies are observed and are ready to harvest.
6. Best results are obtained when the cultures are fed with CHANG Medium C the day before the harvest.

NOTE: For closed systems, flush each culture flask with 5% CO<sub>2</sub> - 95% air for 20 seconds. Tighten the caps to the culture flasks and incubate them at 37°C. (It is recommended that a sterile plugged pasteur pipet be attached to the CO<sub>2</sub> source to ensure sterility of the incoming gas.)

###### Use of CHANG Medium C for Growing Passaged Amniotic Fluid Cells:

To passage the cells, treat the cultures with trypsin (or pronase, etc.) as you would normally do when cells are grown in conventional medium. However, protease treatment should be carefully monitored. Amniotic fluid cells grown in CHANG Medium C tend to be more sensitive to protease treatment than amniotic fluid cells grown in conventional medium. It may be necessary to modify your protocol to take this into account.

###### STORAGE AND STABILITY

Store CHANG Medium C Frozen Supplement at below -10°C, CHANG Medium C Lyophilized Supplement at 2°C to 8°C, CHANG Medium B Basal at 2°C to 8°C, and complete CHANG Medium C at 2°C to 8°C. The CHANG Medium B Basal must not be frozen.

Protect from fluorescent light.

See individual component bottle labels for specific expiration dates. Complete CHANG Medium C may be stored at 2°C to 8°C for 10 days before use without affecting its function. Storage for longer than 10 days is not recommended.

DO NOT FREEZE COMPLETE CHANG MEDIUM C.

###### PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in procedures that include the indicated application for which the device is intended.

Do not use CHANG Medium C components beyond the expiration date indicated on the individual labels.

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