CHANG Medium C For Human Amniotic Fluid Cells

Catalog No.

T101-019 kit includes (C100 & C107)

T101-059 kit includes (C101 & C109),

C100 (CHANG Medium Basal)

+C106 (CHANG Medium C Frozen Supplement-14mL),

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Sterilized using aseptic processing techniques (filtration)

Expiration:
Year - Month - Day

Caution, consult accompanying documents

Consult instructions for use

Storage Temperature 2-8°C for
T101-019, T101-059,
C100, and C101

Storage Temperature -10°C for
C106 and C108

Do Not Re-Sterilize

Do not use if package is damaged

Manufacturer

CE Mark

Emergo Europe - Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices - Symbols to be used with medical device labels, labeling.

ENGLISH

INDICATION FOR USE

CHANG Medium C may be used for the following applications:

1. the primary culture of amniotic fluid cells
2. growing passaged amniotic fluid cells
3. the culture of bone marrow cells
4. solid amnion tissue from chorionic villi sampling.

THIS MEDIUM HAS BEEN DESIGNED FOR USE IN BOTH OPEN (CULTURES EQUILIBRATED WITH A 5% CO₂ ATMOSPHERE) AND CLOSED CULTURE SYSTEMS (CULTURES FLOODED WITH 5% CO₂ AND INCUBATED WITH THE CULTURE VESSELS TIGHTLY CAPPED).

DEVICE DESCRIPTION

CHANG Medium C was developed for the primary culture of human amniotic fluid cells for use in karyotyping and other antenatal genetic testing. THIS FORMULA HAS BEEN OPTIMIZED FOR BOTH OPEN AND CLOSED SYSTEMS.

COMPONENTS

<u>Salts & Ions</u>	<u>Proteins,</u>	<u>pH Indicator</u>
Sodium Chloride	<u>Hormones, and</u>	Phenol red
Sodium Selenite	<u>Growth Factors</u>	<u>Vitamins and trace</u>
Calcium Chloride	Fetal bovine serum	<u>elements</u>
Choline Chloride	(FBS)	Ascorbic acid
Potassium Chloride	Newborn bovine	Folic acid
Potassium	serum	Riboflavin
Phosphate	Human transferrin	Nicotinamide
Magnesium Chloride	Fibroblast growth	Thiamine
Magnesium Sulfate	factor (FGF)	Pantothenic acid
Sodium Phosphate	Insulin	Cobalamin
Ferrous Sulfate	Progesterone	Pyridoxal
Cupric Sulfate	Testosterone	Pyridoxine
Zinc Sulfate	Beta estradiol	
<u>Buffer</u>	Hydrocortisone	
Sodium Bicarbonate	<u>Water</u>	
<u>Amino Acids</u>	WFI Quality	
Alanine	Other	
Arginine	Hypoxanthine	
Asparagine	Biotin	
Aspartic Acid	Putrescine	
Cysteine	Ethyl alcohol	
Cystine	<u>Antioxidant</u>	
Glutamic Acid	Thioctic acid	
Glutamine	<u>Nucleic acids</u>	
Glycine	Deoxyadenosine	
Histidine	Deoxycytidine	
Isoleucine	Deoxyguanosine	
Leucine	Adenosine	
Lysine	Cytidine	
Methionine	Guanosine	
Phenylalanine	Thymidine	
Proline	Uridine	
Serine	<u>Energy Substrates</u>	
Threonine	Glucose	
Tryptophan	Pyruvate	
Tyrosine	Inositol	
Valine		

QUALITY ASSURANCE

STERILITY

Serum used in the production of CHANG Medium C Supplement has been tested for viral contamination per CFR Title 9 Part 113.53. It has also been screened for mycoplasma contamination. Both CHANG Medium B Basal and CHANG Medium C Supplement are sterilized by filtration through a 0.1 micron filter. Representative samples of CHANG Medium B Basal and CHANG Medium C Supplement are tested for possible bacteriological contamination following the sterility testing protocol described in the current USP Sterility test <71>.

PREPARATION FOR USE

Lyophilized Supplement

1. Allow CHANG Medium C Lyophilized Supplement to equilibrate to room temperature.
2. Using aseptic techniques, add 10 mL or 50 mL of sterile distilled water to the lyophilizate. The 10 mL vial is for 90 mL CHANG Medium B Basal; the 50 mL vial is for 450 mL CHANG Medium B Basal.
3. Swirl the vial gently to effect complete solution. (NOTE: It may be necessary to incubate the vial at 37°C for 5 to 10 minutes to completely dissolve the supplement).
4. Aseptically transfer the entire contents of the reconstituted CHANG Medium C Supplement to the

bottle of CHANG Medium B Basal.

5. Mix the complete CHANG Medium C well by swirling the bottle.
6. Add L-Glutamine, 10.0 mL/L (200 mM). Antibiotics may be added if desired.

Frozen Supplement

1. Thaw CHANG Medium C Frozen Supplement rapidly by swirling vial in a 37°C water bath.
2. Aseptically transfer the entire contents of the thawed CHANG Medium C Supplement to the bottle of CHANG Medium B Basal.
3. Mix the complete CHANG Medium C well by swirling the bottle.
4. Add L-Glutamine, 10.0 mL/L (200 mM). Antibiotics may be added if desired.

ALIQUOTING CHANG MEDIUM C SUPPLEMENT

Lyophilized Supplement

If you are unable to use the complete CHANG Medium C within 10 days you may wish to make up smaller amounts by following this modification of procedure:

1. Reconstitute CHANG Medium C Lyophilized Supplement.
2. Distribute aseptically into convenient-sized aliquots and freeze.
3. Also aseptically dispense the CHANG Medium B Basal into a corresponding number of aliquots.
4. To prepare the complete CHANG Medium C, add one aliquot of thawed CHANG Medium C Supplement to one aliquot of CHANG Medium B Basal.
5. Mix well.

Frozen Supplement

1. Thaw CHANG Medium C Frozen Supplement.
2. Distribute aseptically into convenient-sized aliquots and refreeze. (NOTE: The supplement vial contains 14 mL or 70 mL. Distribute equally.)
3. Also aseptically dispense the CHANG Medium B Basal into a corresponding number of aliquots.
4. To prepare the complete CHANG Medium C, add one aliquot of thawed CHANG Medium C Supplement to one aliquot of CHANG Medium B Basal.
5. Mix well.

DIRECTIONS FOR USE

THE pH OF THE MEDIUM USED TO FEED THE CULTURES MUST BE BETWEEN 6.8 - 7.2 (i.e. the medium must be slightly yellowish-salmon color). pH can easily be adjusted by placing the medium in a 5% - 8% CO₂ incubator with the cap slightly loosened. The final pH must be 6.8 - 7.2.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

Use of CHANG Medium C for Primary Cultures: In Situ Methodologies

1. Centrifuge amniotic fluid at low speed to concentrate the cells.
2. Resuspend the cell pellet in a small volume of the patient's own amniotic fluid. For example, aspirate the supernate of 10 mL of spun amniotic fluid to 0.5 mL above the cell pellet and resuspend. Add sufficient CHANG Medium C to the concentrated cell suspension to allow for final plating volume of 0.5 mL per cover slip (total of 4 coverslips) or 2 mL per flaskette.
3. Incubate cultures undisturbed at 37°C 5% CO₂ atmosphere.
4. Flood cultures on day 2 by adding 2 mL of CHANG Medium C.
5. After 4 to 5 days, the cultures should be checked for growth. Cultures should be fed once growth has been observed. Feed cultures by removing all of the culture supernatant and replacing with 2 mL of fresh CHANG Medium C. It is recommended that cultures be fed every 2 days thereafter.
6. Check cultures for growth on/or after day 5, and harvest when sufficient colonies are observed.
7. Best results are obtained when the cultures are fed with CHANG Medium C the day before the harvest.

Use of CHANG Medium C for Primary Cultures: Flask Methodologies

1. Centrifuge amniotic fluid at low speed to concentrate cells.
2. Resuspend the cell pellet in a small volume of the patient's own amniotic fluid. For example, aspirate the supernate of 10 mL of spun amniotic fluid to 1 mL above the cell pellet and resuspend. Add 4 mL of CHANG Medium C for a total volume of 5 mL per flask.
3. Incubate cultures undisturbed at 37°C 5% CO₂ atmosphere.
4. Check for growth on day 5. Change medium with fresh CHANG Medium C and harvest if sufficient cell growth is observed.
5. Check cultures for growth and completely change medium every other day thereafter until sufficient colonies are observed and are ready to harvest.
6. Best results are obtained when the cultures are fed with CHANG Medium C the day before the harvest.

NOTE: For closed systems, flush each culture flask with 5% CO₂ - 95% air for 20 seconds. Tighten the caps to the culture flasks and incubate them at 37°C. (It is recommended that a sterile plugged pasteur pipet be attached to the CO₂ source to ensure sterility of the incoming gas.)

Use of CHANG Medium C for Growing Passaged Amniotic Fluid Cells:

To passage the cells, treat the cultures with trypsin (or pronase, etc.) as you would normally do when cells are grown in conventional medium. However, protease treatment should be carefully monitored. Amniotic fluid cells grown in CHANG Medium C tend to be more sensitive to protease treatment than amniotic fluid cells grown in conventional medium. It may be necessary to modify your protocol to take this into account.

STORAGE AND STABILITY

Store CHANG Medium C Frozen Supplement at below -10°C, CHANG Medium C Lyophilized Supplement at 2°C to 8°C, CHANG Medium B Basal at 2°C to 8°C, and complete CHANG Medium C at 2°C to 8°C. The CHANG Medium B Basal must not be frozen. Protect from fluorescent light.

See individual component bottle labels for specific expiration dates. Complete CHANG Medium C may be stored at 2°C to 8°C for 10 days before use without affecting its function. Storage for longer than 10 days is not recommended.

DO NOT FREEZE COMPLETE CHANG MEDIUM C.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in procedures that include the indicated application for which the device is intended.

Do not use CHANG Medium C components beyond the expiration date indicated on the individual labels.

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KOMPONENTI	Imluha u Joni	Proteini, Ormoni u Fakturi ta' Tkabbir	Indikator tal-pH
	Sodium Chloride	Fetal bovine serum (FBS)	Phenol Red
	Sodium Selenite	Serum ta' bovin tal-twelid	Vitamini u mikroelementi
	Calcium Chloride	Potassium Chloride	Ascorbic acid
	Choline Chloride	Potassium Chloride	Folic acid
	Potassium Chloride	Transferin uman	Nicotinamide
	Phosphate	Fattur ta' tkabir tal-fibroblast (FGF)	Riboflavin
	Magnesium Chloride	Insulina	Thiamine
	Magnesium Sulfate	Progesterone	Pantothenic acid
	Sodium Phosphate	Testosterone	Cobalamin
	Ferrous Sulfate	Beta estradiol	Pyridoxin
	Cupric Sulfate	Hydrocortisone	Pyridoxine
	Zinc Sulfate	Bafer	
	Bafer	Sodium Bicarbonate	
	Alanine	Amino Acids	
	Arginine	Putrescine	
	Asparagine	Cysteine	
	Aspartic Acid	Cystine	
	Cysteine	Glutamic Acid	
	Glutamic Acid	Glutamine	
	Glutamine	Glycine	
	Glycine	Histidine	
	Histidine	Isoleucine	
	Isoleucine	Leucine	
	Leucine	Lysine	
	Lysine	Methionine	
	Methionine	Phenylalanine	
	Phenylalanine	Proline	
	Proline	Serine	
	Serine	Threonine	
	Threonine	Tryptophan	
	Tryptophan	Tyrosine	
	Tyrosine	Valine	
	Valine	Inostol	

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