CHANG Medium *In Situ*

For Human Amniotic Fluid Cells

Catalog No. T104

For *in vitro* diagnostic use.

Zur *In-vitro-Diagnostik*.

Solo per uso diagnostico *in vitro*.

Para uso diagnóstico *in vitro*.

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Sterilized using aseptic processing techniques (filtration)

Expiration: Year - Month - Day

Caution, consult accompanying documents

Consult instructions for use

Storage Temperature below -10°C

Do not resterilize

Do not use if package is damaged

Manufacturer

CE Mark

Emergo Europe - Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

ENGLISH

INDICATION FOR USE

CHANG Medium *In Situ* may be used for the following applications:

1. the primary culture of amniotic fluid cells
2. growing passaged amniotic fluid cells
3. solid amnion tissue from chorionic villi sampling.

This medium has been designed for use in CO_2 incubators (cultures equilibrated with 5%-8% CO_2 atmosphere).

DEVICE DESCRIPTION

CHANG Medium *In Situ* was developed for the primary culture of human amniotic fluid cells for use in karyotyping and other antenatal genetic testing. This formula has been optimized for *In Situ* methodologies.

COMPONENTS

Amino Acids	Potassium Phosphate
Alanine	Calcium Chloride
Arginine	Magnesium Sulfate
Asparagine	Choline Chloride
Aspartic Acid	Magnesium Chloride
Cysteine	Sodium Selenite
Cystine	Cupric Sulfate
Glutamic Acid	Ferrous Sulfate
Glutamine	Zinc Sulfate
Glycine	Proteins, Hormones, and Growth Factors
Histidine	Fetal Bovine Serum
Isoleucine	Newborn Bovine Serum
Leucine	Fibroblast Growth Factor (FGF)
Lysine	Transferrin
Methionine	Insulin
Phenylalanine	Progesterone
Proline	Testosterone
Serine	B-Estradiol
Threonine	Hydrocortisone
Tryptophan	Nucleic acids
Tyrosine	Cytidine
Valine	Sodium Bicarbonate
Buffer	Deoxyadenosine
	Vitamins and trace elements
	Deoxyguanosine
Biotin	Biotin
Riboflavin	Deoxyuridine
Ascorbic Acid	Guanosine
Folic Acid	Thymidine
Nicotinic Acid	Uridine
Pantothenic Acid	Adenosine
Pyridoxal	Hypoxanthine
Pyridoxine	Antioxidant
Thiamine	Thioclic Acid
Vitamin B12	Others
Energy Substrates	Ethyl Alcohol
Glucose	Putrescine
Inositol	
Pyruvate	
pH Indicator	
Phenol Red	
Salts & Ions	
Sodium Chloride	
Potassium Chloride	
Sodium Phosphate	

QUALITY ASSURANCE

STERILITY

Serum used in the production of CHANG Medium *In Situ* has been tested for viral contamination per CFR Title 9 Part 113.53. It has also been screened for mycoplasma contamination. CHANG Medium *In Situ* is sterilized by filtration through a 0.1 micron filter. Samples of CHANG Medium *In Situ* are tested for possible bacteriological contamination following the sterility testing protocol described in the current USP Sterility test <71>.

PREPARATION FOR USE

1. Thaw CHANG Medium *In Situ* rapidly by swirling bottle in a 37°C water bath.
2. Antibiotics may be added if desired.

ALIQUOTING CHANG MEDIUM *In Situ*

1. Thaw CHANG Medium *In Situ* according to instructions.
2. Distribute aseptically into convenient sized aliquots and refreeze.
3. Thaw aliquots in 37°C water bath when ready to use.

DIRECTIONS FOR USE

The pH of the medium used to feed the cultures must be between 6.8-7.2 (i.e. the medium must be slightly yellowish salmon color). pH can easily be adjusted by placing the medium in a 5%-8% CO_2 incubator with the cap slightly loosened for about 30 minutes.

The final pH must be 6.8-7.2.

Use of CHANG Medium *In Situ* for Primary Cultures: in situ Methodologies

1. Centrifuge amniotic fluid at low speed to concentrate the cells.
2. Resuspend the cell pellet in a small volume of the patient's own amniotic fluid. For example, aspirate the supernate of 10 mL of spun amniotic fluid to 0.5 mL above the cell pellet and resuspend. Add sufficient CHANG Medium *In Situ* to the concentrated cell suspension to allow for final plating volume of 0.5 mL per cover slip (total of 4 coverslips) or 2 mL per flaskette.
3. Incubate cultures undisturbed at 37°C 5%-8% CO_2 atmosphere.
4. Flood cultures on day 2 by adding 2 mL of CHANG Medium *In Situ*.
5. After 4 to 5 days, the cultures should be checked for growth. Cultures should be fed once growth has been observed. Feed cultures by removing all of the culture supernatant and replacing with 2 mL fresh CHANG Medium *In Situ*. It is recommended that cultures be fed every 2 days thereafter.
6. Check cultures for growth on/or after day 5, and harvest when sufficient colonies are observed.
7. Best results are obtained when the cultures are fed with CHANG Medium *In Situ* the day before the harvest.

Use of CHANG Medium *In Situ* for Growing Passaged Amniotic Fluid Cells:

To passage the cells, treat the cultures with trypsin (or pronase, etc.) as you would normally do when cells are grown in conventional medium. However, protease treatment should be carefully monitored. Amniotic fluid cells grown in CHANG Medium *In Situ* tend to be more sensitive to protease treatment than amniotic fluid cells grown in conventional medium. It may be necessary to modify your protocol to take this into account.

Note: Calcium Oxalate crystals commonly form in CHANG Medium *In Situ*. The presence of these crystals has not been shown to cause any detrimental effect on product performance.

STORAGE AND STABILITY

Store CHANG Medium *In Situ* frozen at -10°C. Unused CHANG Medium *In Situ* can be refrozen or stored at 2°C to 8°C.

Protect from fluorescent light.

See bottle label for specific expiration date. CHANG Medium *In Situ* may be refrozen a maximum of 2 times and stored thawed at 2°C to 8°C for 14 days without affecting its function. Storage for longer than 14 days is not recommended.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in procedures that include the indicated application for which the device is intended.

Do not use any bottle in which the sterile packaging has been compromised

Do not use CHANG Medium *In Situ* beyond the expiration date indicated on the label.

FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706 • Fax: 1 949 261 6522 • www.irvinesci.com

© 2023 FUJIFILM Irvine Scientific, Inc. All rights reserved. The FUJIFILM Irvine Scientific logo and CHANG Medium are trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions.

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