FUJIFILM

Irvine Scientific

CHANG Medium *In Situ* For Human Amniotic Fluid Cells

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Za diagnostično uporabo *in vitro*.

Glossary of Symbols*:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature below -10°C
	Do not resterilize
	Do not use if package is damaged
	Manufacturer
	CE Mark
	Emergo Europe - Prinsessegracht 20 2514 AP The Hague The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

ENGLISH

INDICATION FOR USE

CHANG Medium *In Situ* may be used for the following applications:

1. the primary culture of amniotic fluid cells
2. growing passaged amniotic fluid cells
3. solid amnion tissue from chorionic villi sampling.

This medium has been designed for use in CO₂ incubators (cultures equilibrated with 5%-8% CO₂ atmosphere).

DEVICE DESCRIPTION

CHANG Medium *In Situ* was developed for the primary culture of human amniotic fluid cells for use in karyotyping and other antenatal genetic testing. This formula has been optimized for *In Situ* methodologies.

COMPONENTS

<u>Amino Acids</u>	Potassium Phosphate
Alanine	Calcium Chloride
Arginine	Magnesium Sulfate
Asparagine	Choline Chloride
Aspartic Acid	Magnesium Chloride
Cysteine	Sodium Selenite
Cystine	Cupric Sulfate
Glutamic Acid	Ferrous Sulfate
Glutamine	Zinc Sulfate
Glycine	<u>Proteins, Hormones, and Growth Factors</u>
Histidine	Fetal Bovine Serum
Isoleucine	Newborn Bovine Serum
Leucine	Fibroblast Growth Factor
Lysine	(FGF)
Methionine	Transferrin
Phenylalanine	Insulin
Proline	Progesterone
Serine	Testosterone
Threonine	B-Estradiol
Tryptophan	Hydrocortisone
Tyrosine	<u>Nucleic acids</u>
Valine	Cytidine
<u>Buffer</u>	Deoxyadenosine
Sodium Bicarbonate	Deoxycytidine
<u>Vitamins and trace elements</u>	Deoxyguanosine
Biotin	Guanosine
Riboflavin	Thymidine
Ascorbic Acid	Uridine
Folic Acid	Adenosine
Nicotinic Acid	Hypoxanthine
Pantothenic Acid	<u>Antioxidant</u>
Pyridoxal	Thioctic Acid
Pyridoxine	<u>Others</u>
Thiamine	Ethyl Alcohol
Vitamin B12	Putrescine
<u>Energy Substrates</u>	
Glucose	
Inositol	
Pyruvate	
<u>pH indicator</u>	
Phenol Red	
<u>Salts & Ions</u>	
Sodium Chloride	
Potassium Chloride	
Sodium Phosphate	

QUALITY ASSURANCE

STERILITY

Serum used in the production of CHANG Medium *In Situ* has been tested for viral contamination per CFR Title 9 Part 113.53. It has also been screened for mycoplasma contamination. CHANG Medium *In Situ* is sterilized by filtration through a 0.1 micron filter. Samples of CHANG Medium *In Situ* are tested for possible bacteriological contamination following the sterility testing protocol described in the current USP Sterility test <71>.

PREPARATION FOR USE

1. Thaw CHANG Medium *In Situ* rapidly by swirling bottle in a 37°C water bath.
2. Antibiotics may be added if desired.

ALIQUOTING CHANG MEDIUM *In Situ*

1. Thaw CHANG Medium *In Situ* according to instructions.
2. Distribute aseptically into convenient sized aliquots and refreeze.
3. Thaw aliquots in 37°C water bath when ready to use.

DIRECTIONS FOR USE

The pH of the medium used to feed the cultures must be between 6.8-7.2 (i.e. the medium must be slightly yellowish salmon color). pH can easily be adjusted by placing the medium in a 5%-8% CO₂ incubator with the cap slightly loosened for about 30 minutes.

The final pH must be 6.8-7.2.

Use of CHANG Medium *In Situ* for Primary Cultures: *in situ* Methodologies

1. Centrifuge amniotic fluid at low speed to concentrate the cells.
2. Resuspend the cell pellet in a small volume of the patient's own amniotic fluid. For example, aspirate the supernate of 10 mL of spun amniotic fluid to 0.5 mL above the cell pellet and resuspend. Add sufficient CHANG Medium *In Situ* to the concentrated cell suspension to allow for final plating volume of 0.5 mL per cover slip (total of 4 coverslips) or 2 mL per flaskette.
3. Incubate cultures undisturbed at 37°C 5%-8% CO₂ atmosphere.
4. Flood cultures on day 2 by adding 2 mL of CHANG Medium *In Situ*.
5. After 4 to 5 days, the cultures should be checked for growth. Cultures should be fed once growth has been observed. Feed cultures by removing all of the culture supernatant and replacing with 2 mL fresh CHANG Medium *In Situ*. It is recommended that cultures be fed every 2 days thereafter.
6. Check cultures for growth on or after day 5, and harvest when sufficient colonies are observed.
7. Best results are obtained when the cultures are fed with CHANG Medium *In Situ* the day before the harvest.

Use of CHANG Medium *In Situ* for Growing Passaged Amniotic Fluid Cells:

To passage the cells, treat the cultures with trypsin (or pronase, etc.) as you would normally do when cells are grown in conventional medium. However, protease treatment should be carefully monitored. Amniotic fluid cells grown in CHANG Medium *In Situ* tend to be more sensitive to protease treatment than amniotic fluid cells grown in conventional medium. It may be necessary to modify your protocol to take this into account.

Note: Calcium Oxalate crystals commonly form in CHANG Medium *In Situ*. The presence of these crystals has not been shown to cause any detrimental effect on product performance.

STORAGE AND STABILITY

Store CHANG Medium *In Situ* frozen at -10°C. Unused CHANG Medium *In Situ* can be refrozen or stored at 2°C to 8°C.

Protect from fluorescent light.

See bottle label for specific expiration date. CHANG Medium *In Situ* may be refrozen a maximum of 2 times and stored thawed at 2°C to 8°C for 14 days without affecting its function. Storage for longer than 14 days is not recommended.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in procedures that include the indicated application for which the device is intended.

Do not use any bottle in which the sterile packaging has been compromised

Do not use CHANG Medium *In Situ* beyond the expiration date indicated on the label.

FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706 • Fax: 1 949 261 6522 • www.irvinesci.com

© 2019 FUJIFILM Irvine Scientific, Inc. All rights reserved. The FUJIFILM Irvine Scientific logo and CHANG Medium are trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions. PN 40280 Rev.8

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