FUJIFILM

IrvineScientific

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STERILE A Sterilized using aseptic processing techniques (filtration)

Expiration:
Year - Month - Day

Caution, consult accompanying documents

Consult instructions for use

Storage Temperature
2-8°C

Do Not Re-Sterilize

Do not use if package is damaged

Manufacturer

Rx Only U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

CE 0050 CE Mark

EC REP Emergo Europe - Prinsessegracht 20
2514 AP The Hague
The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

ENGLISH

EU CAUTION: For Professional Use Only.

INDICATION FOR USE

Modified HTF (mHTF) is intended for use in assisted reproductive procedures which include retrieval, handling, and transfer of human gametes and embryos.

DEVICE DESCRIPTION

mHTF is a synthetic HEPES-based solution for procedures including the retrieval, handling and transfer of human gametes and embryos. This buffering system does not require the use of a CO₂ incubator. This medium contains the antibiotic Gentamicin Sulfate (10 µg/mL) and requires protein supplement.

COMPOSITIONS:

Salts and Ions	Energy Sources
Sodium Chloride	Glucose
Potassium Chloride	Sodium Pyruvate
Magnesium Sulfate	Sodium Lactate
Potassium Phosphate	Antibiotic
Calcium Chloride	Gentamicin Sulfate
Buffers	pH Indicator
HEPES	Phenol Red
Sodium Bicarbonate	Water
	WFI Quality

QUALITY ASSURANCE

mHTF is a culture medium which is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10⁻³.

Each lot of mHTF is tested for:

- Endotoxin by Limulus Amebocyte Lysate (LAL) methodology
- Biocompatibility by Mouse Embryo Assay (one-cell)
- Sterility by the current USP Sterility Test <71>

All results are reported on a lot specific Certificate of Analysis which is available upon request.

BUFFER SYSTEM

mHTF uses a buffering system composed of a 21 mM HEPES (N-2-Hydroxyethylpiperazine-N'-2-ethanesulfonic acid) and 4 mM Sodium Bicarbonate combination. This buffering system provides optimum pH maintenance over the physiologic range (7.2 to 7.4) and does not require the use of a CO₂ incubator.

PROTEIN SUPPLEMENTATION

mHTF does not contain protein components. General laboratory practice includes supplementation when using this medium. The amount of protein supplementation may vary among laboratories and is dependent on the phase of processing/growing the gametes and embryos. Consult your individual laboratory protocols.

DIRECTIONS FOR USE

mHTF Medium should be supplemented with protein, as appropriate, prior to use as a transport medium (prewarmed at 37°C if desired) for the collection of human gametes and for embryo transfer to the patient. This medium is also used for procedures at ambient atmosphere (due to the HEPES buffering capabilities). These procedures include sperm washing and IntraCytoplasmic Sperm Injection (ICSI).

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

STORAGE INSTRUCTIONS AND STABILITY

Store the unopened bottles refrigerated at 2° to 8°C.

Do not freeze or expose to temperatures greater than 39°C.

Duration Following Bottle Opening:

Product should be used within (8) weeks from opening when stored under the recommended conditions of 2° to 8°C.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any medium that is not reddish orange in color, or that shows evidence of particulate matter or cloudiness. Discard the product in accordance with applicable regulations.

mHTF should be tightly capped when used in a CO₂ incubator to avoid pH levels of 7.0 or less.

To avoid problems with contamination, handle using aseptic techniques and discard any excess medium that shows any evidence of contamination after opening.

CONTRAINDICATION

Product contains Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706 Fax: 1 949 261 6522 • www.irvinesci.com

© 2019 FUJIFILM Irvine Scientific Inc. All rights reserved. The FUJIFILM Irvine Scientific logo are trademarks of

FUJIFILM Irvine Scientific, Inc. in various jurisdictions.

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