## Modified Ham's F-10 Basal Medium - HEPES

**Catalog No. 99168**

**100 mL**

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### Glossary of Symbols\*:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature 2-8°C
	Do not re-use
	Do not re-sterilize
	Do not use if package is damaged
	U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	Manufacturer
	CE Mark
	Emergo Europe - Prinsessegracht 20 2514 AP The Hague The Netherlands

\*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

### ENGLISH

**EU Caution:** For Professional Use Only.

#### INDICATION FOR USE

Modified Ham's F-10 Basal Medium - HEPES is intended for use in assisted reproductive procedures which include retrieval, handling and transfer of gametes and embryos.

#### DEVICE DESCRIPTION

Modified Ham's F-10 Basal Medium - HEPES consists of salts, amino acids, 20 mM HEPES, 25 mM bicarbonate buffers and other stable trace components. Supplementation with additional components (sodium lactate, sodium pyruvate, glutamine and gentamicin) are required prior to use. The recombination of the components is performed by the end user and results in a medium which has the same characteristics of pH and osmolality, as it would if made fresh.

#### COMPOSITION

The composition of the complete medium is a modification of the classical tissue culture Ham's F-10 formulation. This is a complex medium of inorganic salts, amino acids and vitamins.

Modified Ham's F-10 Basal Medium - HEPES (catalog #99168) has been specifically formulated without hypoxanthine, and with HEPES for buffering at ambient atmosphere (see section describing Buffer System).

#### COMPOSITION:

<u>Salts and Ions</u>	<u>Amino Acids</u>	<u>Energy Substrates</u>
Sodium Chloride	Alanine	Glucose
Potassium Chloride	Arginine	Inositol
Calcium Chloride	Aspartic Acid	<u>pH Indicator</u>
Magnesium Sulfate	Glutamic Acid	Phenol Red
Ferrous Sulfate	Glycine	<u>Buffer</u>
Zinc Sulfate	Histidine	HEPES
Choline Chloride	Isoleucine	Potassium Bicarbonate
Cupric Sulfate	Leucine	Sodium Bicarbonate
Potassium Phosphate	Lysine	<u>Antioxidant</u>
Sodium Phosphate	Methionine	Thioctic Acid
<u>Vitamins and Minerals</u>	Phenylalanine	<u>Water</u>
Vitamin B-12	Proline	WFI Quality
Folic Acid	Serine	
Riboflavin	Threonine	
Pantothenic Acid	Tryptophan	
Pyridoxine	Tyrosine	
Thiamine	Valine	
Biotin	Cysteine	
Nicotinic Acid Amide	Thymidine	
	Asparagine	

#### QUALITY ASSURANCE

Modified Ham's F-10 Basal Medium - HEPES is a culture medium which is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10<sup>-3</sup>.

Each lot of Modified Ham's F-10 Basal Medium - HEPES is tested for:

- Endotoxin by Limulus Amebocyte Lysate (LAL) methodology
- Biocompatibility by Mouse Embryo Assay (one-cell)
- Sterility by the current USP Sterility Test <71>

All results are reported on a lot specific Certificate of Analysis which is available upon request.

#### BUFFER SYSTEM

Modified Ham's F-10 Basal Medium - HEPES uses a buffering system composed of a 20 mM HEPES (N-2-Hydroxyethylpiperazine-N-2-ethanesulfonic acid) and 25 mM bicarbonate combination. This buffering system provides optimum pH maintenance over the physiologic range (7.2 to 7.4) and does not require the use of a CO<sub>2</sub> incubator.

#### PROTEIN SUPPLEMENTATION

Modified Ham's F-10 Basal Medium-HEPES does not contain protein components. General laboratory practice includes protein supplementation when using this medium. The amount of protein supplementation may vary among laboratories and is dependent on the phase of processing/growing the gametes and embryos. Consult your individual laboratory protocols.

#### DIRECTIONS FOR USE

PREPARATION OF COMPLETE MODIFIED HAM'S F-10 MEDIUM - HEPES

1. Allow bottle of Modified Ham's F-10 Basal Medium - HEPES (#99168) to equilibrate to room temperature.
2. Supplement basal medium aseptically with sterile solutions of sodium lactate, sodium pyruvate, glutamine and gentamicin to obtain final concentrations of 1 mM of each sodium lactate, sodium pyruvate and glutamine, as well as 10 µg/mL of gentamicin.
3. Recap the bottle and mix well by swirling the bottle containing the complete medium.
4. Optional components such as albumin or serum may be added at this point.
5. Note the date of supplementation on the bottle.

Note: Complete Modified Ham's F-10 Medium - HEPES should be tightly capped when used in a CO<sub>2</sub> incubator to avoid pH levels of 7.0 or less.

The complete liquid medium is then ready for use as a transport medium (pre-warmed at 37°C if desired), for the collection of human gametes and for embryo transfer to the patient. This medium is also used for procedures at ambient atmosphere (due to the HEPES buffering capabilities). These procedures include sperm washing and IntraCytoplasmic Sperm Injection (ICSI).

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

#### STORAGE INSTRUCTIONS AND STABILITY

Store the unopened bottles of basal medium at 2° to 8° C. Do not freeze basal medium. The complete medium is stable until 57 days from date of supplementation as a "complete" medium.

Protect from fluorescent light.

#### PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any bottle of medium which shows evidence of particulate matter, cloudiness, or is not pink in color.

To avoid problems with contamination, handle using aseptic techniques, and discard any excess medium that remains in the bottle after the procedure is completed.

Information on known characteristics and technical factors that could pose a risk if the product were to be re-used have not been identified therefore the product is not to be used following the initial use of the container.

If basal medium is supplemented with antibiotic, appropriate precautions should be taken to assure patient is not sensitized to the antibiotic.

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