## Hyaluronidase Solution

**Catalog No. 90101**
**5 x 1 mL**

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
**Glossary of Symbols\*:**

	REF	Catalog Number
	LOT	Lot Number
	STERILE A	Sterilized using aseptic processing techniques (filtration)
		Expiration: Year - Month - Day
		Caution, consult accompanying documents
		Consult instructions for use
		Storage Temperature 2-8°C
		Do not re-use
		Do not re-sterilize
		Do not use if package is damaged
	PHT DEHP DBP	Contains or Presence of Phthalate
		Manufacturer
	Rx Only	U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	CE	CE Mark
	EC REP	Emergo Europe - Prinsessegracht 20 2514 AP The Hague The Netherlands

\*Symbol Reference - EN ISO 15223-1, Medical devices - Symbols to be used with medical device labels, labeling.

**REFERENCES**

Elder, K. and Dale, B. Micromanipulation techniques, in *In Vitro Fertilization*, Cambridge University Press, Cambridge, UK, 2000, p237.

Rombauts, L. and Wood, C.E., Oocyte retrieval, in *Handbook of In Vitro Fertilization*, Trounson, A.O., and Gardner, D.K., Eds., 2<sup>nd</sup> Ed., CRC Press, Boca Raton, FL, 2000, p.78.

**FUJIFILM Irvine Scientific, Inc.**

 2511 Daimler Street, Santa Ana, California 92705 USA  
 Telephone: 1 949 261 7800 • 1 800 437 5706 Fax: 1 949 261 6522 • www.irvinesci.com

© 2019 FUJIFILM Irvine Scientific, Inc. All rights reserved. The FUJIFILM Irvine Scientific logo, ECM, and P-1 are trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions. PN 40694 Rev.14

**ENGLISH**
**EU CAUTION:** For Professional Use Only

**INDICATION FOR USE**

For use in removing cumulus cells surrounding human oocytes (denudation) in preparation for ICSI or other ART procedures.

**DEVICE DESCRIPTION**

An enzymatic solution containing 80 IU/mL of bovine derived hyaluronidase in a HEPES buffered HTF medium supplemented with 5.0 mg/mL (or 0.5 %) therapeutic grade Human Serum Albumin (HSA)\*. Gentamicin Sulfate (10 µg/mL) is included as an antibiotic.

**COMPOSITION:**

<u>Enzyme</u>	<u>Protein Source</u>
Hyaluronidase	Human Serum Albumin
<u>Buffer</u>	<u>pH indicator</u>
HEPES	Phenol Red
<u>Salts and ions</u>	<u>Energy substrate</u>
Sodium Chloride	Sodium Pyruvate
Potassium Chloride	Glucose
Magnesium Sulfate	Sodium Lactate
Potassium Phosphate	<u>Water</u>
Calcium Chloride	WFI Quality
<u>Antibiotic</u>	
Gentamicin Sulfate	

**QUALITY ASSURANCE**

 Hyaluronidase Solution is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10<sup>-3</sup>.

Each lot of Hyaluronidase Solution is tested for:

Endotoxin by Limulus Amebocyte Lysate (LAL) methodology  
 Biocompatibility by Mouse Embryo Assay (one-cell)  
 Enzymatic Activity by the current USP Test  
 Sterility by the current USP Sterility Test <71>  
 All results are reported on a lot specific Certificate of Analysis which is available upon request.

**BUFFER SYSTEM**

 Hyaluronidase Solution uses a buffering system composed of a 21 mM HEPES (N-2-Hydroxyethylpiperazine-N'-2-ethanesulfonic acid) and 4 mM Sodium Bicarbonate combination. This buffering system provides optimum pH maintenance over the physiologic range (7.2 to 7.4) and does not require the use of a CO<sub>2</sub> incubator.

**DIRECTIONS FOR USE**

Pre-equilibrate the Hyaluronidase Solution to 37°C before use.

 (NOTE: Hyaluronidase Solution should be tightly capped if warmed in a CO<sub>2</sub> incubator to avoid pH levels of 7.0 or less.)

Dispense contents of the vial into an appropriate culture dish and place the oocytes in the Hyaluronidase Solution for no more than 30 seconds while stripping off the surrounding cumulus cells with continuous pipetting.

Immediately remove the oocytes with a minimal volume of solution (leaving behind as much cumulus as possible) and transfer to wash medium (Modified HTF + 5 mg/mL HSA). Aspirate the oocytes with a pipette to remove any remaining cumulus and corona cells.

Wash the denuded oocytes 2 more times until completely denuded.

Transfer the oocytes to Complete HTF (Cat. #9922), Complete ECM (Cat. #90140 and #90142), Complete P-1 (Cat. #9926 or #9910), HTF (Cat. #90125) supplemented with protein, P-1 (Cat. #99242) supplemented with protein or Early Cleavage Medium (Cat. #90138) supplemented with protein and identify which oocytes will be used for ICSI (or other) procedures. Denuded oocytes should

be cultured for a minimum of 1 hour before beginning ICSI procedures.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

**STORAGE INSTRUCTIONS AND STABILITY**

Store the unopened bottles refrigerated at 2-8°C.

Do not expose to temperatures greater than 39°C.

Hyaluronidase Solution is stable until the expiration date shown on the bottle label when stored as directed.

**PRECAUTIONS AND WARNINGS**

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any vial of solution which shows evidence of damage, particulate matter, cloudiness or is not reddish-orange in color. Discard the product in accordance with applicable regulations.

To avoid problems with contamination, handle using aseptic techniques and discard any excess medium that remains in the bottle or vial after the procedure is completed.

Do not expose oocytes to the Hyaluronidase Solution for greater than 30 seconds

Information on known characteristics and technical factors that could pose a risk if the product were to be re-used have not been identified therefore the product is not to be used following the initial use of the container.

**EU:** Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/ removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. It is strongly recommended that every time FUJIFILM Irvine Scientific, Inc. Reproductive Media Products culture media are administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

**US:** Hyaluronidase Solution contains Human Serum Albumin (HSA). \*Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal pre-cautions. Donors of the source material have also been screened for CJD.

Hyaluronidase Solution contains hyaluronidase derived from bovine testes that is currently not identified as high risk material with regards to Transmissible Spongiform Encephalopathies (TSE).

**CONTRAINDICATION**

Product contains Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

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