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

Za postopke asistirane reprodukcije.

Glossary of Symbols:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature 2-8°C
	Do Not Re-Sterilize
	Do Not Use If Package Is Damaged
	Phthalate, DBP, DEHP
	Manufacturer
	U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	CE Mark
	Emergo Europe - Westervoortsedijk 60 6827 AT Arnhem The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

 **FUJIFILM Irvine Scientific, Inc.**

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706 • Fax: 1 949 261 6522 • www.irvinesci.com

© 2023 FUJIFILM Irvine Scientific, Inc. All rights reserved. FUJIFILM Irvine Scientific and its logo, Vit Kit, and CryoTip, are registered trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions. All other trademarks are the property of their respective owners.

PN 40695 Rev.25

Effective Date: 31-JUL-2023

Figure 1:

H = Modified HTF with HEPES (90126) + protein
ES = Equilibration Solution (90131)


Figure 2:


Figure 3:

Prepare CryoTip for loading


Figure 4:

Prepare HSV straw for loading


Figure 5:

Prepare Cryolock for loading

Figure 6:


Figure 7a:


Figure 7b:



Figure 8:

ENGLISH

EU CAUTION: For Professional Use Only.

INDICATIONS FOR USE

Vit Kit-Freeze is intended for use in assisted reproductive procedures for vitrification and storage of human oocytes (MII), pronuclear (PN) zygotes through day 3 cleavage stage embryos and blastocyst stage embryos. This kit is designed for use with CryoTip (Catalog #40709), and Vitrification Thaw Kit (Vit Kit-Thaw) for optimal recovery of specimens.

DEVICE DESCRIPTION

Equilibration Solution-ES is a HEPES buffered solution of Medium-199 containing gentamicin sulfate, 7.5% (v/v) of each DMSO and ethylene glycol and 20% (v/v) Dextran Serum Supplement (DSS).

Vitrification Solution-VS is a HEPES buffered solution of Medium-199 containing gentamicin sulfate, 15% (v/v) of each DMSO and ethylene glycol, 20% (v/v) DSS and 0.5 M sucrose.

DSS is a protein supplement consisting of 50 mg/mL therapeutic grade Human Serum Albumin (HSA) and 20 mg/mL Dextran. DSS is used at 20% (v/v) in Vit Kit-Freeze for a final concentration of 10 mg/mL HSA and 4 mg/mL Dextran.

These two solutions are to be used in sequence according to the step-wise microdrop vitrification protocol.

COMPOSITION

Salts & Ions

Sodium Chloride
Sodium Phosphate
Potassium Chloride
Magnesium Sulfate
Sodium Acetate
Calcium Chloride
Choline Chloride
Ferric Nitrate

Buffer

Sodium Bicarbonate
HEPES

pH Indicator

Phenol Red

Amino Acids

Arginine
Glycine
Histidine
Lysine

Proline
Tyrosine
Alanine
Aspartic Acid
Glutamic Acid
Isoleucine
Leucine
Methionine
Phenylalanine
Serine
Threonine
Tryptophan
Valine
Hydroxyproline
Cystine
Cysteine

Antioxidant

Glutathione

Other

Adenine Sulfate
Deoxyribose
Ribose
Guanine
Uracil
Xanthine
Thymine
Hypoxanthine
Adenosine

Vitamins And Minerals

Calciferol
Ascorbic Acid
Aminobenzoic Acid
Nicotinic Acid
Nicotinic Acid Amide
Pantothenic Acid
Riboflavin
Thiamine
Biotin

Pyridoxine
Sodium Bisulfite
Folic Acid
Alpha-Tocopherol

Antibiotics

Gentamicin Sulfate

Energy Substrates

Glucose
Inositol

Protein

Human Serum Albumin

Cryoprotectant

Dextran
Sucrose
Ethylene Glycol
Dimethylsulfoxide

Water

WFI Quality

QUALITY ASSURANCE

The solutions in Vit Kit-Freeze are membrane filtered and aseptically processed according to manufacturing procedures which have been validated.

Each lot of Vit Kit-Freeze receives the following tests:

Solutions and CryoTips.

Endotoxin by Limulus Amebocyte Lysate (LAL) methodology (≤ 0.6 EU/mL)

Mouse Embryo Assay (one-cell) ($\geq 80\%$ expanded blastocyst)

Sterility by the current USP Sterility Test <71> (Pass)

All results are reported on a lot specific Certificate of Analysis which is available upon request.

MATERIALS REQUIRED BUT NOT INCLUDED

- FUJIFILM Irvine Scientific Inc. CryoTip (Catalog #40709) or HSV Straw (Catalog #25246-25251) or Cryolock™ (Catalog #CL R-CT-B, CL-R-CT-G, CL-R-CT-C, CL-R-CT-Y, CL-R-CT-O)
- FUJIFILM Irvine Scientific Inc. Connector (Catalog #40736) or other adaptor
- Sterile petri dishes (50 X 9 mm, Falcon 351006 or equivalent)
- Cryotubes (4.5 mL) or goblets and cryocanes
- Modified HTF - HEPES (Catalog #90126) culture medium supplemented with protein
- Hyaluronidase (Catalog #90101)
- Disposable gloves
- Hamilton GASTIGHT® Syringe (50 μ L, Catalog #80901) or other aspiration tool
- Transfer pipettes (pulled glass pipettes or micropipette tips with an inner tip diameter of $\sim 200\mu$ m)
- Tweezers or forceps
- Impulse Heat Sealer
- SYMS Sealer for HSV Straw
- Stopwatch or timer
- Liquid nitrogen reservoir (dewar or styrofoam container with lid, 1-2 L volume)
- Liquid nitrogen (sufficient volume to achieve 4 inch depth in reservoir)

DIRECTIONS FOR USE

Vit Kit-Freeze component requirements (per application):

- Equilibration Solution (ES):
60 µL for Oocyte Vitrification Protocol
Or
50 µL for Embryo Vitrification Protocol
- Vitrification Solution (VS):
50 µL for Vitrification Protocol
- 1 CryoTip or HSV Straw or Cryolock (stores up to 2 specimens)
- 1 Connector

VITRIFICATION PROTOCOL:

NOTE: Procedures are to be done at room temperature (20-27°C). DO NOT use heated microscope stage for the following procedures. **CAUTION:** Minimize exposure of specimen to light during equilibration in ES and VS solutions.

1. Bring the quantity to be used of ES and VS to room temperature (20-27°C). **NOTE:** Avoid bringing the entire vials of ES and VS to room temperature repeatedly when a partial of the solution is needed each time. It is better to aliquot the quantity to be used and return the vials to 2-8°C right after aliquoting. Modified HTF (HEPES) with protein is also required for oocyte vitrification protocol.
2. Fill the liquid nitrogen reservoir with liquid nitrogen (sufficient to achieve a depth of 4 inches or to completely submerge cryotube on cane) and place close to microscope. Attach a cryotube or goblet (uncapped) to the bottom clamp of a cryocane and submerge in the liquid nitrogen in preparation for storage of the vitrified specimens.
3. Determine the number of specimens to be vitrified.
4. Label each sterile petri dish (or lid) and Cryo storage device with necessary information.
5. Gently invert each vial of ES and VS twice to mix contents before use.
6. Prepare dish with droplets of solutions for Vitrification Procedure as follows:

A. OOCYTE (MII) Vitrification Protocol:

NOTE 1: Retrieved oocytes are to be denuded with Hyaluronidase to confirm they are MII.

NOTE 2: Refer to Section B for embryo vitrification protocol.

1. Aseptically dispense 20 µL drop of culture medium, Modified HTF - HEPES with protein, and ES in close proximity on an inverted lid of sterile petri dish as shown in Figure 1, and place the dish on the microscope stage:
 - one- 20 µL drop of Modified HTF (HEPES with protein)
 - three- 20 µL drops (60 µL total) of ES (ES1, ES2, ES3)
2. Remove the culture dish containing MII oocytes from the incubator and check the quality of the specimens under microscope. Where possible, select only the best quality MII stage oocyte(s).
CAUTION: Minimize the exposure of the specimen(s) to light during equilibration in the H, ES and VS drops.
3. Transfer the oocyte (up to 2 at a time) with minimal volume of medium from the culture dish (in incubator) into the 20 µL drop of H for one minute.
4. Merge the drop of H to ES1 (See Fig. 1, arrow 1) with the tip of the transfer pipette and allow spontaneous mixing of the two solutions to occur for 2 minutes.
5. Then merge the drop of ES2 (arrow 2) to the previously merged drops and leave for 2 minutes.
6. Transfer the oocyte(s) with minimal volume of solution from merged drop to ES3 drop for 6-10 minutes. **Note:** equilibration of oocyte(s) in ES3 is complete when the thickness of the zona pellucida and perivitelline space is equal. The oocyte(s) will settle to the bottom of the drop within 3 minutes.
7. During the equilibration time in ES3:
Aseptically dispense one (1) 50 µL drop of VS 2 minutes prior to complete equilibration and prepare the CryoTip (Fig. 3), HSV Straw (Fig. 4), or Cryolock (Fig. 5) for loading:
NOTE: Carefully examine the vitrification device and tip before starting procedure
 - CryoTip: connect to the Hamilton syringe or appropriate aspiration tool using a connector or adaptor to assure a tight seal. **NOTE:** Keep the metal cover sleeve over the fine pulled tip to protect it until ready to load specimens.
 - HSV Straw: connect the longer end of the blue plastic insertion device to the colored end of the handling rod.
 - Cryolock: detach the cap.
8. The following steps (9-13) should be completed in 80-110 seconds. **CAUTION:** Exposure of specimens to VS should be limited to prevent cytotoxicity. Specimen(s) tend to float in VS, so adjust the focus through the microscope to maintain continuous visualization during exposure and keep the tip of the transfer pipette nearby to assure rapid transfer between VS drops. Refer to Figure 6.
9. After equilibration in ES is complete, draw up some ES into the transfer pipet and transfer the specimen(s) with minimal volume from the drop of ES into the drop of VS for 30 seconds.
10. Load and heat seal the CryoTip as follows (See Figure 7a):
 - Slide the metal cover sleeve up along the CryoTip to expose the fine fragile tip end.
 - Handling the CryoTip and Hamilton syringe while observing under the microscope, carefully aspirate a small volume of VS to the Mark #1 on the CryoTip.
 - Continue observation under the microscope and gently aspirate the specimen with VS to the Mark #2 on the CryoTip.
 - Now observe the CryoTip directly and aspirate more VS to the Mark #3.
 - Specimen must be located between Mark #2 and Mark #3.

- Heat seal (Seal #1) the CryoTip on (or just below) Mark #1, and slide the cover sleeve back down to cover and protect the fine fragile tip.
- Carefully remove the CryoTip from the aspiration tool and adaptor and then heat seal (Seal #2) at the thick end of the CryoTip above the Mark #4.
- Plunge the covered CryoTip directly into liquid nitrogen (cooling at a rate of $-12,000^{\circ}$ C/min) (See Figure 7b).

Load and seal the HSV straw as follows:

- Using a micropipette, carefully deposit the specimen(s) into the gutter of the capillary rod at 1 mm from the end. The drop holding the specimen(s) must be under 0.5 μ l. Maximum of 2 oocytes or embryos for each capillary rod.
- Immediately place the capillary rod and handler into the straw and push until the rectangular portion of the handler comes in contact with the flared end of the straw.
- Slightly pinch the straw between your thumb and finger and remove the insertion device.
- While still holding the straw in place, seal the open end using a SYMS sealer.
- Hold the straw using tweezers in the area of the handling rod.
- Quickly plunge the entire straw into LN₂ vertically. Gently stir the straw in LN₂ for a few seconds so as to avoid formation of an isolating air bubble layer around the straw.

Load the Cryolock as follows:

- Using a micropipette, carefully load a maximum of 2 specimen(s) on the concave surface of the tip (same side of Cryolock logo) about 3 mm (1/8") from the edge of tip (use black mark as a reference), removing any excess of cryoprotectant solution leaving as minimum volume of vitrification media as possible ($\leq 1 \mu$ l).
- Option A: Immediately and before immersing the Cryolock into LN₂, carefully insert the tip into the cap twisting tightly until secure
- Option B: Immediately immerse tip and cap under LN₂. Wait for bubbling to stop to allow for equilibration. Carefully insert the tip into the cap, twisting tightly enough until secure.
NOTE: Option B not cleared for use in the U.S.
- Quickly plunge Cryolock into liquid nitrogen.
NOTE: Always store the Cryolock with the cap facing down.

11. Place the vitrified CryoTip, HSV straw or Cryolock into the submerged LN₂ filled cryotube or goblet (on the cryocane). Cap the cryotube (or goblet) or attach up-side-down with another uncapped cryotube in order to secure the vitrified device in liquid nitrogen.
12. Move the LN₂ reservoir close to the LN₂ cryo-freezer and transfer the cryocane with contents to the cryo-freezer for long-term storage.

B. EMBRYOS (PN to Blastocyst)

Vitrification Protocol:

1. Aseptically dispense one- 50 μ l drop of ES on an inverted lid of Petri dish.
2. Remove the culture dish with embryo(s) from the incubator and check the quality of the specimen(s) under microscope. Where possible, select only the best quality embryo(s) for vitrification.
3. Carefully transfer the specimen (up to two at one time) with a minimal volume of medium from the culture dish to the drop of ES and start the timer.
Embryos should equilibrate in the ES drop slowly by free-fall for 6-10 minutes.
Note: The specimen will shrink and then gradually return to its original size, which indicates that the equilibration is complete.
CAUTION: Minimize the exposure of specimen(s) to light during equilibration in ES and VS drops.
4. During this equilibration time in ES:
 - set up one-50 μ l drop of VS solution as shown in Fig 8 and prepare the CryoTip, HSV Straw, or Cryolock for loading.

Follow protocol as written above (Section A - Oocyte [MI]) Vitrification Protocol) from steps 9 through 12 for exposure to VS solutions, loading of CryoTip, HSV Straw, or Cryolock, plunging in LN₂ and long term storage.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

STORAGE INSTRUCTIONS AND STABILITY

Store the unopened vials refrigerated at 2°C to 8°C. When stored as directed, Vitrification Freeze Kit Solutions are stable until the expiration date shown on the vial labels.

Do not use media for more than eight (8) weeks once containers have been opened.

As human source material is present in the product it may develop some particulate matter during storage. This type of particulate matter is not known to have an effect on product performance.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

As an added precaution during the preparation procedure, we recommend that each CryoTip be carefully examined when taken out of the package. Prior to use, the CryoTips should be examined under a suitable magnification (40x power) for possible damage (such as tip breakages or cracks) that may have occurred during transport.

Do not use any vial of solution which shows evidence of damage, leaking, particulate matter, cloudiness or has shifted color. Discard the product in accordance with applicable regulations.

To avoid problems with contamination, handle using aseptic techniques.

Currently, research literature indicates the long-term effects of vitrification on oocytes and embryos remains unknown.

Do not use any bottle in which the sterile packaging has been compromised.

EU: Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopeia specifications by established processes. It is strongly recommended that every time FUJIFILM Irvine Scientific, Inc. Reproductive Media Products culture media are administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

US: This product contains Human Serum Albumin (HSA). Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal pre-cautions. Donors of the source material have also been screened for CJD.

CONTRAINDICATION

Product contains Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

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pH-Indikator

Phenolrot

Valin
Hydroxyprolin
Cystin
Cystein

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