CHANG Medium BMC

Bone Marrow Culture Medium

Catalog No. 91004
100mL, 500mL

 For *in vitro* diagnostic use.

 Zur *In-Vitro*-Diagnostik.

 Solo per uso diagnostico *in vitro*.

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 Para uso diagnóstico *in vitro*.

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 Za diagnostično uporabo *in vitro*.

REFERENCES

Tijo, JH, and Whang-Peng, J: Direct Preparation of Bone Marrow Cells. Human Chromosome Methodology (JJ Yunis, ed.), Academic Press, New York, 1974.

Hozier, JC, and Lindquist, L: Banded Karyotypes from Bone Marrow: A Clinically Useful Approach. Human Genetics, 53:205-209, 1980.

Williams, DL, et al: A Direct Bone Marrow Chromosome Technique for Acute Lymphoblastic Leukemia. Cancer Genetics and Cytogenetics, 13:239-257, 1984.

Babor, B, and Stossel, T: Hematology, A Pathophysiological Approach, Churchill-Livingstone, Inc., New York 1990.

LeBeau, M: Cytogenetic Analysis of Hematologic Malignant Diseases. ACT Cytogenetics Laboratory Manual, Raven Press, New York, 1991.

 Mitelman, F: Catalog of Chromosome Aberrations in Cancer (4th ed.), Alan Liss, New York, 1991.

Kaplan, B, and Dale, K (eds): The ACT Cytogenetic Symposia, CA 1994.

Mitelman, F, and Heim, S: Cancer Cytogenetics, Wiley-Liss, New York, 1995.

FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706 • Fax: 1 949 261 6522 • www.irvinesci.com

© 2023 FUJIFILM Irvine Scientific, Inc. All rights reserved. The FUJIFILM Irvine Scientific logo and CHANG Medium are trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions.

PN 40746 Rev. 09

Effective Date: 31-JUL-2023

ENGLISH

INDICATION FOR USE

CHANG Medium BMC is intended for use in primary culture of clinical Human Bone Marrow Cultures for karyotyping and other genetic testing of various hematological disorders.

DEVICE DESCRIPTION

CHANG Medium BMC is a ready-to-use medium consisting of RPMI Medium 1640, with FBS, HEPES buffer, L-glutamine, Giant Cell Tumor (GCT) Conditioned Medium and Gentamicin Sulfate. CHANG Medium BMC has been optimized to support efficient cell attachment and growth of bone marrow cells for cytogenetic analysis. No addition of any components prior to culturing bone marrow is required.

COMPONENTS

<u>Amino Acid</u>	<u>Proteins,</u>	<u>Vitamins and trace</u>
Arginine	<u>Horomones, and</u>	<u>elements</u>
Asparagine	<u>Growth Factors</u>	Folic acid
Aspartic Acid	Fetal bovine serum	Nicotinamide
Cystine	(FBS)	Riboflavin
Glutamine		Thiamine
Glutamic Acid		Pantothenic acid
Glycine	<u>pH Indicator</u>	Cobalamin
Histidine	Phenol red	Pyridoxine
Hydroxyproline		Aminobenzoic acid
Isoleucine	<u>Salts & Ions</u>	
Leucine	Sodium chloride	<u>Other</u>
Lysine	Choline chloride	Giant cell tumor
Methionine	Potassium chloride	conditioned medium
Phenylalanine	Magnesium sulfate	(GCT-CM)
Proline	Sodium phosphate	Glutathione
Serine	Calcium Nitrate	Biotin
Threonine		
Tryptophan	<u>Buffers</u>	<u>Energy Substrates</u>
Tyrosine	Sodium bicarbonate	Glucose
Valine	HEPES	Inositol
	<u>Antibiotic</u>	
	Gentamicin Sulfate	
<u>Water</u>		
WFI Quality		

MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

1. Plastic Sterile Centrifuge Tubes and Culture Flasks
2. CO₂ Incubator at 37°C
3. Bench Centrifuge
4. Vortex Mixer
5. Colcemid Stock Solution, 10 µg/mL
6. Potassium Chloride Solution, 0.075 M
7. Fixative Solution, Methanol:Acetic Acid (3:1)

QUALITY ASSURANCE

Several factors including source of specimens, culture conditions and selection of reagents can influence the result obtained. Users are advised to run each new batch of reagent in parallel with reference material of known suitable activity before adoption in routine use. Each lot of CHANG Medium BMC has been performance tested on Clinical Bone Marrow Cultures at an independent Clinical Cytogenetics Laboratory compared to a control medium. Results are reported on a lot specific Certificate of Analysis.

PREPARATION FOR USE

CHANG Medium BMC should be thawed overnight in the refrigerator (2-8°C) then gently mixed to assure homogeneity. Aseptically dispense 10 mL of medium into sterile culture flasks and equilibrate to 37°C for immediate use for bone marrow cultures.

Note: Calcium carbonate crystals commonly form in CHANG Medium BMC. The presence of these crystals has not been shown to cause any detrimental effect on product performance.

DIRECTIONS FOR USE

Sample Preparation:

Use 0.5 to 1.0 mL of sodium heparinized bone marrow aspirate. Lithium heparin, EDTA, or citrate anticoagulants are unsuitable for cytogenetic studies.

- If more than 5 mL of bone marrow aspirate is received, the sample may be hemodilute. Spin the specimen down to isolate the bone marrow fraction.

- If the specimen arrives in transport medium, spin the sample down and remove the transport medium (supernatant). Inoculate using the remaining aspirate fraction.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

Bone Marrow Culture:

Label all culture vessels with patient name, specimen number, and culture type. For each specimen prepare a flask containing:

1. 10.0 mL CHANG Medium BMC.
2. Equilibrate flask to 37°C before inoculation of specimen.
3. Inoculate 0.5 mL (500 µL) of specimen, or the appropriate amount depending upon the white blood cell (WBC) count, into each flask containing 10.0 mL pre-equilibrated CHANG Medium BMC. Add less specimen if WBC is high (> 30,000) or more specimen if WBC is low (< 5,000).
4. Incubate flask at 37°C for 1-2 days.

Harvesting the Cultures:

1. Remove cultures from incubator and gently swirl to resuspend cells.
2. Transfer the contents of the flask to a 15 mL centrifuge tube.
3. Add 100 µL of stock Colcemid (10 µg/mL) to each tube.
4. Cap tubes and mix by inverting.
5. Incubate tubes at 37°C for 20 minutes.
6. After incubation, centrifuge tubes for 8 minutes at 1200 rpm (300 x g).
7. Carefully aspirate supernatant from each tube.
8. Resuspend the cell pellet by gently mixing, or flicking the bottom of the tube with forefinger.
9. VERY SLOWLY add 10 mL of hypotonic solution (0.075 M Potassium Chloride) to each tube while vortexing (on the lowest setting).
10. Let tubes stand at room temperature for 20 minutes (hypotonic treatment).
11. Centrifuge tubes for 8 minutes at 1200 rpm (300 x g).
12. Aspirate supernatant leaving about 1.0 mL of hypotonic solution above cell pellet.
NOTE: Be cautious of fibrous material that may extend from the cell pellet up into the supernatant after centrifugation. The last few mL of supernatant may need to be removed by hand with a Pasteur pipette (not using vacuum aspiration) to avoid aspirating the entire cell pellet into the waste container.
13. Resuspend cell pellet as described in step 8.
14. VERY SLOWLY add 10 mL of 3:1 Methanol:Acetic acid fixative to each tube while vortexing (on the lowest setting).
15. Let tubes stand at room temperature for 20 minutes (first fix).
16. Repeat steps 11 - 13.
17. Add 5 mL of fixative as in step 14.
18. Let tubes stand at room temperature for 10 minutes (second fix).
19. Repeat steps 16-18 (third fix).
20. At this point, fixed cell pellets can be used immediately for slide preparation according to the laboratory's standard protocol or stored in the refrigerator (2-8°C) for future use.

STORAGE AND STABILITY

CHANG Medium BMC should be stored frozen below -10°C until ready to use. CHANG Medium BMC is stable until the expiration date shown on the bottle label when stored frozen. After thawing, any unused product can be dispensed into working aliquots and refrozen for later use, or lightly capped and stored at 2°C to 8°C for up to 30 days. Protect from fluorescent light.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in procedures that include the indicated application for which the device is intended.

CHANG Medium BMC contains FBS and GCT conditioned medium and should be handled with universal laboratory precautions. The medium contains an antibiotic (gentamicin) to reduce the potential for bacterial contamination, but aseptic techniques should always be used when dispensing the medium. Do not use any medium that is not red in color.

Glossary of Symbols*:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature below -10°C
	Do not resterilize
	Do not use if package is damaged
	Manufacturer
	CE Mark
	Emergo Europe - Westervoorsdijk 60 6827 AT Arnhem The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices - Symbols to be used with medical device labels, labeling.

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