CHANG Medium MF Mitogen-Free

Catalog No. 91005

100 mL, 500 mL

For *in vitro* diagnostic use.

Zur *In-vitro*-Diagnostik.

Solo per uso diagnostico *in vitro*.

Para uso diagnóstico *in vitro*.

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
Za diagnostično uporabo *in vitro*.

Glossary of Symbols*:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature below -10°C
	Do not resterilize
	Do not use if package is damaged
	Manufacturer
	CE Mark
	Emergo Europe - Westervoortsedijk 60 6827 AT Arnhem The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices -
Symbols to be used with medical device labels, labeling.

 **FUJIFILM Irvine Scientific, Inc.**

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706 • Fax: 1 949 261 6522 • www.irvinesci.com

REFERENCES

1. Arakaki, D.T. and Sparkes, R.S. (1963), *Cytogenetics*, 2, 57-60.
2. Li, J.G. et Osgood, E.E. (1949). *Blood*, 4, 670
3. Maluish, A.E. and Strong, D. M. (1986). In *Manual of Clinical Laboratory Immunology*, pp. 274-281, 3rd Edition, Eds. N.R. Rose, H. Friedman and J.L. Fahey. *American Society for Microbiology*, Washington, D.C.
4. Nowell, P. C. (1960). *Cancer Research* 20, 462-466.
5. Waithe, W. I. And Hirschron, K. (1978). In *Handbook of Experimental Immunology*, 3rd Edition, D. M. Weir, Ed. Blackwell Scientific Publications, Oxford, Chapter 26- Lymphocyte responses to activators.
6. Watt, J.L. and Stephen G.S. (1986). In *Human Cytogenetics: a practical approach*, pp. 39-55. Eds.D. E. Rooney and B.H. Czepulkowski, IRL Press Ltd., Oxford.
7. A proposed Standard System of Nomenclature of Human Mitotic Chromosomes. (1960) *Lancet*, I, 1063.

ENGLISH

INDICATION FOR USE

CHANG Medium MF is a mitogen-free, ready to use medium for use in culturing peripheral blood and other specimens for purposes of cytogenetic analysis.

DEVICE DESCRIPTION

CHANG Medium MF consists of RPMI containing 20% FBS, 2 mM glutamine, 20 mM HEPES buffer and the antibiotic Gentamicin Sulfate. It may require the addition of mitogenic agents, such as phytohemagglutinin (PHA) to optimize the growth of peripheral blood and other cells. The required concentration of PHA (or other mitogens) should be determined by the individual laboratory.

COMPONENTS

<u>Amino Acid</u>	<u>Water</u>	<u>Salts & Ions</u>
Arginine	WFI Quality	Sodium chloride
Asparagine	<u>Proteins, Hormones, and Growth Factors</u>	Choline chloride
Aspartic Acid	Fetal bovine serum (FBS)	Calcium nitrate
Cystine	<u>pH Indicator</u>	Potassium chloride
Glutamic Acid	Phenol red	Magnesium sulfate
Glutamine	<u>Other</u>	Sodium phosphate
Glycine	Biotin	<u>Antibiotic</u>
Histidine	Hydroxyproline	Gentamicin Sulfate
Isoleucine	Glutathione	<u>Vitamins and trace elements</u>
Leucine	<u>Buffers</u>	Folic acid
Lysine	Sodium bicarbonate	Nicotinamide
Methionine	HEPES	Riboflavin
Phenylalanine	<u>Energy Substrates</u>	Thiamine
Proline	Glucose	Pantothenic acid
Serine	Inositol	Cobalamin
Threonine		Pyridoxine
Tryptophan		Aminobenzoic acid
Tyrosine		
Valine		

QUALITY ASSURANCE

STERILITY

Serum used in the production of CHANG Medium MF has been tested for viral contamination per CFR Title 9 Part 113.53. It has also been screened for mycoplasma contamination. CHANG Medium MF is sterilized by filtration through a 0.1 micron filter. Samples of CHANG Medium MF are tested for possible bacteriological contamination following the sterility testing protocol described in the current USP Sterility test <71>.

Several factors including source of sample, culture conditions and selection of reagents can influence the result obtained. Users are advised to run each new batch of reagent in parallel with reference material of known suitable activity before adoption in routine use.

Each Lot of CHANG Medium MF is evaluated using peripheral blood for mitotic index, chromosomal length and quality compared to a control medium. Results are reported on a lot specific Certificate of Analysis.

MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

1. Phytohemagglutinin (PHA) Solution (9 mg/mL), Sterile
2. Phenol-Free Heparin
3. Colcemid 25 µg/mL Solution
4. Potassium Chloride Solution, 75 mM
5. Acetic Alcohol, 1 part Glacial Acetic Acid: 3 parts Methanol (Analytical Reagent Grade)
6. Giemsa or 2% Acetic Acid-orcein
7. Chromic acid
8. Mountant
9. Glass Slides and Coverslips
10. Plastic Centrifuge Tubes
11. CO₂ Incubator
12. Bench Centrifuge
13. Vortex Mixer
14. Light Microscope

SPECIMEN COLLECTION AND HANDLING

The success of whole blood cultures for cytogenetic analysis is influenced by the level of lymphocytes with normal function at the time of sampling. As this level can be affected by infection and drugs, wherever possible subjects for cytogenetic studies should not have taken drugs for 7 days prior to collection of blood for tests. Similarly, the mitotic index may be greatly reduced during the anergic phases of certain diseases (e.g. Hodgkin's disease, sarcoidosis, etc.) and, to a lesser degree, in normal individuals during the later stages of pregnancy. Preservative must not be added to blood samples for lymphocyte culture. Aseptic techniques are essential.

Blood samples should be tested without delay whenever possible. If absolutely necessary, they may be stored at 2°C to 8°C for no longer than 48 hours.

PREPARATION FOR USE

Thaw overnight in a refrigerator (2–8°C) then gently mix to assure homogeneity. Aseptically dispense 10 mL of medium into sterile culture flasks and equilibrate to 37°C for immediate use.

Note: Calcium Carbonate crystals commonly form in CHANG Medium MF. The presence of these crystals has not been shown to cause any detrimental effect on product performance.

DIRECTIONS FOR USE

Whole blood or separated leukocytes have been used for cytogenetic studies, but the former is the simplest and most widely used in routine studies. As with any cell culture procedure, optimal results are dependent on establishing adequate culture conditions. Since the relative content of active PHA may vary slightly between different batches, it may be beneficial to test two concentrations of PHA.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

I. Sample Preparation:

Obtain 5-10 mL for adults, of fresh blood in sodium heparin, and 2-3 mL for pediatrics.

- Whole blood should be received at room temperature and mixed by inversion.
- Lymphocytes from whole blood can be stimulated with Phytohemagglutinin (PHA) and may be cultured with a synchronization technique.

II. Whole Blood Culture

Label all culture vessels with patient name, specimen number, and culture type.

1. Before inoculation of specimen bring CHANG Medium MF to ambient temperature.
2. Reconstitute PHA by adding 5 mL of sterile distilled water, using a sterile syringe.
3. Prepare required 5 mL of CHANG Medium MF necessary for each culture vessel, using aseptic technique.
4. Add 0.1 mL of reconstituted PHA per culture vessel.
5. Inoculate 0.3 mL of sample per culture. If the patient is an infant (<1mo old), inoculate with 0.2 mL of sample.
6. Each individual laboratory should determine the number of cultures to set up depending on the clinical indication and age of the patient. An additional 48 hour culture is to be initiated on specimens from newborns (<1mo old) without synchronization.
7. Incubate cultures at 35 - 39°C, 5 - 8% CO₂ atmosphere until ready for harvest.
8. For synchronization: After 48 hours of incubation add 50µL per 5mL of working solution Methotrexate to each culture to be harvested at 72 hours. Add 100µL per 5 mL, working solution Thymidine to each culture 18- 19 hours after adding Methotrexate (5 - 6 hours before harvest).

III. Harvesting the Cultures:

1. Remove culture ready for harvest from the incubator and gently swirl flask to re-suspend cells.
2. Transfer the contents of each flask to a 15 mL centrifuge tube.
3. Add 40 µL of stock Colcemid (10 µg/mL) to each culture tube. Tightly cap tubes and mix gently by inverting.
4. Incubate tubes at 35 - 39°C, for 45 minutes.
5. After incubation, centrifuge tubes for 8 minutes at 1,000 rpm.
6. Carefully aspirate supernatant from each tube using a vacuum aspirator, with solvent trap. Be careful to not aspirate pellet.
7. Resuspend pellet by tapping bottom or side of each tube with finger.
8. Initiate a 20 minute timer.
9. Add 3 - 4 mL dropwise, of prewarmed (35 - 37°C) hypotonic solution (0.075M Potassium Chloride).
10. Tightly cap tube and mix gently by tapping bottom or side of the tube with finger.
11. Add 5 - 6 mL, drop-wise, of pre-warmed (35 - 37°C) hypotonic solution. Tightly cap tube and invert tube.
12. Repeat Steps 8 -10 for each tube.
13. Using a water bath, allow tubes to stand at 35 - 37°C. Invert tubes once at midpoint of 20 minute timer.
14. Add 1 mL of fresh 3:1 Carnoy's Fixative to each tube, at the end of the 20 minute timer. Tightly cap and invert each tube. (This is the Pre-Fixative step)
15. Centrifuge tubes for 8 minutes at 1,000 rpm.
16. Aspirate supernatant from each tube, leaving about 1 mL above the cell pellet. Be careful to not aspirate pellet. Be cautious of fibrous material that may extend from the cell pellet up into the supernatant after centrifugation. The last few mL of supernatant may need to be removed by hand with a Pasteur pipette (not using vacuum aspiration) to avoid aspirating the entire cell pellet into the waste container.
17. Resuspend cell pellet, as described in Step 7.

18. Add 3 - 4 mL dropwise of fresh 3:1 Carnoy's Fixative.
19. Add remaining fixative up to 7 mL.
20. Repeat steps 16 - 18 for each tube.
21. Let stand for 10 minutes at room temperature. (This is the First fixative step).
22. Centrifuge tubes for 8 minutes at 1,000 rpm.
23. Aspirate supernatant leaving about 1 mL above pellet. Resuspend cell pellet.
24. Add fix up to 7 mL. Centrifuge tubes for 8 minutes at 1,000 rpm. (Second fixative step).
25. Repeat steps 23 - 24. (Third fixative step).
26. At this point, fixed cell pellets can be used immediately for slide preparation according to the laboratory standard protocol or stored in the refrigerator (2 - 8°C) or freezer for future use.

CLINICAL APPLICATIONS

The modal human chromosome number is 46. Human chromosomes have been classified according to their length and the position of the centromere (Denver Classification). Aberrations of chromosomal constitution have been associated with a number of congenital disorders such as Down syndrome (typically with an additional small autosome) and syndromes associated with indeterminate sexuality (Turner's syndrome, Klinefelter's syndrome and others, where the sex chromosomes are found to be abnormal). An acquired chromosomal abnormality can be detected in a proportion of the leukocytes in chronic myelocytic leukemia (the "Philadelphia" chromosome) and the progress of treatment may be assessed by following this marker. As the limit of tolerance is approached in radiation therapy, there is a marked increase in the proportion of cells with bizarre chromosomal constitution and the appearance of these cells has been used as a guide to dosage.

STORAGE AND STABILITY

Store CHANG Medium MF frozen, below -10°C, until ready to use. CHANG Medium MF is stable until the expiration date shown on the bottle label when stored frozen. After thawing, any unused product can be dispensed into working aliquots and refrozen (maximum of two times) for later use, or tightly capped and stored at 2-8°C for up to 30 days. Protect from fluorescent light.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in procedures that include the indicated application for which the device is intended.

CHANG Medium MF contains FBS and should be handled with universal laboratory precautions. The medium contains an antibiotic (gentamicin sulfate), to reduce the potential of bacterial contamination, but aseptic techniques should always be used when dispensing the media. Do not use any medium that is not red in color.

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