## 10% PVP Solution with HSA

Catalog No. 90123

5 x 0.5 mL Kits

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
Za postopke asistirane reprodukcije.

### Glossary of Symbols\*:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature 2-8°C
	Do not re-use
	Do not re-sterilize
	Do not use if package is damaged
	Contains or Presence of Phthalate
	Manufacturer
	U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	CE Mark
	Emergo Europe - Prinsessegracht 20 2514 AP The Hague The Netherlands

\*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

### REFERENCES

Attee S, Pool T, Martin J: A Simple Approach to Intracytoplasmic Sperm Injection. Fertil Steril 63:3, 652-665, 1995.

Van Steirteghen A.C, et. al.: High Fertilization and Implantation Rates After Intracytoplasmic Sperm Injection. Human Reproduction 8:7, 1061-1066, 1993.

Palermo, G, et. al.: Intracytoplasmic sperm injection: a novel treatment for all forms of male factor infertility. Fertil Steril 63:1231-1240, 1995.

### ENGLISH

**EU Caution:** For Professional Use Only

#### INDICATION FOR USE

10% PVP Solution with HSA\* is intended for use in assisted reproductive procedures which include human gamete and embryo manipulation. These procedures include the use of PVP solution for immobilizing sperm for ICSI procedures.

#### DEVICE DESCRIPTION

Polyvinylpyrrolidone (PVP) dissolved in mHTF (an isotonic HEPES buffered medium) to make a 10% (w/v) solution. Contains 5 mg/mL Human Serum Albumin.

#### COMPOSITION

Salts and Ions	Energy Substrates
Sodium Chloride	Glucose
Potassium Chloride	Sodium Pyruvate
Magnesium Sulfate	Sodium Lactate
Potassium Phosphate	Protein Source
Calcium Chloride	Human Serum Albumin
Buffers	Polymer
HEPES	Povidone
Sodium Bicarbonate	Water
	WFI Quality

#### QUALITY ASSURANCE

PVP solutions are aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10<sup>-3</sup>.

Each lot of PVP Solution is tested for:

- Endotoxin by Limulus Amebocyte Lysate (LAL) methodology
- Sterility by the current USP Sterility Test <71>

All results are reported on a lot specific Certificate of Analysis which is available upon request.

#### DIRECTIONS FOR USE

The 10% PVP Solution is recommended for normal and high motility sperm specimens.

Place aliquot of sperm into liquid PVP solution, 10%. Fill ICSI pipette with PVP solution and then capture single, immobilized sperm according to standard laboratory protocols. Sperm is now ready for ICSI procedure.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

#### STORAGE INSTRUCTIONS AND STABILITY

Store the unopened vials refrigerated at 2°C to 8°C.

Do not freeze or expose to temperatures greater than 39°C.

PVP solutions are stable in unopened vials until the expiration date shown on the kit label when stored as directed.

#### PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures that include the indicated application for which the device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any vial of medium which shows evidence of particulate matter, cloudiness or is not colorless.

Do not use any bottle in which the sterile packaging has been compromised.

To avoid problems with contamination, handle using aseptic techniques.

Information on known characteristics and technical factors that could pose a risk if the product were to be re-used have not been identified therefore the product is not to be used following the initial use of the container.

Sperm exposed to PVP solution should be used as soon as possible. Extended exposure to PVP may significantly impact sperm viability.

EU: Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. It is strongly recommended that every time this product is administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

US: This product contains Human Serum Albumin (HSA). Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal pre-cautions. Donors of the source material have also been screened for CJD.

<b>DEUTSCH</b>	<b>ITALIANO</b>	<b>ESPAÑOL</b>
----------------	-----------------	----------------

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