

## Ultrafiltered Yeast Hydrolysate

---

Catalog #	Product	Size
96863	Ultrafiltered Yeast Hydrolysate	1 kg

### Intended Use

For research or further manufacturing use only. Ultrafiltered Yeast Hydrolysate may be used as a supplement to improve cell growth and protein production in a variety of cell cultures.

### Quality Assurance

All quality control test results are reported on a lot specific Certificate of Analysis (COA), which is available at [www.irvinesci.com](http://www.irvinesci.com) or upon request.

### Storage Instructions and Stability

Handle using aseptic techniques after sterile filtration to avoid contamination. Store dry at 15-30°C. Do not use after the assigned expiration date.

### Suggested Instructions for Use

---

1. Add 2-6 g/L of hydrolysate supplement to 1 L of culture medium (optimal amount may vary depending on cell line and chosen hydrolysate).
2. Dissolve completely.
3. Sterile filter the supplemented medium by 0.2 µm membrane filtration.

## Related Products

---

Catalog #	Product	Size
91128	BalanCD CHO Growth A, Liquid	1L
94120	BalanCD CHO Growth A, Powder	10L
94110	IS CHO-CD XP, Powder	10L
96850	Optimized Hydrolysate Mixture for CHO Cells	made-to-order
96857	Ultrafiltered Soy Hydrolysate	made-to-order
91104	IS MAB-CD	1L

## Technical support

---

### CONTACT US

For more information or assistance, contact Customer Service at:

- Email: [tmrequest@irvinesci.com](mailto:tmrequest@irvinesci.com)
- Direct line: +1 800 577 6097

### WEBSITE RESOURCES

Visit the website at [www.irvinesci.com](http://www.irvinesci.com) for technical resources and information including:

- Safety Data Sheets (SDS)
- COAs (when available)
- FAQs
- Product literature
- Complete list of offices and contact information by country

---

#### **FUJIFILM Irvine Scientific, Inc.**

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706

Fax: 1 949 261 6522 • [www.irvinesci.com](http://www.irvinesci.com)