Continuous Single Culture (CSCM) with Gentamicin

Catalog No. 90164

60 mL

For assisted reproductive procedures.
Für assistierte Reproduktionsverfahren.
Per tecniche di riproduzione assistita.
Para utilización en técnicas de reproducción asistida.
Pour les techniques de procréation médicalement assistée.
Para técnicas de reprodução assistida.
Για διαδικασίες υποβοηθούμενης αναπαραγωγής.
Pro postupy asistované reprodukce.
Til assisteret reproduktionsbehandling.
Avusteisiin lisäätymismenetelmiin.
Ar palīdzītkēļiem veicamām reprodūktīvām procedūram.
Voor geassisteerde voortplantingsprocedures.
Do procedur wspomagane go rozrodu.
Pentru proceduri de reproducere asistată.
För procedurer för assisterad befruktning.
Kasutamiseks abistatud viljastamisprotseduurides.
Asszisztált reprodukciós eljárásokhoz.
Skirta pagalbinio apvaisinimo procedūroms.
Yardımcı üreme işlemleri içindir.
Na postupy asistovanej reprodukcie.
За процедури за асистирана репродукција.
За postupke potpomognute oplodnje.
Għal proceduri ta' riproduzzjoni assistita.
За postopke asistirane reprodukcije.

REFERENCES

Biggers, JD. and Racowsky, C. *The development of fertilized human ova to the blastocyst stage in KSOM[®] medium: is a two-step protocol necessary?* RBMOnline, 5:133-140, 2002.
Pool, TB. *Recent advances in the production of viable human embryos in vitro.* RBMOnline, 4:294-302, 2002.
Biggers, JD. *Thoughts on embryo culture conditions.* RBMOnline, 4 (suppl.1):30-38, 2001.
Lane, M., Hooper, K., and Gardner, DK. *Effect of essential amino acids on mouse embryo viability and ammonium production.* J. Asst. Reprod. Genet. 18: 519-525, 2001
Biggers, JD. and McGinnis, LK. *Evidence that glucose is not always an inhibitor of mouse preimplantation development in vitro.* Hum. Reprod 16:153-163, 2001.
Devreker, F., Van den Bergh, M., Biramane, J., Winston, RML., Englert, Y., and Hardy, K. *Effects of taurine on human embryo development in vitro.* Hum. Reprod. 14: 2350-2356, 1999.

FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705 USA
Telephone: 1 949 261 7800 • 1 800 437 5706 Fax: 1 949 261 6522 • www.irvinesci.com

© 2019 FUJIFILM Irvine Scientific, Inc. All rights reserved. The FUJIFILM Irvine Scientific logo, and Continuous Single Culture are trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions. PN 40944 Rev.8

ENGLISH

EU CAUTION: For Professional Use Only.

INDICATION FOR USE

Continuous Single Culture (CSCM) is intended for use in assisted reproductive procedures which include gamete and embryo manipulation. These procedures include the use of CSCM as a culture medium from fertilization through day 5/6 of development.

DEVICE DESCRIPTION

Continuous Single Culture (CSCM) is optimized to be used in an uninterrupted culture system without dish change or medium renewal. This product contains 10 µg/mL of Gentamicin Sulfate

COMPOSITION:

<u>Salts & Ions</u>	<u>Amino Acids</u>
Sodium Chloride	Alanine
Potassium Chloride	Asparagine
Potassium Phosphate	Aspartic Acid
Calcium Chloride	Glutamic Acid
Magnesium Sulfate	Glycine
<u>Buffer</u>	Proline
Sodium Bicarbonate	Serine
<u>Energy Substrates</u>	Arginine
Sodium Pyruvate	Cystine
Dextrose	Histidine
Sodium Lactate	Isoleucine
<u>Antioxidant</u>	Leucine
EDTA	Lysine
Sodium Citrate	Methionine
<u>Dipeptide</u>	Phenylalanine
Alanyl-glutamine	Threonine
<u>pH Indicator</u>	Tryptophan
Phenol Red	Tyrosine
<u>Antibiotic</u>	<u>Water</u>
Gentamicin Sulfate	WFI Quality

QUALITY ASSURANCE

CSCM is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10⁻³.

Each lot of CSCM is tested for:

- Endotoxin by Limulus Amebocyte Lysate (LAL) methodology
- Biocompatibility by Mouse Embryo Assay (one-cell to expand blastocysts after 96 hours)
- Sterility by the current USP Sterility Test <71>
- Human sperm survival assay (HSSA) at 24 hours

All results are reported on a lot specific Certificate of Analysis which is available upon request.

BUFFER SYSTEM

CSCM uses sodium bicarbonate as a buffering system. This is specifically designed for use in a CO₂ incubator.

DIRECTIONS FOR USE

PROTEIN SUPPLEMENTATION

CSCM does not contain protein components. General laboratory practice includes protein supplementation when using this medium. The amount of protein supplementation may vary among laboratories and is dependent on the phase of processing/growing the gametes and embryos. Consult your individual laboratory protocols.

The following are recommendations for protein supplementation based upon the indications for use of the CSCM:

For Fertilization and Embryo Culture:

When using FUJIFILM Irvine Scientific, Inc. Human Serum Albumin (HSA), a 100 mg/mL solution, use at 5 mg/mL. For 10 mL of medium, add 0.5 mL of HSA solution to 9.5 mL of the medium. When using FUJIFILM Irvine Scientific, Inc. Serum Substitute Supplement (SSS), a 50mg/mL protein Solution, use at 10% (v/v). For 10 mL of medium, add 1.0 mL SSS to 9.0 mL of medium.

EQUILIBRATION

CSCM (supplemented with protein) should be pre-warmed to 37°C and equilibrated to the desired pH overnight in a 5-6% CO₂ incubator prior to use. A sufficient volume of protein supplemented medium is required so that oocyte recovery, insemination and embryo culture dishes can be prepared.

The following are general procedures for the indications for use of CSCM.

Fertilization:

On the day before oocyte retrieval, prepare oocyte collection and insemination dishes with pre-supplemented CSCM overlaid with oil and pre-equilibrated overnight to 37°C in a CO₂ incubator. Immediately upon oocyte collection and identification, place oocytes into the oocyte collection dish with pre-equilibrated medium and return to the incubator for the desired period (1-4 hours) prior to insemination by conventional IVF or ICSI.

Conventional IVF (use insemination dishes):

- It is recommended to aseptically dispense 50,000-100,000/mL motile sperm per microdroplet containing 1-3 oocytes.
- Return the insemination dish to the incubator and check for normal fertilization 16-20 hours post insemination.

Intracytoplasmic Sperm Injection (ICSI):

- Following at least 1 hour post oocyte denuding (and no more than 4 hours following oocyte retrieval), remove denuded oocytes from incubator and inseminate with sperm per standard ICSI protocol for your individual laboratory.
- Immediately following insemination, place 1-3 inseminated oocytes into a fresh drop of the pre-equilibrated insemination dish, return dish to the incubator and check for normal fertilization 16-20 hours post insemination.

Embryo Culture:

On the day of fertilization (one day prior to fertilization assessment), prepare embryo culture dishes with pre-supplemented CSCM overlaid with oil and pre-equilibrate overnight to 37°C in a 5-6% CO₂ incubator.

Following fertilization assessments with the identification of the presence of normal fertilization (two pronuclei and two polar bodies), transfer 2PN zygotes into the pre-equilibrated CSCM culture dish previously prepared. It is recommended to allow the embryos to grow in a continuous, uninterrupted culture system without changing medium, until the desired developmental stage is reached (up to day 5/6 of development).

If medium change is desired for embryo culture beyond day 3, after 48 hours of embryo culture (of the fertilized embryos), the embryos should be transferred into a new dish of fresh pre-equilibrated CSCM (pre-supplemented with protein).

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

STORAGE INSTRUCTIONS AND STABILITY

Store the unopened bottles refrigerated at 2° to 8°C.

Do not freeze or expose to temperatures greater than 39°C.

Duration Following Bottle Opening:

The product without protein supplement should be used within eight (8) weeks from opening.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any bottle of medium which shows evidence of particulate matter, cloudiness or is not pale orange in color.

CAUTION: To avoid problems with contamination, handle using aseptic techniques and discard any excess medium that shows any evidence of contamination after opening.

Not for injection use.

CONTRAINDICATION

CSCM contains the antibiotic Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

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