**EU CAUTION:** For Professional Use Only.

**INDICATION FOR USE**

Multipurpose Handling Medium - Complete (MHM-C) is intended for use in assisted reproductive procedures which involve the manipulation of human gametes or embryos. Specifically, MHM-C is indicated for use as an oocyte retrieval medium during ovarian follicle aspiration procedures (not for flushing ovarian follicles), washing sperm prior to IVF and ICSI fertilization procedures, and for transport of the embryo to the uterus during embryo transfer procedures.

**DEVICE DESCRIPTION**

MHM-C is dual buffered solution (HEPES and MOPS). MHM contains Human Serum Albumin (HSA) for a final total protein concentration of 5 mg/mL and the antibiotic Gentamicin Sulfate (9.5 µg/mL).

**COMPOSITION:**

<u>Salts &amp; Ions</u>	<u>Energy Substrates</u>
Potassium Phosphate	Sodium Lactate
Sodium Chloride	Sodium Pyruvate
Potassium Chloride	Glucose
Calcium Chloride	<u>Protein Source</u>
Magnesium Sulfate	Human Serum Albumin
<u>Amino Acids</u>	<u>Antibiotic</u>
Glycine	Gentamicin Sulfate
Taurine	<u>pH Indicator</u>
<u>Buffer</u>	Phenol Red
HEPES	<u>Water</u>
Sodium Bicarbonate	WFI Quality
MOPS	

**QUALITY ASSURANCE**

MHM-C is a handling medium which is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10<sup>-3</sup>.

Each lot of MHM-C is tested for:

- Endotoxin by Limulus Amebocyte Lysate (LAL) methodology (≤ 1.00 EU/mL)
- Biocompatibility by Mouse Embryo Assay (one-cell at ≥80% expanded blastocyst 96h).
- Sterility by the current USP Sterility Test <71>
- Human Sperm Survival Assay (HSSA) (≥70% motility at 24h).

All results are reported on a lot specific Certificate of Analysis which is available upon request.

**BUFFER SYSTEM**

MHM-C uses a buffering system composed of a HEPES (N-2-Hydroxyethylpiperazine-N'-2-ethanesulfonic acid), MOPS (3 Morpholinopropane-1-sulfonic acid) and Sodium Bicarbonate combination. This buffering system provides pH maintenance over the physiologic range (7.2 to 7.4) and does not require the use of a CO<sub>2</sub> incubator.

**PROTEIN SUPPLEMENTATION**

MHM-C contains 5 mg/mL HSA. No protein supplementation is required prior to use.

**DIRECTIONS FOR USE**

The following are general procedures for the indications of use for the MHM-C.

**Sperm Washing:**

The general procedure for washing sperm from its surrounding seminal fluid includes:

1. Bring medium to room temperature or 37°C.
2. Allow the semen to liquefy at room temperature for 20 to 30 minutes.
3. Using aseptic techniques, transfer the liquefied semen into a sterile 10 mL conical centrifuge tube and add 2 to 3 volumes of room temperature MHM-C (for example, a 2 mL semen sample requires 4 to 6 mL of medium). Should the volume of the sperm medium mixture be greater than 5 mL, divide into two sterile conical

centrifuge tubes, minimizing the volume per tube to 4 – 6 mL, the recovery of sperm will be maximized. Samples having high viscosity may require a further processing to ensure total sperm recovery. (See the Special Processing Considerations section).

4. Centrifuge the tubes at ambient temperature for 10 minutes using a g-force of 200 - 300 x g.
5. Using a sterile pipette, remove and discard the supernatant above the "sperm pellet" by aspiration. The sperm should then be resuspended by gently flicking the tube externally with the index finger. (Note: Do not use a vortex mixer for this step). Resuspend the sperm in 1 to 2 mL of fresh medium, recap and gently mix by inversion. Samples which were fractionated for the first centrifugation step should now be recombined into one tube.
6. Recentrifuge as in Step 4.
7. Using a sterile pipette, remove and discard the supernatant and resuspend the sperm pellet gently by manual agitation. Add fresh medium to a final volume of 0.5 mL. The sperm are ready for assisted reproductive procedures. (Note: The total volume of the nonpregnant uterus is 15 - 56 mL).

**SPECIAL PROCESSING CONSIDERATIONS Processing the highly viscous semen sample:**

Some samples are naturally highly viscous even after liquefaction. These samples have the consistency of heavy syrup and may be among the most difficult to process.

1. After the medium is added to an ejaculate, aspirate and expel the mixture gently using an 18 gauge needle and syringe. This will "shear" some of the viscous mucus.
2. Limit the amount of medium-sperm mixture from Step 3 to 5 mL per centrifuge tube for the first centrifugation step.
3. If after preprocessing the sample with the needle and syringe (Step 1), the sperm do not "pellet" in a normal manner (the sperm will appear as a "cloudy fiber" attached to the bottom of the centrifuge tube), carefully aspirate as much of the supernatant as possible without disrupting the "cloudy sperm fiber" using a sterile needle and syringe. This can be done by keeping the beveled edge of the needle firmly against the wall of the centrifuge tube and slowly start aspiration from the top of the tube downward. When as much of the supernatant as possible has been removed, add 2 or 3 mL of fresh medium. Repeat the process of drawing the mixture through the 18 gauge needle and syringe. Recentrifuge the mixture. The sperm should pellet normally after the second processing.
4. On subsequent sample collections, the patient should be requested to produce a split ejaculate which will minimize the viscosity in the sperm rich portion of the sample.

**Oocyte Retrieval (not for flushing ovarian follicles):**

MHM-C may be supplemented with quality tested pharmaceutical grade heparin (2.5 - 10 units/mL) to reduce clotting of the follicular aspirates containing blood.

1. Bring MHM-C to room temperature or 37°C.
2. The collected follicle aspirates should be transferred into an empty sterile dish.
3. Identify the oocytes and remove them from follicular fluid and possible blood contamination using sterile pipettes pre-rinsed with MHM-C.
4. Rinse the oocytes in pre-warmed MHM-C.
5. Place oocytes into an equilibrated culture medium for further handling.

**Embryo Transfer:**

Transfer of embryos from culture medium on day 3 or day 5:

1. On day 3 or day 5 following assessment of embryos for development, bring MHM-C to room temperature or 37°C.

2. Set up one sterile wash dish containing pre-warmed MHM-C for each set of embryos.
3. Place 1.0 mL of the pre-warmed MHM-C into the well of a sterile 1-well dish.
4. Place the wash dish onto a heated stage.
5. Wash the embryos in the wash dish by picking up the embryos 2 - 3 times and moving them around in minimal volume of the pre-warmed MHM-C within the well.
6. After washing the embryos are ready for transfer into the patient.

**STORAGE INSTRUCTIONS AND STABILITY**

Store the unopened bottles refrigerated at 2° to 8°C. Do not freeze or expose to temperatures greater than 39°C. Duration following bottle opening: 100 mL & 500 mL bottles- product should be used within five (5) weeks of opening when stored under the recommended conditions of 2° to 8°C. 12 mL bottle- product should be used within one (1) week of opening when stored under the recommended conditions of 2° to 8°C.

**PRECAUTIONS AND WARNINGS**

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended. This device is not for use in ovarian follicles flushing procedure. This medium is not for use in oocyte flushing procedures.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any bottle of medium which shows evidence of particulate matter or cloudiness.

Multipurpose Handling Medium - Complete (MHM-C) should be tightly capped when used in a CO<sub>2</sub> incubator to avoid pH levels of 7.0 or less.

To avoid problems with contamination, handle using aseptic techniques and once removed from the bottle discard any excess medium that remains after the procedure is completed.

Not for injection use.

**EU:** Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applied to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. It is strongly recommended that every time MHM-C are administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of product.

**US:** MHM-C contains Human Serum Albumin (HSA). Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal precautions. Donors of the source material have also been screened for Creutzfeldt-Jakob Disease (CJD).

**CONTRAINDICATION**

Product contains Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

IrvineScientific

**Multipurpose Handling Medium-Complete (MHM-C)**

with Gentamicin and HSA

Catalog No. 90166

100 mL, 500 mL, 12 x 12 mL

For assisted reproductive procedures.

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Sterilized using aseptic processing techniques (filtration)

Expiration: Year - Month - Day

Caution, consult accompanying documents

Consult instructions for use

Storage Temperature 2-8°C

Do not resterilize

Do not use if package is damaged

Manufacturer:

U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

CE Mark

Emurgo Europe - Westervoortseidijk 60 6827 AT Arnhem The Netherlands

\*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

**FUJIFILM Irvine Scientific, Inc.**  
2511 Daimler Street, Santa Ana, California 92705 USA  
Telephone: 1 949 261 7800 • 1 800 437 5706 Fax: 1 949 261 6522 • www.irvinesci.com

© 2023 FUJIFILM Irvine Scientific, Inc. All rights reserved. The FUJIFILM Irvine Scientific logo, Multipurpose Handling Medium and MHM are registered trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions.

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