

FUJIFILM



IrvineScientific

CHANG Amnio with Gentamicin and L-Glutamine

Catalog No. 99473

100 mL, 500 mL

For *in vitro* diagnostic use.

Zur *In-vitro*-Diagnostik.

Solo per uso diagnostico *in vitro*.

Para uso diagnóstico *in vitro*.

Pour diagnostics *in vitro*.

Para utilização em diagnóstico *in vitro*.

Για *in vitro* διαγνωστική χρήση.

Pro diagnostické použití *in vitro*.

Til *in vitro*-diagnostik.

In vitro -diagnostiikkaan.

Lietošanai *in vitro* diagnostikā.

Uitsluitend voor *in vitro* diagnostisch gebruik.

Do diagnostyki *in vitro*.

Pentru uz diagnostic *in vitro*.

För *in vitro*-diagnostik.

In vitro diagnostiliseks kasutamiseks.

In vitro diagnostikai alkalmazáshoz.

Skirta *in vitro* diagnostikai.

In vitro diagnostik kullanim için.

Na diagnostické použitie *in vitro*.

За *in vitro* диагностична употреба.

За употребу *in vitro* dijagnostici.

Ghal uz dijanjostiku *in vitro*.

За дијагностично uporabo *in vitro*.

Glossary of Symbols*:



Catalog Number



Lot Number



Sterilized using aseptic processing techniques (filtration)



Expiration:
Year - Month - Day



Caution, consult accompanying documents



Consult instructions for use



Storage Temperature
store below -10°C



Do Not Re-Sterilize



Do not use if package is damaged



Manufacturer



CE Mark



Emergo Europe - Prinsessegracht 20
2514 AP The Hague
The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

ENGLISH

INDICATION FOR USE

CHANG Amnio with Gentamicin and L-Glutamine may be used for the following applications:

1. the primary culture of amniotic fluid cells
2. growing passaged amniotic fluid cells
3. solid amnion tissue from chorionic villi sampling.

This medium has been designed for use in CO₂ incubators (cultures equilibrated with 5%-8% CO₂ atmosphere).

The final pH must be 6.65-7.44 Please see Directions for Use.

DEVICE DESCRIPTION

CHANG Amnio is a complete, ready-to-use medium for the primary culture of human amniotic fluid cells (AFC), chorionic villus sampling (CVS), and products of conception (POC) for use in karyotyping and other prenatal genetic testing. It has been optimized for both flask and *in situ* methodologies. This product contains the antibiotic Gentamicin Sulfate (50 µg/mL).

Buffers

Sodium bicarbonate

Antioxidant
Thioctic acid

Antibiotic
Gentamicin Sulfate

Amino Acid

Alanine
Arginine
Asparagine
Aspartic Acid
Cysteine
Glutamic Acid
Glutamine
Glycine
Histidine
Isoleucine
Leucine
Lysine
Methionine
Phenylalanine
Proline
Serine
Threonine
Tryptophan
Tyrosine
Valine

Salts & Ions

Sodium chloride
Sodium selenite
Calcium chloride
Choline chloride
Potassium chloride
Potassium phosphate

QUALITY ASSURANCE

STERILITY

Serum used in the production of CHANG Amnio has been tested for viral contamination per CFR Title 9 Part 113.53. It has also been screened for mycoplasma contamination. CHANG Amnio is sterilized by filtration through a 0.1 micron filter. Samples of CHANG Amnio are tested for possible bacteriological contamination following the sterility testing protocol described in the current USP Sterility test <71>.

PREPARATION FOR USE

Thaw on a sterile countertop at room temperature or by placing bottle in a 37°C water bath.

CHANG Amnio contains gentamicin (50 mg/L). Additional antibiotics may be added if desired.

ALIQUOTING CHANG Amnio

1. Thaw CHANG Amnio according to the above instructions.
2. Distribute aseptically into convenient sized aliquots and refreeze.
3. Thaw aliquots in 37°C water bath when ready to use.

DIRECTIONS FOR USE

Use of CHANG Amnio for Primary Cultures: *in situ* Methodologies

1. Centrifuge amniotic fluid at approx. 1,200 rpm for 10 minutes to concentrate the cells.
2. Aspirate supernatant from the centrifuged tube, leaving approx. 0.5 mL above cell pellet (or about 2x volume of pellet) of spun amniotic fluid. Aliquot supernatant

COMPONENTS

Proteins, Hormones, and Growth Factors
Fetal bovine serum (FBS)

Newborn bovine serum
Human transferrin
Fibroblast growth factor (FGF)
Insulin
Progesterone
Testosterone
Beta estradiol
Hydrocortisone

Nucleic acids

Deoxyadenosine
Deoxycytidine
Deoxyguanosine
Adenosine
Cytidine
Guanosine
Thymidine
Uridine

Other

Ethyl alcohol
Thyronine

Vitamins and trace elements

Ascorbic acid
Folic acid
Nicotinamide
Riboflavin
Thiamine
Pantothenic acid
Cobalamin
Pyridoxine
Biotin

Water

WFI Quality

(at least 1 mL, if possible) for alpha-fetoprotein (AFP) and acetyl cholinesterase assays, if necessary. If specimen is bloody, prepare an additional aliquot for further testing. Resuspend cell pellet in small volume of patient's own amniotic fluid. Add sufficient CHANG Amnio to the concentrated cell suspension to allow for final plating volume of 0.5 mL per coverslip (total of 4 coverslips, depending on size of cell pellet) or 2 mL per flaskette. If specimen is received from a patient in the third trimester of pregnancy, the pellet may be larger but contain less viable cells, thus requiring heavier seeding (less media than normal).

3. Incubate cultures undisturbed at 37°C 5-8% CO₂ atmosphere.
4. Flood cultures on day 2 by adding 2 mL of CHANG Amnio.
5. After 4 to 5 days, cultures should be checked for growth. Cultures should be fed once growth has been observed. Feed cultures by removing all of the culture supernatant and replacing with 2 mL fresh CHANG Amnio. It is recommended that cultures be fed every 2 days thereafter. For bloody specimens, cultures may require more frequent media changes.
6. Check cultures for growth on, or after, day 5, and harvest when sufficient colonies are observed.
7. Best results obtained when cultures are fed with CHANG Amnio the day before the harvest.

Use of CHANG Amnio for Primary Cultures: Flask Methodologies

1. Centrifuge amniotic fluid at approx. 1,200 rpm for 10 minutes to concentrate the cells.
2. Aspirate supernatant from the centrifuged tube, leaving approx. 0.5 mL above cell pellet (or about 2x volume of pellet) of spun amniotic fluid. Aliquot supernatant (at least 1 mL, if possible) for alpha-fetoprotein (AFP) and acetyl cholinesterase assays, if necessary. If specimen is bloody, prepare an additional aliquot for further testing. Resuspend cell pellet in small volume of patient's own amniotic fluid. Add 4 mL of CHANG Amnio for a total volume of 5 mL per flask. If specimen is received from a patient in the third trimester of pregnancy, the pellet may be larger but contain less viable cells, thus requiring heavier seeding (less media than normal).
3. Incubate cultures undisturbed at 37°C, 5%-8% CO₂ atmosphere.
4. Check for growth on day 5. Change medium with 2 mL of fresh CHANG Amnio and harvest if sufficient cell growth is observed.
5. Check cultures for growth and completely change medium every day thereafter until sufficient colonies are observed and are ready to harvest. For bloody specimens, cultures may require more frequent media changes.
6. Best results obtained when cultures are fed with CHANG Amnio the day before the harvest.

Use of CHANG Amnio for Growing Passaged Amniotic Fluid Cells:

To passage the cells, treat cultures with trypsin (or pronase, etc.) as you would normally do when cells are grown in conventional medium. However, protease treatment should be carefully monitored. Amniotic fluid cells grown in CHANG Amnio tend to be more sensitive to protease treatment than when grown in conventional medium. It may be necessary to modify your protocol to take this into account.

Note: The pH of the medium used to feed the cultures must be between 6.65 - 7.44 (i.e. the medium must be slightly yellowish salmon color). pH can easily be adjusted by placing the medium in a 5-8% CO₂ incubator with the cap slightly loosened for about 30 minutes.

Note: Calcium Oxalate crystals commonly form in CHANG Amnio. The presence of these crystals has not been shown to cause any detrimental effect on product performance.

Chorionic Villi Sample Preparation:

1. Label a Petri dish for each specimen received, and transfer contents of sample aseptically. Add 5 mL of the prewarmed CHANG Amnio complete media to

- the dish, and place in an incubator for 30 minutes, at 35 - 39°C, 5 - 8% CO₂ atmosphere, to allow the blood to settle.
2. Clean specimen using a dissecting microscope, initially under 1.5X magnification, then adjusting to about 3X magnification. Note: Approx. 20 - 40 mg of chorionic villi sample is required.
3. Using two pairs of sterile forceps, remove blood clots and any maternal decidua from the villus material, while working within the dish, under the dissecting microscope. Villus material is light colored, tubular and/ or lumpy with visible branches and veins.
4. Transfer clean villi to another Petri dish containing prewarmed CHANG Amnio complete media. Perform final cleaning, using forceps to grasp villi and agitate gently. Note: 5 mg is the ideal amount to use per culture. Be careful to avoid damaging the fragile villi.
5. Transfer villi and media to a 15 mL centrifuge tube, and add 4 drops of antibiotic (ie. Gentamicin Sulfate, 50 µL/mL) to the centrifuge tube. Let sit for 30 minutes.
6. Centrifuge villi at approx. 1,400 rpm for 5 minutes.

Chorionic Villi Sample Culture:

1. Aspirate supernatant from the centrifuged tube, making sure to leave 0.5 mL of media above cell pellet (or about 2X volume of pellet).
2. Gently resuspend pellet. Add 2 mL of prewarmed CHANG Amnio media to centrifuge tube.
3. Add 2 mL of Trypsin EDTA and incubate culture undisturbed at 35 - 39°C, 5 - 8% CO₂ atmosphere for 10 minutes. Remove tube from incubator, resuspend pellet and place in incubator for 10 additional minutes.
4. Remove centrifuge tube from incubator, resuspend pellet and centrifuge for 8 - 10 minutes at 1,400 rpm.
5. Aspirate supernatant from the centrifuged tube. Resuspend pellet, then add 1 mL of collagenase to tube and place in incubator for 5 minutes.
6. Remove from incubator and visually check to see if pellet is cloudy and no distinct pieces of individual villi can be seen. If pellet is not cloudy, place back in incubator for 5 more minutes.
7. Add 3 mL of prewarmed CHANG Amnio to centrifuge tube to stop the action of the collagenase.
8. Centrifuge tube for 8 - 10 minutes at 1,400 rpm.
9. Aspirate supernatant from the centrifuged tube, leaving 0.5 mL of media above cell pellet. Resuspend pellet before adding CHANG Amnio used for set up.
10. Set up optimal number of cultures using 0.5 mL of CHANG Amnio per culture for each petri dish that contains a coverslip.
11. Incubate cultures undisturbed at 35 - 39°C, 5 - 8% CO₂ atmosphere.
12. Flood cultures on day 2 by adding 1.5 mL of prewarmed CHANG Amnio.
13. At 4 days, cultures should be checked for growth. If growth is observed, remove media and add 2 mL of fresh prewarmed CHANG Amnio to each coverslip. Cultures should be fed every 2 days thereafter. For bloody specimens, cultures may require more frequent media changes.
14. Check cultures for growth on day 5, and harvest when sufficient colonies are observed. Best results obtained when cultures are fed with CHANG Amnio the day before the harvest.

STORAGE AND STABILITY

Store frozen below -10°C. Product is stable until the expiration date on the bottle label when stored frozen. Unused product can be dispensed into working aliquots and refrozen for later use, or tightly capped and stored at 2°C to 8°C for up to 30 days; it may be frozen a maximum of two times. Protect from fluorescent light.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in procedures that include the indicated application for which the device is intended. Do not use any bottle in which the sterile packaging has been compromised. Do not use CHANG Amnio beyond the expiration date indicated on the label.

FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706 • Fax: 1 949 261 6522 • www.irvinesci.com

© 2020 FUJIFILM Irvine Scientific, Inc. All rights reserved. The Irvine Scientific logo and CHANG Amnio are trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions. PN 40965 Rev.4

