FUJIFILM

IrvineScientific

Continuous Single Culture - NX Complete

(CSCM-NXC)
with Gentamicin Sulfate
with Human Serum Albumin (HSA)

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REFERENCE

Biggers, JD. and Racowsky, C. *The development of fertilized human ova to the blastocyst stage in KSOM[®] medium: is a two-step protocol necessary?* RBMOnline, 5:133-140, 2002.
Pool, TB. *Recent advances in the production of viable human embryos in vitro.* RBMOnline, 4:294-302, 2002.
Biggers, JD. *Thoughts on embryo culture conditions.* RBMOnline, 4 (suppl.1):30-38, 2001.
Lane, M., Hooper, K., and Gardner, DK. *Effect of essential amino acids on mouse embryo viability and ammonium production.* J. Asst. Reprod. Genet. 18: 519-525, 2001
Biggers, JD. and McGinnis, LK. *Evidence that glucose is not always an inhibitor of mouse preimplantation development in vitro.* Hum. Reprod 16:153-163, 2001.
Devreker, F., Van den Bergh, M., Biramane, J., Winston, RML., Englert, Y., and Hardy, K. *Effects of taurine on human embryo development in vitro.* Hum. Reprod. 14: 2350-2356, 1999.

FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705 USA
Telephone: 1 949 261 7800 • 1 800 437 5706 • Fax: 1 949 261 6522 • www.irvinesci.com
© 2023 FUJIFILM Irvine Scientific, Inc. All rights reserved. The FUJIFILM Irvine Scientific logo, Continuous Single Culture and CSCM-NXC are trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions.
PN 41070-EU Rev.02
Effective Date: 23-JUL-2023

ENGLISH

EU CAUTION: For Professional Use Only.

INDICATIONS FOR USE

Continuous Single Culture-NX Complete (CSCM-NXC) is intended for use in assisted reproductive procedures which include human gamete and embryo manipulation. These procedures include the use of CSCM-NXC as a culture medium from fertilization through day 5/6 of development in vitro.

DEVICE DESCRIPTION

CSCM-NXC is a single step culture medium that has been optimized to be used in an uninterrupted culture system, without dish change or medium renewal. CSCM-NXC contains Human Serum Albumin (HSA) for final protein concentration of 5 mg/mL and the antibiotic Gentamicin Sulfate 10 µg/mL

COMPOSITION:

Amino Acids	Salts & Ions
Alanine	Sodium Chloride
Arginine	Potassium Chloride
Asparagine	Potassium Phosphate
Aspartic Acid	Calcium Chloride
Cystine	Magnesium Sulfate
Glutamic Acid	Sodium Citrate
AlanylGlutamine	Energy Substrates
Glycine	Glucose
Histidine	Pyruvate
Isoleucine	Lactate
Leucine	Buffers
Lysine	Sodium Bicarbonate
Methionine	Antibiotic
Phenylalanine	Gentamicin Sulfate
Proline	Others
Serine	EDTA
Threonine	Hydrochloric Acid
Tryptophan	Proteins
Tyrosine	Human Serum Albumin
Valine	Water
	WFI Quality

QUALITY ASSURANCE

CSCM-NXC is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10⁻³.

Each lot of CSCM-NXC is tested for:

- Endotoxin (LAL): ≤ 0.25 EU/mL
- One-Cell MEA ≥80% expanded blastocyst at 96 hours
- Sterility by the current USP Sterility Test <71>
- Human Sperm Survival Assay ≥70% of original motility at 24 hours

All results are reported on a lot specific Certificate of Analysis which is available upon request.

BUFFER SYSTEM

CSCM-NXC uses sodium bicarbonate as a buffering system. This is specifically designed for use in a CO₂ incubator.

DIRECTIONS FOR USE

CSCM-NXC is a complete, ready-to-use medium containing Gentamicin and Human Serum Albumin (HSA). It is not necessary to add any protein before use.

EQUILIBRATION

CSCM-NXC should be pre-warmed to 37°C and equilibrated to the desired pH overnight in a 5-6% CO₂ incubator prior to use. A sufficient volume of medium is required so that oocyte recovery, insemination, embryo culture and embryo transfer dishes can be prepared.

The following are general procedures for the indications for use of CSCM-NXC.

Fertilization:

On the day before oocyte retrieval, prepare insemination

dishes with CSCM-NXC overlaid with oil and pre-equilibrate overnight to 37°C in a CO₂ incubator. Immediately upon oocyte collection and identification, place oocytes into the pre-equilibrated medium and return to the incubator for the desired period (1-4 hours) prior to insemination by conventional IVF or ICSI.

Conventional IVF:

- It is recommended to aseptically dispense 50,000-100,000/mL motile sperm per microdroplet containing 1-3 oocytes.
- Return the insemination dish to the incubator and check for normal fertilization 16-20 hours post insemination.

Intracytoplasmic Sperm Injection (ICSI):

- Following at least 1 hour post oocyte denuding (and no more than 4 hours following oocyte retrieval), remove denuded oocytes from incubator and inseminate with sperm per standard ICSI protocol for your individual laboratory.
- Immediately following insemination, place 1-3 inseminated oocytes into a fresh drop of the pre-equilibrated insemination dish, return dish to the incubator and check for normal fertilization 16-20 hours post insemination.

Embryo Culture:

On the day of fertilization (one day prior to fertilization assessment), prepare embryo culture dishes with CSCM-NXC overlaid with oil and pre-equilibrate overnight to 37°C in a 5-6% CO₂ incubator.

Following fertilization assessments with the identification of the presence of normal fertilization (two pronuclei and two polar bodies), transfer 2PN zygotes into the pre-equilibrated CSCM-NXC culture dish previously prepared. It is recommended to allow the embryos to grow in a continuous, uninterrupted culture system, without changing medium, until the desired developmental stage is reached (up to day 5/6 of development).

Embryo Transfer:

On the day prior to embryo transfer CSCM-NXC should be pre-warmed to 37°C and equilibrated to the desired pH overnight in a 5-6% CO₂ incubator prior to use. A sufficient volume of medium is required for the embryo transfer. On the day of transfer, remove the cultured embryo that has been selected for transfer to the equilibrated volume of CSCM-NXC. Using a syringe and embryo transfer catheter aseptically draw-up 5-8 µl of the equilibrated CSCM-NXC followed by the selected embryo and lastly another 5-8 µl of more equilibrated CSCM-NXC for a total of 15-20 µl. After the embryo is loaded into the catheter it is ready for immediate transfer to the patient.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

STORAGE INSTRUCTIONS AND STABILITY

Store the unopened bottles refrigerated at 2° to 8°C.

Do not freeze or expose to temperatures greater than 39°C.

Duration Following Bottle Opening:

Product should be used within four (4) weeks from opening.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any bottle of medium which shows evidence of particulate matter, or cloudiness.

To avoid problems with contamination, handle using aseptic techniques and discard any excess medium that remains in the bottle or vial after the procedure is completed.

Not for injection use.

EU: Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. It is strongly recommended that every time FUJIFILM Irvine Scientific, Inc. Reproductive Media Products culture media are administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

US: CSCM-NXC contains Human Serum Albumin (HSA). Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are non-infectious. Handle all human source material as if it were capable of transmitting infection, using universal pre-cautions. Donors of the source material have also been screened for CJD.

CONTRAINDICATION

CSCM-NXC contains the antibiotic Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

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