ISolate Concentrate

Catalog No. 99306

100 mL

REFERENCES


Adaniya GK, Jackson KV: *Comparison of Percoll® and ISolate in the Preparation of Semen for use in Assisted Reproductive Technologies*. 53rd Annual Meeting of the American Society for Reproductive Medicine, Cincinnati, Ohio; 0-113, 1997.

Tarchala SM, Fahy MM, Hauserman HM, Kallmann K, Volentine K, Radwanska E, Binor Z, Molo MW, Rawlins RG: *ISolate As a New Method of Sperm Separation*. 53rd Annual Meeting of the American Society for Reproductive Medicine, Cincinnati, Ohio; P-116, 1998.

Miller, KF, Fry KL, Arciaga RL, Falcone T: *Semen preparation for in-vitro fertilization using ISolate results in sperm recovery, fertilization and pregnancy outcomes that are indistinguishable from those obtained with Percoll®*. 14th Embryology, Goteborg, Sweden; P-116, 1998.

Tarchala SM, Volentine KK, Rawlins RG: *A comparison of sperm processing using ISolate and Pure Sperm™*. 14th Annual Meeting of the European Society of Human Reproduction and Embryology, Goteborg, Sweden; R-041, 1998.

Glossary of Symbols*:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature 2-8°C
	Do not re-sterilize
	Do not re-use
	Do not use if package is damaged
	Manufacturer
	U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	CE Mark
	Emergo Europe - Prinsessegracht 20 2514 AP The Hague The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices - Symbols to be used with medical device labels, labeling.

ENGLISH

EU Caution: For Professional Use Only.

INDICATION FOR USE:

ISolate is intended for assisted reproductive procedures that involve the manipulation of human sperm. ISolate is intended for the separation of the motile fraction of sperm from seminal fluid.

DEVICE DESCRIPTION

ISolate is a density gradient medium designed to separate the motile fraction of sperm from seminal fluid. As a two layer gradient system, it effectively reduces cellular contaminants such as dead sperm, white blood cells and miscellaneous debris. The resulting sample contains predominantly motile sperm.

COMPOSITION

Salts and Ions	Energy Substrate
Sodium Chloride	Glucose
Potassium Chloride	Sodium Pyruvate
Magnesium Sulfate	Sodium Lactate
Potassium Phosphate	Other
Calcium Chloride	Colloidal Suspension of Silica Particles
Buffer	
Sodium Bicarbonate	
HEPES	

QUALITY ASSURANCE

ISolate is a membrane filtered, aseptically processed colloidal suspension of silica particles stabilized with covalently bound hydrophilic silane in HEPES-buffered HTF. The sterility assurance level (SAL) is 10⁻³. ISolate is tested and found negative for pyrogens by rabbit pyrogen test.

DIRECTIONS FOR USE

Catalog #99306 - CONCENTRATE

Colloidal silica concentrate for further dilution into density gradients; formulated in an isotonic buffer and approximates the density of undiluted Percoll®. ISolate CONCENTRATE was evaluated according to this package insert at 50% and 90% gradients. Each laboratory may choose to make gradients at different concentrations which should be evaluated for effectiveness.

A. Materials Suggested:

- ISolate Concentrate
- FUJIFILM Irvine Scientific, Inc. Catalog #99306 Modified Human Tubal Fluid
- FUJIFILM Irvine Scientific, Inc. Catalog #90126 Sperm Washing Medium
- FUJIFILM Irvine Scientific, Inc. Catalog #9983 Centrifuge
- FUJIFILM Irvine Scientific, Inc. Catalog #IS-300 Sterile, disposable conical centrifuge tubes
- Pasteur Pipettes, sterile
- Incubator, 37°C

B. Gradient Preparation:

Two solutions are required to perform the gradient sperm separation method: An upper layer and a lower layer.

1. Upper Layer solution - 50%:
 - a) Take CONCENTRATE (Catalog #99306) and pipette 4.5 mL into a sterile centrifuge tube
 - b) Add 5.5 mL of Modified Human Tubal Fluid (Catalog #90126)
 - c) Cap tube and invert gently 3-4 times to allow for complete mixing
 - d) Label tube as "UPPER LAYER"
2. Lower Layer solution - 90%:
 - a) Take CONCENTRATE (Catalog #99306) and pipette 8.1 mL into a sterile centrifuge tube
 - b) Add 1.9 mL of Modified Human Tubal Fluid (Catalog #90126)
 - c) Cap tube and invert gently 3-4 times to allow for complete mixing
 - d) Label tube as "LOWER LAYER"
3. Continue to GENERAL PROCEDURE section of product insert.

Note: The terms 50% and 90% do not represent an actual concentration. The 50% and 90% are relative concentrations based on the previous nomenclature that defined a 1:9 dilution of Percoll® as equal to a 100% working, isotonic solution.

GENERAL PROCEDURE

The following is a general procedure for a two-step gradient separation of motile sperm from semen. The volume and concentration can be modified to meet each laboratory's preference (i.e. mini-gradient, cryopreserved semen, three-layer gradient).

1. Bring all media components to room temperature or 37°C.
2. Using a sterile pipette, transfer 1.5 - 2.0 mL of the "LOWER LAYER" into a sterile, disposable, conical centrifuge tube.
3. Using a new sterile pipette, transfer an equal volume of "UPPER LAYER" on top of the "LOWER LAYER". This is done by contacting the surface of the "LOWER LAYER" at the side of the tube with the tip of the pipette. Carefully dispense the "UPPER LAYER" by spiraling the pipette tip around the circumference of the tube in an upward motion as the level of the "UPPER LAYER" rises.
4. Gently place 1.5 - 2.0 mL of liquefied semen onto the "UPPER LAYER" using a new sterile pipette.
5. Centrifuge for 10-20 minutes at approximately 200 - 300 x g.* Carefully expose the pellet by either aspirating off the "UPPER" and "LOWER LAYERS", or directly removing the pellet and transferring to a new sterile centrifuge tube.

Note: The quality of the sperm sample should be evaluated and taken into consideration when determining the proper centrifugation speed and time. These should be adjusted according to the individual specimen quality for optimization of the procedure.

6. Using a new sterile pipette, add 2.0 - 3.0 mL of appropriate washing medium such as Sperm Washing Medium (Catalog #9983) or Modified HTF (Catalog #90126) with protein supplementation. Resuspend the isolated pellet. Centrifuge (~200xg) for 8 - 10 minutes and remove the supernatant. Repeat this step for a second wash. Discard the supernatant and resuspend the pellet using a suitable volume of appropriate medium. The sample is now ready for analysis.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

STORAGE INSTRUCTIONS AND STABILITY

Store the unopened containers refrigerated at 2° to 8°C.

Warm to ambient or incubator (37°C) temperature prior to use.

Do not freeze or expose to high temperatures.

ISolate is stable until the expiration date shown on the box and bottle labels when stored as directed.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

ISolate will appear opaque. This is normal for this product. Do not use any vial of medium which shows evidence of particulate matter or contamination (nonuniform cloudiness).

ISolate should be tightly capped when used in a CO₂ incubator to avoid pH changes.

Information on known characteristics and technical factors that could pose a risk if the product were to be re-used have not been identified therefore the product is not to be used following the initial use of the container.

In case of damage, do not use. Please discard or return to Manufacturer for replacement.

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© 2012 FUJIFILM Irvine Scientific, Inc. Kaikki oikeudet pidätetään.