

CTGrade GMP rh IL-4

Catalog #	Product	Size
500-05	CTGrade GMP rh IL-4	50µg, 100µg, 1mg lyophilized

Intended Use

This product is for research or further manufacturing use only. Not for injection or diagnostic procedures. The safety and efficacy of this product in diagnostic or other clinical procedures has not been established.

Product Description

This product is produced from *E. coli* and is manufactured in a facility that does not use or process beta-lactam containing materials. No animal- or human-derived materials were used during manufacturing or as ingredients. As such, the risk for BSE/TSE contamination can be considered negligible. This product is manufactured, tested, and released in an ISO 9001:2015 certified facility and follows cGMP practices. USP chapter <1043> for ancillary materials has been considered in the manufacture of this product.

Synonyms: B-cell stimulatory factor-1, BSF-1, BCDF, BCGF

NIH Accession Number: P05112

Background: Interleukin 4 (IL-4) is an immunomodulatory cytokine that stimulates B cell proliferation and activates eosinophils, basophils, and mast cells (1, 2). IL-4 is also involved in collagen production by fibroblasts and can induce vascular cell adhesion molecule (VCAM)-1 expression on endothelial cells [3]. IL-4 acts a key role in promotion of Th2 differentiation (4, 5, 6, 7) from naïve helper T cells. Functionally, IL-4 is best known for defining the Th2 phenotype of lymphocytes and for regulating cell proliferation, apoptosis, and expression of numerous genes in various cell types, including lymphocytes, macrophages, and fibroblasts, as well as epithelial and endothelial cells (8).

Specifications

Formulation:	CTGrade GMP rh IL-4 is lyophilized from a 0.2 µm filtered solution containing 20 mM Sodium Phosphate, pH 7.5
Protein Purity:	≥98% determined by reducing and non-reducing SDS-PAGE analysis.
Endotoxin:	<0.05 EU/µg using USP <85>/ EP 2.6.14
Bioactivity:	ED50 is determined by the dose-dependent TF-1 cell proliferation cells. The ED50 is typically less than 0.25 ng/mL. The international units of CTGrade GMP rh IL-4 is approximately 1.4 x 10 ⁴ IU/µg, which is calibrated against recombinant Human Interleukin 4 WHO International Standard (NIBSC code 88/656).
Quality:	Carrier-free and no animal- or human-derived materials were used during manufacturing.

Quality Assurance

All quality control test results are reported on a lot specific Certificate of Analysis, which is available at www.irvinesci.com or upon request.

Shipping

This product is shipped at ambient temperature. Immediately upon receipt, store at the recommended temperature below.

Storage Instructions and Stability

Upon receipt, store the lyophilized protein at -20°C in a manual defrost freezer. Unopened vials are stable for 36 months from the date of manufacture. Reconstituted material should be apportioned in working volumes and stored at or below -20°C in manual defrost freezer.

For short term storage reconstituted material is stable for 4-6 weeks when stored at 2-8°C. Stability can be increased by adding at least 0.1% carrier protein.

Precautions

For *ex vivo* use only. Not for injection or diagnostic procedures. The safety and efficacy of this product in diagnostic or other clinical uses has not been established. Please refer to the Safety Data Sheet for information regarding hazards and safe handling practices.

Directions for Use

1. Reconstitution

Allow the vial and sterile water (e.g. FUJIFILM Irvine Scientific, Inc. P/N 9309 Water for Injection) to equilibrate to room temperature. Draw up desired volume of reconstitution buffer. Aseptically puncturing through rubber stopper with sterile needle, inject the buffer to achieve the desired concentration (0.1-0.5 mg/mL). Swirl the vial gently, **do not vortex**. Allow protein to rehydrate for 10-15 minutes at room temperature with occasional gentle mixing.

2. Optimum Concentration

The optimum concentration varies depending on cell type and culture conditions. Working concentration should be determined for each specific application.

Related Products

Catalog #	Product	Size
91146	PRIME-XV Dendritic Cell Maturation CDM	500 mL
91211	PRIME-XV Hematopoietic Cell Basal XSFM	500 mL
9160	RPMI Medium 1640 w/o L-Glutamine - Liquid	100mL, 500mL, 1L
9161	RPMI Medium 1640 with L-Glutamine	100mL, 500mL, 1L
9309	Water for Injection	1L

References

1. Paul, W.E. Interleukin-4: A prototypic immunoregulatory lymphokine. *Blood* 1991, 77, 1859–1870.
2. Yanagihara, Y.; Ikizawa, K.; Kajiwara, K.; Koshio, T.; Basaki, Y.; Akiyama, K. Functional significance of IL-4 receptor on B cells in IL-4-induced human IgE production. *J. Allergy Clin. Immunol.* 1995, 96 Pt 2, 1145–1151
3. Schnyder, B.; Lugli, S.; Feng, N.; Etter, H.; Lutz, R.A.; Ryffel, B.; Sugamura, K.; Wunderli-Allenspach, H.; Moser, R. Interleukin-4 (IL-4) and IL-13 bind to a shared heterodimeric complex on endothelial cells mediating vascular cell adhesion molecule-1 induction in the absence of the common gamma chain. *Blood* 1996, 87, 4286–4295.
4. Barner, M.; Mohrs, M.; Brombacher, F.; Kopf, M. Differences between IL-4R α -deficient and IL-4-deficient mice reveal a role for IL-13 in the regulation of Th2 responses. *Curr. Biol.* 1998, 8, 669–672.
5. Haas, H.; Falcone, F.H.; Holland, M.J.; Schramm, G.; Haisch, K.; Gibbs, B.F.; Bufe, A.; Schlaak, M. Early interleukin-4: Its role in the switch towards a Th2 response and IgE-mediated allergy. *Int. Arch. Allergy Immunol.* 1999, 119, 86–94.
6. Mosmann, T.R.; Coffman, R.L. TH1 and TH2 Cells: Different Patterns of Lymphokine Secretion Lead to Different Functional Properties. *Annu. Rev. Immunol.* 1989, 7, 145–173.
7. Paul, W.E. History of interleukin-4. *Cytokine* 2015, 75, 3–7.
8. Kelly-Welch A. E., Hanson E. M., Boothby M. R., Keegan A. D. (2003) Interleukin-4 and interleukin-13 signaling connections maps. *Science* 300, 1527–1528

Technical Support

CONTACT US

For more information or assistance contact Customer Service at:

- Email: fisitmrequest@fujifilm.com
- Direct line: +1 800 577 6097

WEBSITE RESOURCES

Visit the website at www.irvinesci.com for technical resources and information including:

- Safety Data Sheets (SDS)
- Certificate of Analysis (CoA) (when available)
- FAQs
- Product literature

FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705-5588 USA

Telephone: 1 949 261 7800 • 1 800 437 5706

Fax: 1 949 261 6522 • www.irvinesci.com

© 2022 FUJIFILM Irvine Scientific, Inc. All rights reserved. PRIME-XV, FUJIFILM Irvine Scientific and its logo are registered trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions. Rev.00

 		Category: CofA, Label, PI Title: CTGrade rh IL-4 PI	
Version 03	State Effective	Effective Date 01-APR-2024	Document ID 883010
Department Category QA	Protein Name N/A		Training Type N/A

Printed by anindita.rudra@fujifilm.com from app.zenqms.com on 01-May-2024 at 6:08:45 PM UTC • Page 6 of 6

REVISION HISTORY

Version 01 Effective on 13-Jun-2023

Initial

Version 02 Effective on 08-Aug-2023

Updated Related products, added approximate IU/ug in activity section

Version 03 Effective on 01-Apr-2024

Added GMP in name, add BSE/TSE statement

DOCUMENT ELECTRONIC SIGNATURES

DOCUMENT APPROVAL WORKFLOW

Author Approval

Iane Charles
Quality Control Supervisor
ianec@shenandoah-bt.com

I am the author of this document.
Signed 5:12:54 PM UTC 28-Feb-2024

Required Workflow Steps for this Category

Iane Charles
Quality Control Supervisor
ianec@shenandoah-bt.com

Technical Reviewer/Approval
I have reviewed and approve this document.
Signed 5:41:44 PM UTC 28-Feb-2024

Amanda Cinquin
amanda.cinquin@fujifilm.com

Compliance Reviewer/Approval
I have reviewed and approve this document.
Signed 8:41:30 PM UTC 28-Feb-2024

Anindita Rudra
Product Manager
anindita.rudra@fujifilm.com

Marketing Reviewer/Approval
I have reviewed and approve this document.
Signed 8:58:28 PM UTC 28-Feb-2024

Christopher Unger
QA Supervisor
christopher.unger@fujifilm.com

QA Approval
I have reviewed and approve this document.
Signed 8:38:21 PM UTC 29-Feb-2024