



BILL OF MATERIALS

PART#	QTY	DESCRIPTION
-	100	20 LITER BAG
3C-1021	100	20 LITER DIE
PL-01026	-	ASI THERM FILM
PK-01098	-	ASI OUTER BAGGER FILM
CN-01249	100	6.325" DIA. HANGING BAR 12" LING
TB-01131	100	BOAT-FITMENT WITH THREE 3/8" BARS
CL-01011	100	3/8" ID. X 5/8" OD. C-FLEX CLEAR TUBING
CL-01010	200	SMALL PINCH CLAMP
CN-01057	100	LARGE PINCH CLAMP
CP-01055	100	MALE PVC TUBE WITH 3/8" BARS
TB-01023	100	FEMALE PVC CAP
CH-01061	0.50 FT.	1/4" ID. X 3/8" OD. POKES TUBING
CH-01025	100	FEMALE LUER WITH 1/4" BAR
CP-01025	100	SWABBLE VALVE
CP-01025	100	MALE COVER CAP

REV	DATE	DESCRIPTION
A	03-07-06	NEW PRODUCT.
B	07-09-10	CHANGED PACKAGING NOTES.
C	10-17-12	CHANGED CUSTOMER NAME FROM MICROBIX TO IRVINE SCIENTIFIC. DRAWING FORMAT UPDATE.



ASI

163 Research Lane
Millersburg, PA 17061-9402
717-692-2104 - Phone
717-692-2197 - Fax
www.advancedscientifics.com

CONFIDENTIAL: THIS IS PROPRIETARY INFORMATION AND IS NOT TO BE SHARED WITH ANY PARTY WITHOUT WRITTEN CONSENT OF ADVANCED SCIENTIFICS INC.

CUSTOMER
IRVINE SCIENTIFIC

TITLE
15 X 20 LITER MANIFOLD

SHEET
1 OF 2

SCALE
NTS

DRAWING NO.
B100452-1

REVISION
C

APPROVALS
ENGINEERING
AES

DATE
10-17-12