



**BILL OF MATERIALS**

PART#	QTY	DESCRIPTION
-	100	20 LITER BAG
BC-10261	100	20 LITER BIE
PL-01026	-	ASI TINNER FILM
PL-01029	-	ASI OUTER BARRIER FILM
PK-01098	100	6325' DIA. HANGING BAR 12' LING
CN-01249	100	BOAT-FITMENT WITH THREE 3/8" BARS
TB-01131	667 FT.	3/8" ID. X 5/8" OD. C-FLEX CLEAR TUBING
CL-01011	100	SMALL PINCH CLAMP
CL-01010	200	LARGE PINCH CLAMP
CN-01057	100	MALE PVC TUBE WITH 3/8" BARS
CP-01068	100	FEMALE PVC CAP
TB-01023	0.50 FT.	1/4" ID. X 3/8" OD. POKES TUBING
CN-01061	100	FEMALE LUER WITH 1/4" BARB
CN-01225	100	SWIVABLE VALVE
CP-01055	100	MALE COVER CAP

REV	DATE	DESCRIPTION
A	03-07-06	NEW PRODUCT.
B	07-09-10	CHANGED PACKAGING NOTES.
C	10-17-12	CHANGED CUSTOMER NAME FROM MICROBIX TO IRVINE SCIENTIFIC. DRAWING FORMAT UPDATE.



**ASI**

163 Research Lane  
 Millersburg, PA 17061-9402  
 717-692-2104 - Phone  
 717-692-2197 - Fax  
[www.advancedscientifics.com](http://www.advancedscientifics.com)  
 THIS IS PROPRIETARY INFORMATION AND IS NOT TO BE SHARED WITH ANY PARTY WITHOUT WRITTEN CONSENT OF ADVANCED SCIENTIFICS INC.

**CUSTOMER**  
 IRVINE SCIENTIFIC

**TITLE**  
 15 X 20 LITER MANIFOLD

**SHEET**  
 1 OF 2

**SCALE**  
 NTS

**DRAWING NO.**  
 B100452-1

**REVISION**  
 C

**APPROVALS**

ENGINEERING: AES  
 CUSTOMER: ZUSTWICK

**DATE**  
 10-17-12