



Irvine**Scientific**[®]

Anti-Clumping Supplement



Irvine**Scientific**[®]

2511 Daimler Street, Santa Ana, California 92705-5588 USA

Telephone: 1 949 261 7800 • 1 800 437 5706

Fax: 1 949 261 6522 • www.irvinesci.com

PN 41062 Rev. 0

Anti-Clumping Supplement

Catalog # 91150

PRODUCT DESCRIPTION

Anti-Clumping Supplement is designed to reduce cell aggregation for various cell types in suspension culture. This supplement is animal component and enzyme-free. This supplement is not intended to dissociate adherent cells.

INTENDED USE

For research or further manufacturing use only.

QUALITY ASSURANCE

All quality control test results are reported on a lot specific Certificate of Analysis that is available upon request.

STORAGE INSTRUCTIONS AND STABILITY

Handle using aseptic techniques to avoid contamination. Store at room temperature or 2-8°C. Do not use any bottle of medium that shows evidence of particulate matter or cloudiness.

INSTRUCTIONS FOR USE

Anti-Clumping Supplement is designed for use at a 1:200 – 1:1000 dilution (1 – 5mL/L) range. Typically used at 2mL/L. Supplement additional amount if necessary.

Note: Anti-Clumping Supplement will completely inhibit transfection via cationic liposomes or cationic polymers such as polyethylenimine (PEI). Cells should be spun down and re-suspended in media without Anti-Clumping Supplement before transfection.

2mL/L of Anti-Clumping Supplement added to HEK293 cultures

Without supplement



With supplement

