

## PRODUCT SHEET

---

### IDENTIFICATION OF PRODUCT

<b>NAME OF PRODUCT</b>	ASPIGLAIRE/GLAIR SUCKER
<b>REF. NUMBER(S)</b>	007447 - 007448
<b>PACKAGING</b>	PACKAGED PER 1, BOX OF 20. PACKAGED PER 5 UNITS, BOX OF 40
<b>APPLICATION</b>	HÜHNER TEST – MUCUS VISCOSITY EVALUATIONS
<b>COUNTRY OF ORIGIN</b>	FRANCE
<b>MANUFACTURER</b>	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
<b>MEDICAL DEVICE</b>	CE CLASS I, STERILE
<b>STERILIZED</b>	YES (STERILIZATION BY IRRADIATION)
<b>SHELF LIFE</b>	3 YEARS

---

### COMPOSITION / INFORMATION ON THE COMPONENTS

Material = PVC (Polyvinyl Chloride)  
Blue color index, biocompatible, non toxic

---

### SPECIFICATIONS

Length = +/- 280 mm  
Diameter = +/- 3 mm

---

### BIOCOMPATIBILITY

MEA (Mouse Embryo Assay) and LAL (endotoxin) tested

---

### MANIPULATION AND STORAGE / STABILITY

Do not use if the product or the packaging is damaged as sterility and integrity of the product can no longer be guaranteed  
Storage cautions: Preserve away from light, heat and damp

Cryo Bio System  
Groupe IMV Technologies  
27, rue Taitbout  
75009 PARIS - FRANCE  
Tél. : +33 (0)149 240 505  
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY