

PRODUCT SHEET

IDENTIFICATION OF PRODUCT

NAME OF PRODUCT	CBS™ ASPIRATION TUBING
REF. NUMBER	007451
PACKAGING	BLISTER 5X10 UNITS
APPLICATION	STERILE ASPIRATION NOZZLE FOR CBS™ HIGH SECURITY SPERM STRAW ON PACE OR MAPI
COUNTRY OF ORIGIN	FRANCE
MANUFACTURER	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
MEDICAL DEVICE	CE CLASS 2A, STERILE
STERILIZED	YES (STERILIZATION BY IRRADIATION)
SHELF LIFE	3 YEARS

COMPOSITION/INFORMATION ON THE COMPONENTS

PVC (PolyVinyl Chloride) biocompatible, non toxic

SPECIFICATIONS

Length = +/- 160 mm

Diameter = +/- 2 mm

BIOCOMPATIBILITY

NA

MANIPULATION AND STORAGE/STABILITY

Do not use if the product or the packaging is damaged as sterility and integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp

Cryo Bio System
Groupe IMV Technologies
27 rue Taitbout
75009 PARIS - FRANCE
Tél. : +33 (0) 149 240 505
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY