

PRODUCT SHEET

IDENTIFICATION OF PRODUCT

| | |
|--------------------------|-------------------------------------------------------------------------|
| NAME OF PRODUCT | FILLING NOZZLE, YELLOW |
| REF. NUMBER | 008656 |
| PACKAGING | BAG OF 10 UNITS |
| APPLICATION | YELLOW STERILE FILLING NOZZLE FOR CBS™ HIGH SECURITY SPERM STRAW |
| COUNTRY OF ORIGIN | FRANCE |
| MANUFACTURER | CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES |
| MEDICAL DEVICE | CE CLASS 2A, STERILE |
| STERILIZED | YES (STERILIZATION BY IRRADIATION) |
| SHELF LIFE | 3 YEARS |

COMPOSITION/INFORMATION ON THE COMPONENTS

Nozzle = PETG (PolyEthylene Terephthalate), + EVA (Ethylene Vinyl Acetate)
biocompatible, non toxic

SPECIFICATIONS

Length = +/- 50 mm
Diameter = +/- 2 / 7 mm

BIOCOMPATIBILITY

MANIPULATION AND STORAGE/STABILITY

Do not use if the product or the packaging is damaged as sterility and integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp

Cryo Bio System
Groupe IMV Technologies
27 rue Taitbout
75009 PARIS - FRANCE
Tél. : +33 (0) 149 240 505
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY