

PRODUCT SHEET

IDENTIFICATION OF PRODUCT

NAME OF PRODUCT	PETG SPERM STRAW
REF. NUMBER	010261 (0.25 ml) AND 014550 (0.50 ml)
PACKAGING	BAG OF 100 UNITS
APPLICATION	PACKAGING AND CRYOGENIC STORAGE OF SPERM
COUNTRY OF ORIGIN	FRANCE
MANUFACTURER	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
MEDICAL DEVICE	CE CLASS 2A, STERILE
STERILIZED	YES (STERILIZATION BY IRRADIATION)
SHELF LIFE	3 YEARS

COMPOSITION/INFORMATION ON THE COMPONENTS

Tube = PETG (PolyeEthylene Terephtalate), biocompatible, non toxic
Plug = Polyester viscosse + polymerizing powder

SPECIFICATIONS

Length = +/- 133 mm

Diameter (outer) = +/- 1.9 mm for 0.25 ml straw and +/- 2.5 mm for 0.50 ml straw

BIOCOMPATIBILITY

Validated on human sperm

MANIPULATION AND STORAGE/STABILITY

Precautions for use: The powder of the plug should be wet between 1/3 and 1/2

Do not use if the product or the packaging is damaged as sterility and integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp



Cryo Bio System
Groupe IMV Technologies
27, rue Taitbout
75009 PARIS - FRANCE
Tél. +33 (0) 149 240 505
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY