

PRODUCT SHEET

IDENTIFICATION OF PRODUCT

NAME OF PRODUCT	EMPTYING DEVICE FOR CBS™ HIGH SECURITY STRAWS
REF. NUMBER	010284
PACKAGING	BAG 20 UNITS
APPLICATION	STERILE EMPTYING DEVICE COMPOSED OF PROTECTIVE SHEETH AND PISTON FOR CBS™ HIGH SECURITY SPERM STRAW
COUNTRY OF ORIGIN	FRANCE
MANUFACTURER	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
MEDICAL DEVICE	CE CLASS 2A, STERILE
STERILIZED	YES (STERILIZATION BY IRRADIATION)
SHELF LIFE	3 YEARS

COMPOSITION/INFORMATION ON THE COMPONENTS

LDPE (Low Density PolyEthylene) + PVC (PolyVinyl Chloride) biocompatible, non toxic

SPECIFICATIONS

Length = +/- 130 / 125 mm

Diameter = +/- 2.2 / 5 mm

BIOCOMPATIBILITY

Validated on human sperm

MANIPULATION AND STORAGE/STABILITY

Do not use if the product or the packaging is damaged as sterility and integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp

Cryo Bio System
Groupe IMV Technologies
27, rue Taitbout
75009 PARIS - FRANCE
Tél. +33 (0) 149 240 505
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY