

PRODUCT SHEET

IDENTIFICATION OF PRODUCT

NAME OF PRODUCT	INTRACERVICAL INSEMINATION CATHETER
REF. NUMBER	015821
PACKAGING	INDIVIDUALLY PACKAGED, BOX OF 20 UNITS
APPLICATION	ARTIFICIAL INTRACERVICAL INSEMINATION WITH FROZEN SPERM IN CBS™ HIGH SECURITY STRAWS
COUNTRY OF ORIGIN	FRANCE
MANUFACTURER	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
MEDICAL DEVICE	CE CLASS I, STERILE
STERILIZED	YES (STERILIZATION BY IRRADIATION)
SHELF LIFE	3 YEARS

COMPOSITION/INFORMATION ON THE COMPONENTS

PVC (Polyvinyl Chloride)
PETG (PolyEthylene Terephthalate)
LDPE (Low Density PolyEthylene)
White color index, biocompatible, non toxic

SPECIFICATIONS

Length = +/- 280 mm
Diameter = +/- 4 mm


BIOCOMPATIBILITY

MEA (Mouse Embryo Assay) and LAL (endotoxin) tested

MANIPULATION AND STORAGE/STABILITY

Do not use if the product or the packaging is damaged as sterility and integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp


Cryo Bio System
Groupe IMV Technologies
27, rue Taitbout
75009 PARIS - FRANCE
Tél. +33 (0) 149 240 505
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY