

PRODUCT SHEET

IDENTIFICATION OF PRODUCT

NAME OF PRODUCT	CBS™ INJECTION TUBING LONG
REF. NUMBER	018620
PACKAGING	BOX OF 25 UNITS
APPLICATION	STERILE LONG BLUE INJECTION NOZZLE FOR CBS™ HIGH SECURITY SPERM STRAW ON PACE OR MAPI
COUNTRY OF ORIGIN	FRANCE
MANUFACTURER	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
MEDICAL DEVICE	CE CLASS 2A, STERILE
STERILIZED	YES (STERILIZATION BY IRRADIATION)
SHELF LIFE	3 YEARS

COMPOSITION/INFORMATION ON THE COMPONENTS

PVC (PolyVinyl Chloride) biocompatible, non toxic

SPECIFICATIONS

Length = +/- 220 mm

Diameter = +/- 2 mm

BIOCOMPATIBILITY

Validated on human sperm

MANIPULATION AND STORAGE/STABILITY

Do not use if the product or the packaging is damaged as sterility and integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp

Cryo Bio System
Groupe IMV Technologies
27, rue Taitbout
75009 PARIS - FRANCE
Tél. : +33 (0) 149 240 505
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY