

PRODUCT SHEET

IDENTIFICATION OF PRODUCT

NAME OF PRODUCT	CBS™ HIGH SECURITY TISSUE 2 ML STRAW
REF. NUMBER	018960
PACKAGING	INDIVIDUALLY PACKAGED, BAG WITH 5 BLISTERS OF 4 UNITS
APPLICATION	SEALABLE CONTAINER FOR CRYOGENIC STORAGE
COUNTRY OF ORIGIN	FRANCE
MANUFACTURER	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
MEDICAL DEVICE	CE CLASS 2A, STERILE
STERILIZED	YES (STERILIZATION BY IRRADIATION)
SHELF LIFE	3 YEARS

COMPOSITION/INFORMATION ON THE COMPONENTS

Tube = Ionomeric resin, biocompatible, non toxic

SPECIFICATIONS

Length = +/- 133 mm

Diameter (outer) = +/- 6 mm


BIOCOMPATIBILITY

MEA (Mouse Embryo Assay) and LAL (endotoxin) tested

MANIPULATION AND STORAGE/STABILITY

Do not use if the product or the packaging is damaged as sterility and integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp


Cryo Bio System
Groupe IMV Technologies
27 rue Taitbout
75009 PARIS - FRANCE
Tél. : +33 (0) 149 240 505
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY