

PRODUCT SHEET

IDENTIFICATION OF PRODUCT

NAME OF PRODUCT	DAISY GOBLETS
REF. NUMBER	015144, 015146 – 015152, 015154 – 15159, 018141, 018142
PACKAGING	INDIVIDUALLY PACKED OR BAG OF 5
APPLICATION	STORAGE GOBLET FOR ALL TYPES OF STRAWS
COUNTRY OF ORIGIN	FRANCE
MANUFACTURER	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
MEDICAL DEVICE	NO, ACCESSORY
STERILIZED	No
SHELF LIFE	ND

COMPOSITION/INFORMATION ON THE COMPONENTS

LDPE (Low Density PolyEthylene)

SPECIFICATIONS

Height: 140 mm

Diameter: 68 mm


BIOCOMPATIBILITY

ND

MANIPULATION AND STORAGE/STABILITY

Do not use if the product or the packaging is damaged as integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp

 Cryo Bio System
Groupe IMV Technologies
27, rue Taitbout
75009 PARIS - FRANCE
Tél. : +33 (0) 149 240 505
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY