

PRODUCT SHEET

IDENTIFICATION OF PRODUCT

NAME OF PRODUCT	HSV HIGH SECURITY VITRIFICATION KIT
REF. NUMBER	019916 (RED), 019917 (WHITE), 019918 (YELLOW), 019919 (PURPLE), 019920 (BLUE), 019921 (GREEN)
PACKAGING	INDIVIDUALLY PACKAGED, POUCH WITH 1 BLISTER OF 4 UNITS
APPLICATION	VITRIFICATION OF OOCYTES AND EMBRYOS
COUNTRY OF ORIGIN	FRANCE
MANUFACTURER	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
MEDICAL DEVICE	CE CLASS 2A, STERILE
STERILIZED	YES (STERILIZATION BY IRRADIATION)
SHELF LIFE	3 YEARS

COMPOSITION/INFORMATION ON THE COMPONENTS

Tube = Ionomeric resin, biocompatible, non toxic

Spatula = Styrene Butadiene Copolymer + PETG, biocompatible, non toxic

Blue color index for the introduction device only, biocompatible, non toxic

SPECIFICATIONS

Length = +/- 133 mm

Diameter (outer) = +/- 1.9 mm


BIOCOMPATIBILITY

MEA (Mouse Embryo Assay) and LAL (endotoxin) tested

MANIPULATION AND STORAGE/STABILITY

Do not use if the product or the packaging is damaged as sterility and integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp


Cryo Bio System
Groupe IMV Technologies
27, rue Taitbout
75009 PARIS - FRANCE
Tél. : +33 (0) 149 240 505
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485;2003 CERTIFIED COMPANY