

PRODUCT SHEET

IDENTIFICATION OF PRODUCT

NAME OF PRODUCT	ID AND OBTURATING ROD
REF. NUMBER	007437 (yellow), 007438 (purple), 007440 (red), 007442 (orange), 014484 (pink), 014485 (blue), 014486 (grey), 014487 (green), 014488 (royal blue) and 014489 (white)
PACKAGING	BLISTER OF 4 X 5 UNITS
APPLICATION	OBTURATION OF PETG EMBRYO STRAWS
COUNTRY OF ORIGIN	FRANCE
MANUFACTURER	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
MEDICAL DEVICE	CE CLASS 2A, STERILE
STERILIZED	YES (STERILIZATION BY IRRADIATION)
SHELF LIFE	3 YEARS

COMPOSITION/INFORMATION ON THE COMPONENTS
PS (Polystyrene + color index), biocompatible, non toxic

SPECIFICATIONS
Length = 42 mm
Diameter (outer) = +/- 2.5 mm


BIOCOMPATIBILITY
Validated on animal and human semen

MANIPULATION AND STORAGE/STABILITY

Precaution for use: The storage of the straw in LN2 should be cotton plug at the bottom, ID rod at the top

Do not use if the product or the packaging is damaged as sterility and integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp


Cryo Bio System
Groupe IMV Technologies
27, rue Taitbout
75009 PARIS - FRANCE
Tél. : +33 (0)149 240 505
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY