

PRODUCT SHEET

IDENTIFICATION OF PRODUCT

NAME OF PRODUCT	ID ROD
REF. NUMBER	010266 – 010271, 019021 - 019026
PACKAGING	BAGS OF 10
APPLICATION	COLORED WEIGHTED IDENTIFICATION RODS FOR CBS™ EMBRYO AND SPERM STRAWS
COUNTRY OF ORIGIN	FRANCE
MANUFACTURER	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
MEDICAL DEVICE	NO, ACCESSORY
STERILIZED	NO
SHELF LIFE	ND

COMPOSITION/INFORMATION ON THE COMPONENTS

LDPE (Low Density PolyEthylene) + Stainless steel

SPECIFICATIONS

Length: +/- 30 / 40 mm

BIOCOMPATIBILITY

NA

MANIPULATION AND STORAGE/STABILITY

Do not use if the product or the packaging is damaged as integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp

Cryo Bio System
Groupe IMV Technologies
27, rue Taitbout
75009 PARIS, FRANCE
tel. +33 (0) 149 240 505
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY