

## PRODUCT SHEET

### IDENTIFICATION OF PRODUCT

<b>NAME OF PRODUCT</b>	<b>CBST™ HIGH SECURITY SPERM STRAW</b>
<b>REF. NUMBER</b>	010287 – 010288 – 014650 – 014651 017797 to 802 – 016584 – 016611 to 615
<b>PACKAGING</b>	BLISTER OF 20 (4 X 5) OR 75 (5 X 15) OR 100 (5 X 20)
<b>APPLICATION</b>	SPERM PACKAGING AND CRYOGENIC STORAGE
<b>COUNTRY OF ORIGIN</b>	FRANCE
<b>MANUFACTURER</b>	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
<b>MEDICAL DEVICE</b>	CE CLASS 2A, STERILE
<b>STERILIZED</b>	YES (STERILIZATION BY IRRADIATION)
<b>SHELF LIFE</b>	3 YEARS

### COMPOSITION/INFORMATION ON THE COMPONENTS

Tube = Ionomeric resin, biocompatible, non toxic  
Plug = Polyester viscose + polymerizing powder

### SPECIFICATIONS

Length = +/- 133 mm  
Diameter (outer) = +/- 2.5 mm


### BIOCOMPATIBILITY

MEA (Mouse Embryo Assay) and LAL (endotoxin) tested

### MANIPULATION AND STORAGE/STABILITY

Do not use if the product or the packaging is damaged as sterility and integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp

  
Cryo Bio System  
Groupe IMV Technologies  
27, rue Taitbout  
75009 PARIS - FRANCE  
Tél. +33 (0) 149 240 505  
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY