

PRODUCT SHEET

IDENTIFICATION OF PRODUCT

NAME OF PRODUCT	VISOTUBES
REF. NUMBER	020470, -72, -020482 – 020499, 020500 – 020523, 022141
PACKAGING	BAGS OF 10
APPLICATION	STORAGE COMPARTMENTS FOR GOBLETTS FOR ALL TYPES OF STRAWS
COUNTRY OF ORIGIN	FRANCE
MANUFACTURER	CRYO BIO SYSTEM, GROUP IMV TECHNOLOGIES
MEDICAL DEVICE	No, ACCESSORY
STERILIZED	No
SHELF LIFE	ND

COMPOSITION/INFORMATION ON THE COMPONENTS

LDPE (Low Density PolyEthylene)

SPECIFICATIONS

Height: 120 / 132mm

Diameter: 7.1 - 13 mm

BIOCOMPATIBILITY

ND

MANIPULATION AND STORAGE/STABILITY

Do not use if the product or the packaging is damaged as integrity of the product can no longer be guaranteed

Storage cautions: Preserve away from light, heat and damp

Cryo Bio System
Groupe IMV Technologies
27, rue Taitbout
75009 PARIS - FRANCE
Tél. +33 (0) 149 240 505
www.cryobiosystem.com

CRYO BIO SYSTEM IS AN ISO 13485:2003 CERTIFIED COMPANY