

PBS DULBECCO'S	LIQUID		
	PBS 1X with Ca & Mg	PBS 1X without Ca & Mg	PBS 10X without Ca & Mg
<i>Catalog Nos.</i>	9236	9240	9242 (10X)
<i>Component</i>	<i>mg/L</i>	<i>mg/L</i>	<i>mg/L</i>
Hydrochloric Acid 5N	—	—	—
Potassium Chloride	200	200	2000
Potassium Phosphate monobasic KH_2PO_4	200	200	2000
Sodium Chloride	8000	8000	80000
Sodium Phosphate, dibasic Na_2HPO_4 anhyd.	1150	1150	11500
Magnesium Chloride•6H ₂ O $\text{MgCl}_2\cdot 6\text{H}_2\text{O}$	100	—	—
Calcium Chloride, anhyd.	100	—	—